**MELISSA**

**TECHNICAL NOTE 96.1**

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# ***TECHNICAL NOTE 96.1***

## **Results of Chamber Reassembly at UAB**

Prepared by/Préparé par	Stasiak, M. and Peiro. E.
Reference/Référence	MELISSA Pilot Plant Frame Contract 19445/05/NL/CP
Issue/Edition	0
Revision/Révision	0
Date of issue/Date d'édition	30/09/10
Status/Statut	Final

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# MELiSSA

## List of acronyms

MELiSSA: Micro-Ecological Life Support System Alternative

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HPC: Higher plants chamber

PID: Pipelines and instrumentation diagram

VFD: Variable frequency drive

## TABLE OF CONTENTS

1. Introduction .....	5
2. Applicable Documents .....	5
2. Reference Documents .....	5
3. Tasks distribution .....	5
4. Safety aspects .....	7
5. Description of the installation tasks .....	7
5.1 Delivery and visual inspection (November 2008) .....	7
5.2 Reassembly (December 2008 - February 2009) .....	7
5.3 Installation inside the MPP (January 2009) .....	8
5.4 Initial Functional Testing .....	9
5.5 Wiring Cabinet Layout update .....	12
5.6 Review of the Corrosion status .....	13
5.7 Results of Assembly .....	13
6. Additional hardware modifications (July 09) .....	13
7. Comments .....	15
8. Appended documents .....	17
- Appendix 1: UAB inspection report .....	17
- Appendix 2: Layout diagrams .....	17
- Appendix 3: PID diagrams .....	17
- Appendix 4: Electrical drawings of the power cabinet .....	17
- Appendix 5: Lamp and ballasts technical documentation .....	17
- Appendix 6: CIFA Report on hot water circuit installation and hydraulic testing .....	17
- Appendix 7: Wiring Cabinet Layout .....	18
- Appendix 8: Argus Wiring .....	18
- Appendix 9: Gas injection Wiring Layout .....	18
- Appendix 10: Corrosion report .....	18
- Appendix 11: Hardware modification drawings .....	18
- Appendix 12: Gas analyser technical documentation .....	18
- Appendix 13: Hydroponics pump and VFD technical documentation .....	18

### 1. Introduction

The present document describes the activities of delivery, reassembly and installation of the first higher plant chamber installed in the MPP.

The following are the main tasks performed within this installation:

- Preparation of Utilities lines in HPC room of the MELiSSA Pilot Plant (9D area)
- Installation of lamp ballasts and power cabinet
- Delivery and location of HPC modules in the room.
- Re-assembly of modules and airlocks
- Re-assembly on the chamber of parts disassembled before shipment to UAB.
- Connection of HPC interfaces with the different supply lines (all external interfaces): power, chilled and hot water, auxiliary gases, etc.
- Mounting of additional elements in the temperature control loop
- Checking of hardware configuration: internal connections, connections to Utilities, Utilities lines and equipment, and individual equipment installation.
- Checking of software configuration: PLC and Supervision / HMI.
- Hardware additional modifications

### 2. Applicable Documents

Ref.	Title	Reference	Issue	Date
AD1	MPP Quality manual	MPP-QA-07-0001	1	20/03/08
AD2	MPP rules for good lab practices	MPP-QA-07-0003	0	11/12/07

### 2. Reference Documents

Ref.	Title	Reference	Issue	Date
RD1	MELiSSA Pilot Plant: General Resources, Interfaces and Environment	MPP-TN-08-0001	0	01/04/08
RD2	Higher Plant Chamber Prototype for the MELiSSA Pilot Plant: Detailed Design and Verification	TN85.5		02/11/06
RD3	Prototype Operations Manual	TN85.91	1	August 2008
RD4	Prototype Interface Specifications User Manual	TN85.73	1	May 2008

### 3. Tasks distribution

Responsibilities for the development of the tasks are distributed among the different parts involved in this TN, as follows:

This document is confidential property of the MELiSSA partners and shall not be used, duplicated, modified or transmitted without their authorization

### MPP(UAB)

- Inspects visually the hardware delivered and unpacked with UoG
- Defines the location of chamber and Utilities in HPC room of the MELiSSA Pilot Plant.
- Defines the specification of Utilities lines and equipment, and interface connections to the HPC.
- Reconnects the Argus controller with UoG.
- Supervises the checking of the hardware and software installation.
- Supervises the performance of the hardware modifications defined after the Functional Tests.
- Stores and controls all documentation to maintain its integrity.
- Defines hardware additional modifications
- 

### UoG/ANGSTROM

- Unpacks the crated modules
- Inspects visually the hardware delivered and unpacked with UAB
- Carries out the Checking of the hardware and software installation.
- Locates HPC modules into the room and connects them to each other (internal interfaces).
- Re-assembles HPC parts disassembled before shipment to UAB.
- Performs functional checks of the instruments that were dismantled.
- Performs checking of hardware configuration: internal connections, connections to Utilities, Utilities lines and equipment, and individual equipment installation.
- Connects Argus PLC to the chamber with UAB.
- Checking of software configuration: Argus PLC and Supervision / HMI.
- Defines a list of spare parts, items to be procured.
- Trains MPP, SHERPA and NTE staff for HPC operation, HPC maintenance and ancillary equipment use.

### CIFA / SATE

- Prepares Utilities lines in HPC room.
- Installs of lamp ballasts and mounts and wires lamps
- Installs and connects the power interface cabinet (EPIC)
- Connects the Argus control system to power and wiring to the chamber with UoG
- Connects the chilled/hot water lines and proportional valves and recirculation pumps
- Insulates the chilled/hot water lines
- Connects compressed air supply
- Performs the approved hardware modifications

### SHERPA

- Checking of software configuration: Schneider PLC and Supervision / HMI.

### NTE

- Connects Schneider PLC to the chamber with UAB.
- Checking of software configuration: Schneider PLC and Supervision / HMI.

## CARBUROS METALICOS

- Connects auxiliary gases supply: CO<sub>2</sub>, N<sub>2</sub>, calibration mixture for gas analyzer, and airlock purge gas.

## 4. Safety aspects

The tasks reported in this TN were performed following the safety rules established for the operation in the MPP, particularly the following:

- All people involved in the installation activity had to know and observe the rules related to safety required by MPP and UAB (see AD2).
- Any specific precaution and rule that must be observed in the activities related to the execution of these tasks had to be pre-determined by the MPP Safety Manager or, where feasible, detailed in the relevant specification documents.
- A Safety check-list had to be used and compiled prior to the initiation of the installation tasks, considering the acceptance of the critical parameters there indicated as blocking for the installation phase.

## 5. Description of the installation tasks

### 5.1 Delivery and visual inspection (November 2008)

The higher plant chamber (HPC) built by the University of Guelph with Angstrom Engineering was shipped on November 14, 2008 and delivered to the MELISSA Pilot Plant on November 25, 2008. UAB inspected the shipped subsections and accompanying support packages for content and damage. The visual inspection report is appended (Appendix 1).

UoG confirmed that there was no deflector installed inside the plenum between the blower outlet and the louvers air grids. This point could not be visually inspected by UAB in November 2008 (but later on it was indeed checked).

### 5.2 Reassembly (December 2008 - February 2009)

The HPC was reassembled in three one week stages as follows:

December 15 – 19, 2008	Reassembly of main shell components
January 19 – 23,	Argus control system

2009	Hydroponics system Gas sampling and injection lines Gas and purge system solenoid valves Control system wiring (80%) Initial subsystem testing
February 2 – 6, 2009	Completion of control system wiring Complete system testing (initial testing)

The layout of the HPC within the MPP facility after the reassembly is included in Appendix 2.

### 5.3 Installation inside the MPP (January 2009)

The following table summarizes the activities performed in the MPP in order to interface with the HPC1. These tasks were performed mainly by the company CIFA (otherwise specified):

January 2009	Preparation of Hot and chilled water circuits and equipment in the MPP
January 2009	Electrical power supply for the lamps ballasts, including the power cabinet installation
January 2009	Mounting and connection of valves and forced circulation pumps on chilled and hot water loops
January 2009	Hot and chilled water connection to the chamber
January 2009	Insulation of chilled/hot water lines
January 2009	Connection of Argus control system to power and wiring
January 2009	Supply of compressed air
February 2009	Supply of calibration gas, nitrogen, CO2 and airlock purge gas (by CARBUROS METALICOS)

The Layouts of the corresponding Utilities and the ballasts layout are shown in Appendix 2 (HPC1\_Layout\_diagrams.pdf file).

The PID diagrams of the HPC and the specific Utilities PID are shown in Appendix 3 (HPC1\_PIDs.pdf file; HPC1 PID document ref.: MPP-4100-A-001-A1, dated 15/06/09; Utilities document PID ref.: Cooling water loop diagram, dated 30/09/10).

The electrical cabinet drawings are included in Appendix 4 (HPC1\_Electrical\_cabinet\_drawings.pdf file; CIFA document, dated 10/12/08)

The technical documentation of the lamps and ballasts is included in Appendix 5 (HPC1\_Lamps\_ballasts\_info.pdf file).

The installation and hydraulic testing report corresponding to the hot water circuit performed by the company CIFA is annexed to the present document as Appendix 6 (HPC1\_Hot\_water\_install\_testing.pdf, CIFA document, dated 30/09/10).

### 5.4 Initial Functional Testing

The following table documents the initial input/output connection and testing of the HPC.

ZS_4100_01	Airlock A upper door switch	<i>Connected</i>	<i>Functional</i>
ZS_4100_02	Airlock A lower door switch	<i>Connected</i>	<i>Functional</i>
ZI_4100_01	Airlock A door indicator LEDs	<i>Connected</i>	<i>Functional</i>
RV_4100_01	Airlock A passive vent	not used	As in UoG
ZS_4101_01	Airlock C upper door switch	<i>Connected</i>	<i>Functional</i>
ZS_4101_02	Airlock C lower door switch	<i>Connected</i>	<i>Functional</i>
ZI_4101_01	Airlock C door indicator LEDs	<i>Connected</i>	<i>Functional</i>
RV_4101_01	Airlock C passive vent	not used	As in UoG
SV_4102_01	Airlock A gas injection solenoid valve	<i>Connected</i>	<i>Functional</i>
SV_4102_01	Airlock A gas ventilation solenoid valve	<i>Connected</i>	<i>Functional</i>
PT_4102_01	Airlock A pressure gauge	<i>Connected</i>	<i>Functional</i>
PS_4102_01	Airlock A over pressure switch	<i>Connected</i>	<i>Functional</i>
HV_4102_01	Airlock A purge override manual valve	<i>Connected</i>	<i>Functional</i>
SV_4103_01	Airlock C gas injection solenoid valve	<i>Connected</i>	<i>Functional</i>
SV_4103_02	Airlock C gas ventilation solenoid valve	<i>Connected</i>	<i>Functional</i>
PT_4103_01	Airlock C pressure gauge Reassigned to measure external pressure	<i>Connected</i>	<i>Functional</i>
PS_4103_01	Airlock C over pressure switch	<i>Connected</i>	<i>Functional</i>
HV_4103_01	Airlock C purge override manual valve	<i>Connected</i>	<i>Functional</i>
IY_4104_01	Lamp string A	<i>Connected</i>	<i>Functional</i>
IY_4104_02	Lamp string B	<i>Connected</i>	<i>Functional</i>
IY_4104_03	Lamp string C	<i>Connected</i>	<i>Functional</i>
RT_4104_01	Module A PAR sensor	<i>Connected</i>	<i>Functional</i>
RT_4104_02	Module B PAR sensor	<i>Connected</i>	<i>Functional</i>
RT_4104_03	Module C PAR sensor	<i>Connected</i>	<i>Functional</i>

FSL_4105_01	Lamp loft A fan flow signal	<i>Connected</i>	<i>Functional</i>
FSL_4105_02	Lamp loft B fan flow signal	<i>Connected</i>	<i>Functional</i>
FSL_4105_03	Lamp loft C fan flow signal	<i>Connected</i>	<i>Functional</i>
FAN_4105_01	Lamp loft A fans	<i>Connected</i>	<i>Functional</i>
FAN_4105_02	Lamp loft B fans	<i>Connected</i>	<i>Functional</i>
FAN_4105_03	Lamp loft C fans	<i>Connected</i>	<i>Functional</i>
TT_4105_01	Lamp loft A temperature sensor	<i>Connected</i>	<i>Functional</i>
TT_4105_02	Lamp loft B temperature sensor	<i>Connected</i>	<i>Functional</i>
TT_4105_03	Lamp loft C temperature sensor	<i>Connected</i>	<i>Functional</i>
GP_4106_01	Irrigation pump	<i>Connected</i>	<i>Functional</i>
FT_4106_01	Nutrient flow in sensor	<i>Connected</i>	<i>Functional</i>
HV_4106_01	Irrigation manual override valve A	<i>Connected</i>	<i>Functional</i>
HV_4106_02	Irrigation manual override valve B	<i>Connected</i>	<i>Functional</i>
HV_4106_03	Irrigation drain manual override valve	<i>Connected</i>	<i>Functional</i>
HV_4106_04	Irrigation bypass manual override valve	<i>Connected</i>	<i>Functional</i>
AT_4107_01	pH sensor	<i>Connected</i>	<i>Functional</i>
SV_4107_01	Acid tank solenoid valve	<i>Connected</i>	<i>Functional</i>
SV_4107_02	Base tank solenoid valve	<i>Connected</i>	<i>Functional</i>
LSL_4107_01	Acid tank low level switch	<i>Connected</i>	<i>Functional</i>
LSH_4107_01	Acid tank high level switch	not used	As in UoG
LSL_4107_02	Base tank low level switch	<i>Connected</i>	<i>Functional</i>
LSH_4107_02	Base tank high level switch	not used	As in UoG
HV_4107_01	Acid manual override valve	<i>Connected</i>	<i>Functional</i>
HV_4107_02	Base manual override valve	<i>Connected</i>	<i>Functional</i>
AT_4108_01	EC sensor	<i>Connected</i>	<i>Functional</i>
SV_4108_01	Stock A solenoid valve	<i>Connected</i>	<i>Functional</i>
SV_4108_02	Stock B solenoid valve	<i>Connected</i>	<i>Functional</i>
LSL_4108_01	Stock A low level switch	<i>Connected</i>	<i>Functional</i>
LSH_4108_01	Stock A high level switch	not used	As in UoG
LSL_4108_02	Stock B low level switch	<i>Connected</i>	<i>Functional</i>
LSH_4108_02	Stock B high level switch	not used	As in UoG
HV_4108_01	Stock A manual override valve	<i>Connected</i>	<i>Functional</i>
HV_4108_02	Stock B manual override valve	<i>Connected</i>	<i>Functional</i>

TT_4109_01	Nutrient temperature sensor	<i>Connected</i>	<i>Functional</i>
SV_4109_01	Nutrient cooling solenoid valve	not used	As in UoG
LSH_4110_01	Nutrient tank high level switch	<i>Connected</i>	<i>Functional</i>
LSL_4110_01	Nutrient tank low level switch	<i>Connected</i>	<i>Functional</i>
LSH_4110_02	Condensate tank low level switch	<i>Connected</i>	<i>Functional</i>
LSL_4110_02	Condensate tank high level switch	<i>Connected</i>	<i>Functional</i>
GP_4110_01	Condensate pump relay	<i>Connected</i>	<i>Functional</i>
BLWR_4111_01	Blower	<i>Connected</i>	<i>Functional</i>
FT_4111_01	Air velocity sensor	<i>Connected</i>	<i>Functional</i>
MVFD_4111_01	Blower motor VFD	<i>Connected</i>	<i>Functional</i>
TT_4112_01	Module A temperature sensor 1	<i>Connected</i>	<i>Functional</i>
AT_4112_01	Module A relative humidity sensor	<i>Connected</i>	<i>Functional</i>
TT_4112_02	Module B temperature sensor 1	<i>Connected</i>	<i>Functional</i>
AT_4112_02	Module B relative humidity sensor	<i>Connected</i>	<i>Functional</i>
TT_4112_03	Module C temperature sensor 1	<i>Connected</i>	<i>Functional</i>
AT_4112_03	Module C relative humidity sensor	<i>Connected</i>	<i>Functional</i>
TT_4112_04	Module A temperature sensor 2	<i>Connected</i>	<i>Functional</i>
TT_4112_05	Module A temperature sensor 3	<i>Connected</i>	<i>Functional</i>
TT_4112_06	Module A temperature sensor 4 Sensor reassigned to measure external T	<i>Connected</i>	<i>Functional</i>
TT_4112_07	Module B temperature sensor 2	<i>Connected</i>	<i>Functional</i>
TT_4112_08	Module B temperature sensor 3	<i>Connected</i>	<i>Functional</i>
TT_4112_09	Module B temperature sensor 4	<i>Connected</i>	<i>Functional</i>
TT_4112_10	Module C temperature sensor 2	<i>Connected</i>	<i>Functional</i>
TT_4112_11	Module C temperature sensor 3	<i>Connected</i>	<i>Functional</i>
TT_4112_12	Module C temperature sensor 4 Sensor reassigned to measure external T	<i>Connected</i>	<i>Functional</i>
TT_4112_13	Chilled water source temperature	<i>Connected</i>	<i>Functional</i>
TT_4112_14	Hot water source temperature	<i>Connected</i>	<i>Functional</i>
TT_4112_15	Chilled coil surface temperature 1	<i>Connected</i>	<i>Functional</i>
TT_4112_16	Hot coil surface temperature 1	<i>Connected</i>	<i>Functional</i>
TT_4112_17	Chilled water exit temperature	<i>Connected</i>	<i>Functional</i>
TT_4112_18	Hot water exit temperature	<i>Connected</i>	<i>Functional</i>
TT_4112_19	Chilled coil surface temperature 2	<i>Connected</i>	<i>Functional</i>

TT_4112_20	Hot coil surface temperature 2	<i>Connected</i>	<i>Functional</i>
TT_4112_21	Chilled water inlet temperature	<i>Connected</i>	<i>Functional</i>
TT_4112_22	Hot water inlet temperature	<i>Connected</i>	<i>Functional</i>
S3CV_4112_01	Chilled water proportional valve	<i>Connected</i>	<i>Functional</i>
S3CV_4112_02	Hot water proportional valve	<i>Connected</i>	<i>Functional</i>
FS_4114_01	Vent detect flow switch	<i>Connected</i>	<i>Functional</i>
FC_4113_01	CO2 MFM	<i>Connected</i>	<i>Functional</i>
FC_4113_01	CO2 MFM set-point	<i>Connected</i>	<i>Functional</i>
AT_4113_01	CO2 analyser *	Connected	Functional
AT_4113_02	O2 analyser *	Not yet available	Not yet available
SV_4113_01	CO2 injection solenoid valve	<i>Connected</i>	<i>Functional</i>
PT_4114_01	Growing area pressure sensor (module B)	<i>Connected</i>	<i>Functional</i>
RV_4114_01	Chamber pressure relief valve	<i>Connected</i>	<i>Functional</i>
	Hydroponics system plumbing/return	<i>No leaks / balanced</i>	<i>Functional</i>
	CO2 injection plumbing	<i>No leaks</i>	<i>Functional</i>
	Purge plumbing	<i>No leaks</i>	<i>Functional</i>
	Acid/Base A/B reservoirs	<i>No leaks</i>	<i>Functional</i>
	Hydroponics Reservoir	<i>No leaks</i>	<i>Functional</i>
	Condensate Reservoir	<i>No leaks</i>	<i>Functional</i>
	Heat exchangers	<i>No leaks</i>	<i>Functional</i>
	Expansion bladders	<i>Installed</i>	<i>Functional</i>
	Outer airlock doors		<i>Functional</i>
	Inner airlock doors		<i>Functional</i>
	Roller system for trays		<i>Functional</i>
	T/RH control		<i>Functional</i>

\* A VAISALA GMP343 CO<sub>2</sub> analyser was temporarily installed until the definite equipment (combined O<sub>2</sub> and CO<sub>2</sub> analyser) was received (see Section 6).

## 5.5 Wiring Cabinet Layout update

The intermediate box colour coding layout was updated by J. Lawson and is appended as the excel document 'wiring\_cabinet\_jamie.pdf' (UoG document, Appendix 7). The Argus system

wiring layout is also appended as the PDF file UAB\_Argus\_Wiring.pdf (UoG document, dated 06/02/09, Appendix 8). The wiring layout of the gas injection (CO<sub>2</sub> and airlock purge gas) is appended as HPC1\_Gas\_Panel.pdf (UoG document, Appendix 9).

### 5.6 Review of the Corrosion status

Initial inspection by Angstrom and UoGuelph personnel have not noted a recurrence of the corrosion found on the HPC at the UoGuelph facility prior to shipment: the corrosion which had appeared on most interior surfaces during initial testing, had not reappeared after the complete tear down and passivation treatment by Angstrom Engineering.

The report by Bodycote for Angstrom engineering, Angstrom\_Bodycote\_report.pdf is appended (original document "Investigation of corrosion of a 316SS growth chamber, by B. Elliot and C. Fleck, BODYCOTE TESTING GROUP, Ref. G808171, dated 21/10/08, Appendix 10)

Anyway, as indicated in the UAB inspection report (Appendix 1), the previous impact on the stainless steel is still visible on the chamber modules although the current status is OK (evolution to be checked periodically).

It was decided in agreement with ESA that UAB and UoG would define the needs for a service on corrosion follow-up, focused especially on the prevention approach, and ask specialised companies to quote for this service.

### 5.7 Results of Assembly

By the end of the first week of February 2009, the HPC was completely reassembled, tested, and ready for subsequent functional testing.

## 6. Additional hardware modifications (July 09)

Some upgrade activities were implemented in the hardware and some maintenance tasks as well performed before the replacement of the Argus control system by the Schneider one took place, as follows:

- Plastic supports manufactured by CIFA were installed by MPP personnel along the chamber by UoG personnel to both avoiding the movement of the taps, and levelling the same so that guaranteeing the correct flow of hydroponic solution to the trays. However, in order to provide better stability and long-term resistance, they were later on replaced by stainless steel supports (draft design drawings shown in Appendix 11a).
- A new electric cabinet was installed by MPP personnel to fit inside the electro valves for CO<sub>2</sub> injection and airlock purging.

- A new CO<sub>2</sub>/O<sub>2</sub> gas analyser CAI 601P (California Analytical Instruments, supplied in Spain by the company EQUITROL) was installed, including electrical connection and piping. The technical datasheet and User's manual of these equipment (California Analytical Instruments commercial info) are included in Appendix 12
- Panels closing the lower part of the chamber were modified to allow tubing and equipment to be fitted through them.
- Two Temp.-Humidity sensors were replaced by new ones due to potential failure (the failure was later confirmed).
- The hydroponic Nutrient Flow Meter was replaced by a new one 4-20 mA provided by UoG, and connected electrically.
- A leak was repaired (twice) in the intermediate connection of the general polypropylene collector of the trays (caused due to insufficient silicon welding). This failure was recorded in a Non Conformance Report (MPP-NCR-09-4104)
- Substitution of one of the Lighting Loft fans, that failed and didn't send feedback signal to the PLC.
- New stainless steel supports for the PAR sensors were designed and installed, to replace the former temporary supports (draft drawing in Appendix 11b).

## 7. Comments

### *Results of Chamber Reassembly at UAB*

#### *Comments*

#### Detailed comments

Page/paragraph	Comment
7/ Section 5.1	<p>A word is missing after UoG</p> <p>OK, completed: <i>“UoG confirmed ...”</i></p>
8/ Section 5.2	<p>The precise order of perforated panels for air flow distribution is not recalled anywhere. Was it identified at that time or not? Please precise</p> <p><i>At that time, the precise order of perforated panels was not identified.</i></p>
13/ Section 6	<p><i>Part of these tasks were performed along WP96.3 duration, but they're not related to the Argus/ Schneider replacement, so I propose to keep them here</i></p> <p>More or less the need for these tasks has been identified during the re-assembly, so indeed I would keep them here</p>
14/ Section 6, end of first paragraph and beginning of second paragraph	<p><i>“New stainless steel supports for the PAR sensors were designed and installed, to replace the former temporary supports”</i></p> <p><i>This task for example was performed later (September 09) after Argus removal, but nothing to do with the control system replacement</i></p> <p><i>“Some points were highlighted as potential for future improvement in the hardware: ...”</i></p> <p><i>To be kept here or moved to TN96.13?</i></p> <p>If these points were identified during the re-assembly, then you can keep them here and recall them in TN 96.13. If they have been identified during the course of the testing ,then only in 96.13</p> <p><i>They were mainly identified during the tests, so then removed here and to be traced in TN96.13</i></p>
14/ Section 7, Title	<p><i>“Additional installation tasks performed in the frame of the WP-96.3”</i></p> <p><i>In principle we consider these out of this TN as this installation is related to the control system performance, but the second one was done before removing the Argus, so to keep it?</i></p> <p>I would keep 7, 8 and 9 in TN 96.13 only, as they were identified/performed along</p>

	<p>the testing phase</p> <p>OK, removed here and to be traced in TN96.13</p>
14/ Section 8, Title	<p><i>"Additional maintenance tasks due to hardware failure"</i></p> <p>Even if maintenance tasks, they are related to the quality of the initial installation, so I don't know if keeping here or include in the TN96.13</p> <p>To be moved to TN96.13 as agreed in previous comment</p>
14/ Section9, Title	<p><i>"Final hardware modifications"</i></p> <p>Within these tasks some are really new installation even if performed late in the COO6, others are maintenance, so same question as previous comment in Section 8</p> <p>To be moved as well to TN96.13</p>
14/ section 9, last bullet	<p><i>"The PLC capacity was extended to allow keeping all sensors implemented on the HPC by the addition of Remote I/O"</i></p> <p>Probably to be removed (traced in other documents)</p> <p>You can anyway keep it here, even recalling the ref of the documents describing these tasks in detail</p>

## 8. Appended documents

### - Appendix 1: UAB inspection report

Electronic document:

MPP-REP 08-4101\_0\_ HPC1 Visual inspection acceptance Report.pdf

### - Appendix 2: Layout diagrams

Electronic document: HPC1\_Layout\_diagrams.pdf

### - Appendix 3: PID diagrams

Electronic document: HPC1\_PIDs.pdf (HPC1 PID document ref.: MPP-4100-A-001-A1, dated 15/06/09; Utilities document PID ref.: Cooling water loop diagram, dated 30/09/10).

### - Appendix 4: Electrical drawings of the power cabinet

Electronic document: HPC1\_Electrical\_cabinet\_drawings.pdf (CIFA document, dated 10/12/08)

### - Appendix 5: Lamp and ballasts technical documentation

Electronic document: HPC1\_Lamps\_ballasts.pdf

### - Appendix 6: CIFA Report on hot water circuit installation and hydraulic testing

Electronic document: HPC1\_Hot\_water\_instal\_testing.pdf (CIFA document, dated 30/09/10)

## **- Appendix 7: Wiring Cabinet Layout**

Electronic document: HPC1\_wiring\_cabinet\_jamie.pdf (UoG document)

## **- Appendix 8: Argus Wiring**

Electronic document: HPC1\_\_Argus\_Wiring.pdf (UoG document, dated 06/02/09)

## **- Appendix 9: Gas injection Wiring Layout**

Electronic document: HPC1\_Gas\_Panel.pdf (UoG document)

## **- Appendix 10: Corrosion report**

Electronic document: HPC1\_Angstrom\_Bodycote\_report.pdf (original document "Investigation of corrosion of a 316SS growth chamber, by B. Elliot and C. Fleck, BODYCOTE TESTING GROUP, Ref. G808171, dated 21/10/08)

## **- Appendix 11: Hardware modification drawings**

Electronic documents (draft design drawings): HPC1\_PAR\_sensors\_support.pdf (Appendix 11a) and HPC1\_Spigots\_support.pdf (Appendix 11b)

## **- Appendix 12: Gas analyser technical documentation**

Electronic document: HPC1\_CAI600\_technical\_info.pdf (California Analytical Instruments commercial info)

## **- Appendix 13: Hydroponics pump and VFD technical documentation**

Electronic document: HPC1\_Hydroponics\_pump\_VFD.pdf (ESPA commercial info)

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## HPC1 Visual inspection acceptance Report

Approval Loop :

Date	Issue	Prepared by (visa):	Checked by (visa):	Approved by (visa):
01/12/2008	0	E. PEIRO	A. FOSSEN	F. GODIA

Distribution List :

Company	Name	Number of copies
ESA/ESTEC	B. LAMAZE	1
	C. PAILLE	1
UOG	M. DIXON	1

This document is confidential property of the MELiSSA partners and shall not be used, duplicated, modified or transmitted without their prior authorization.

# MELiSS Pilot Plant

Document Identification : HPC1 Visual inspection acceptance Report	Type	Chrono	Issue	Page: 2 / 16
	MPP-REP	08-4101	(0)	

Change log :

Date	Issue	Reason of the change	Modified paragraphs
01/12/2008	0	Creation	

This document is confidential property of the MELiSSA partners and shall not be used, duplicated, modified or transmitted without their prior authorization.

Document Identification : HPC1 Visual inspection acceptance Report	Type	Chrono	Issue	Page: 3 / 16
	MPP-REP	08-4101	(0)	

## TABLE OF CONTENTS

1. Purpose.....	3
2. Delivery and visual inspection.....	3
3. Acceptance of the equipment.....	3
4. Appendix 1: Visual inspection check-list .....	5
5. Appendix 2: Pictures of HPC1 after arrival.....	14

### 1. Purpose

The objective of the present document is to describe the delivery of the Higher plant Chamber HPC1 to the MELiSSA Pilot Plant at UAB and to explain the results of the visual inspection and checking performed by MPP and UoG for acceptance of the delivered items.

### 2. Delivery and visual inspection

The equipment was delivered on Tuesday 25 November 2008 by truck by SERGITRANS, company subcontracted by ALBATRANS, representative of FEDEX. The 5 crates were unloaded by means of a pallet truck and the truck lifting platform. No forklift had been planned by SERGITRANS, therefore it could not be used. The crates were opened on the loading dock and the equipment was then moved inside the Pilot Plant.

All the components were visually inspected and checked jointly by MPP and UoG, except for the box containing Hidroponic system and Chamber shell parts, which were visually inspected in general but not checked individually, because they were not broken down into a detailed packing list.

The separate smaller boxes also included in the crates were also opened to check the status of the hardware.

Nothing was found damaged by the transport. The remarks made during the visual inspection are summarized in the attached table (see 4. Appendix: Visual inspection check-list).

### 3. Acceptance of the equipment

Based in the visual review performed by MPP and UoG, the delivered equipment is accepted, with the following main remarks:

This document is confidential property of the MELiSSA partners and shall not be used, duplicated, modified or transmitted without their prior authorization.

# MELiSS Pilot Plant

Document Identification : HPC1 Visual inspection acceptance Report	Type	Chrono	Issue	Page: 4 / 16
	MPP-REP	08-4101	(0)	

- Crate including Hydroponic system and Chamber shell parts: parts not checked individually, to be checked individually by UoG when proceeding to the remounting
- Chamber modules (A, B and C): treated for removal of corrosion; the previous impact on the stainless steel is still visible although the current status is OK. Evolution to be checked periodically : a procedure for inspection will be requested to UoG to specify the dismounting steps to access inspection points and the control operations (visual inspection, sampling, measurements) to be performed.

# MELiSS Pilot Plant

Document Identification : HPC1 Visual inspection acceptance Report	Type	Chrono	Issue	Page: 5 / 16
	MPP-REP	08-4101	(0)	

## 4. Appendix 1: Visual inspection check-list

HPC1 Shipping Container Contents						
ITEM	DELIVERY	CONDITION	REMARKS			
			Location	Ext./Int.	Appearance	
LAMP LOFT MODULE A	OK	OK				
LAMP LOFT MODULE B	OK	OK				
LAMP LOFT MODULE C	OK	OK				
MODULE A	OK	See remarks	Airlock door and cage	Int.	"Washing" areas different colour (Appendix 2, Fig. 1 and 2)	
			Airlock	Ext.	Shadows and circles different colour	
			Main module	Int.	"Washing" areas different colour (much less in left side) (Appendix 2, Figs. 3 and 4)	
			Main module	Int.	White small pitting (lower part on both sides) (Appendix 2, Fig. 5)	
			Lower platform	Int.	White, grey and yellow spots (Appendix 2, Fig. 6)	
			Lower platform	Ext.	Yellow powder (oxidation) in edge	
			Lower platform	Int.	White powder on feeding valves support	
			Ventilation ducts	Int.	White small pitting (Appendix 2, Fig. 7)	

This document is confidential property of the MELiSSA partners and shall not be used, duplicated, modified or transmitted without their prior authorization.

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	<b>MODULE B</b>	OK	See remarks	Main module	Int.	White small pitting (lower part on both sides) (Appendix 2, Fig. 8)
				Lower module	Int.	White small pitting (platform, ceiling, wall to module C, steel supports of heat exchanges -partially- and front cover) (Appendix 2, Figs. 9, 10, 11 and 12)
				Lower module	Int.	Oxidation in steel elbows of heat exchangers (Appendix 2, Figs. 13 y 14)
				Ventilation ducts	Int.	White small pitting (Appendix 2, Fig. 15)
	<b>MODULE C</b>	OK	See remarks	Airlock door	Ext.	Some yellow spots (upper side)
				Airlock door	Int.	Shadows different colour (lower part)
				Airlock cage	Ext./Int.	Circles different colour (upper and lower parts)
				Airlock cage	Ext.	Shadows/ black spots (both sides)
				Airlock cage	Ext.	Yellow spots (right side)
				Main module	Int.	Grey spots all over, no oxidation (left side and lower part)
				Main module	Int.	White small pitting (right side)
				Main module	Int.	Yellow spots over white pitting (right side)
				Main module	Int.	White small pitting (lower part on both sides)
				Lower platform	Ext.	Yellow powder (oxidation) in edge (Appendix 2, Fig. 16)
				Lower platform	Int.	Oxidation in pump wheel (Appendix 2, Fig. 17)
				Ventilation ducts	Int.	White small pitting (Appendix 2, Fig. 18)
	<b>BOX WITH HYDROPONIC SYS. AND CHAMBER SHELL PARTS</b>	OK	OK			

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<b>BOX 1</b>	OK	OK			
<b>Box 1A</b>	OK	OK			
1 * Nutrient flow sensor and cable	OK	OK			
<b>Box 1B</b>	OK	OK			
1 * Mass flow meter, power supply and cable	OK	OK			
<b>Box 1C</b>	OK	OK			
Air handling system temperature probes	OK	OK			
5 * wet well temperature sensors	OK	OK			
2 * air thermisters	OK	OK			
<b>Box 1D</b>	OK	OK			
3 * PAR sensors and cables	OK	OK			
<b>Box 1E</b>	OK	OK			
3 * Vaisala temperature / relative humidity sensor ends	OK	OK			
<b>Box 1F</b>	OK	OK			
1 * Honeywell pressure transducer	OK	OK			
<b>Box 1G</b>	OK	OK			
1 * Honeywell pressure transducer	OK	OK			
<b>Box 1H</b>	OK	OK			
1 * Honeywell pressure	OK	OK			

This document is confidential property of the MELiSSA partners and shall not be used, duplicated, modified or transmitted without their prior authorization.

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	transducer				
	<b>Box 1I</b>	OK	OK		
	1 * Anemometer sensor end	OK	OK		
	<b>Box 1J</b>	OK	OK		
	3 * growth area temperature sensor assemblies	OK	OK		
	1 * coil temperature thermister assembly	OK	OK		
	<b>Box 1K</b>	OK	OK		
	1 * irrigation pump controller	OK	OK		
	<b>Box 1L</b>	OK	OK		
	Argus control panel spare components	OK	OK		
	termination resistors, fuses, mounting hardware, special tools	OK	OK		
	<b>Box 2</b>	OK	OK		
	CO2 injection valve fixture	OK	OK		
	5 * ASCO solenoid valves	OK	OK		
	1 * stainless steel mounting plate	OK	OK		
	1 * wiring harness	OK	OK		

This document is confidential property of the MELiSSA partners and shall not be used, duplicated, modified or transmitted without their prior authorization.

# MELiSS Pilot Plant

Document Identification : HPC1 Visual inspection acceptance Report	Type	Chrono	Issue	Page: 9 / 16
	MPP-REP	08-4101	(0)	

<b>Box 3</b>	OK	OK			
1 * pH probe and signal conditioner	OK	OK			
1 * EC probe and signal conditioner	OK	OK			
<b>Box 4</b>	OK	OK			
2 * proportional valves	OK	OK			
<b>Box 5</b>	OK	OK			
1 * irrigation controller to pump power cable assembly	OK	OK			
<b>Box 6</b>	OK	OK			
Nutrient delivery manifold assembly	OK	OK			
1 * stainless steel mounting plate	OK	OK			
4 * ASCO solenoid valve assemblies	OK	OK			
4 * manual shutoff valves	OK	OK			
4 * tubing assemblies	OK	OK			
<b>Box 7</b>	OK	OK			
4 * polypropylene stock tanks	OK	OK			
1 * polypropylene	OK	OK			

This document is confidential property of the MELiSSA partners and shall not be used, duplicated, modified or transmitted without their prior authorization.

# MELiSS Pilot Plant

Document Identification : HPC1 Visual inspection acceptance Report	Type	Chrono	Issue	Page: 10 / 16
	MPP-REP	08-4101	(0)	

	condensate tank					
	associated tubing assemblies	OK	OK			
	<b>Box 8</b>	OK	OK			
	2 * spare polypropylene stock tanks	OK	OK			
	1 * box plastic fittings (various)	OK	OK			
	<b>Box 9</b>	OK	OK			
	spare electrical connectors	OK	OK			
	feed through assemblies	OK	OK			
	feed through gaskets	OK	OK			
	connector bodies	OK	OK			
	pins	OK	OK			
	shrink wrap boots	OK	OK			
	1 pack feed through mounting bolts	OK	OK			
	1 * crimping tool	OK	OK			
	1 * stainless steel feed through mounting plate	OK	OK			
	3 sizes shrink wire wrap	OK	OK			
	1 set labels (blank)	OK	OK			
	1 * air velocity sensor	OK	OK			

This document is confidential property of the MELiSSA partners and shall not be used, duplicated, modified or transmitted without their prior authorization.

# MELiSS Pilot Plant

Document Identification : HPC1 Visual inspection acceptance Report	Type	Chrono	Issue	Page: 11 / 16
	MPP-REP	08-4101	(0)	

hardware					
1 * Argus USB cable	OK	OK			
2 * spare screws	OK	OK			
3 * bags tie wraps	OK	OK			
1 * spare condensate pump	OK	OK			
1 set irrigation pump brushes	OK	OK			
1 * irrigation pump controller resistor	OK	OK			
1 * spare pressure sensor gasket and cable	OK	OK			
2 * optional 250ohm resistors for proportional valves	OK	OK			
<b>Box 10</b>	OK	OK			
irrigation pump	OK	OK			
<b>Box 11</b>	OK	OK			
nutrient tank plumbing	OK	OK			
<b>Box 12</b>	OK	OK			
multiple plastic conduit	OK	OK			
1 * irrigation tray delivery manifold	OK	Dirty			
1 * roll plastic film	OK	OK			

This document is confidential property of the MELiSSA partners and shall not be used, duplicated, modified or transmitted without their prior authorization.

# MELiSS Pilot Plant

Document Identification : HPC1 Visual inspection acceptance Report	Type	Chrono	Issue	Page: 12 / 16
	MPP-REP	08-4101	(0)	

<b>Box 13</b>	OK	OK			
2 * plastic conduit	OK	OK			
<b>Box 14</b>	OK	OK			
multiple plastic tubing lengths	OK	OK			
2 * rolls wire	OK	OK			
2 * rolls wire wrap	OK	OK			
<b>Box 15</b>	OK	OK	Includes lettuce seeds		
2 * growth trays	OK	OK			
2 * ESA stick on signs	OK	OK			
<b>Box 16</b>	OK	OK			
1 * computer w/ keyboard and mouse – documentation- software	OK	OK			
<b>Box 17 (non labelled)</b>	OK	OK			
Connectors for the Teflon bags, removed gaskets, and bolts	OK	OK			
<b>Box 18 (non labelled)</b>	OK	OK			
Teflon bag	OK	OK			
<b>Box 19 (non labelled)</b>	OK	OK			
Teflon bag	OK	OK			

This document is confidential property of the MELiSSA partners and shall not be used, duplicated, modified or transmitted without their prior authorization.

# MELiSS Pilot Plant



Document Identification : HPC1 Visual inspection acceptance Report	Type	Chrono	Issue	Page: 13 / 16
	MPP-REP	08-4101	(0)	

<b>Box 20 (non labelled)</b>	OK	OK			
Teflon bag	OK	OK			
<b>Box 21 (non labelled)</b>	OK	OK			
Reassembly parts (ARMSTRONG)	OK	OK			
<b>Box 22 (non labelled)</b>	OK	OK			
Reassembly parts and spares (ARMSTRONG)	OK	OK			
<b>Box 23 (non labelled)</b>	OK	OK			
Reassembly parts: plates for connecting trays	OK	OK			
<b>Box 24 (non labelled)</b>	OK	OK			
Rockwool	OK	OK			

Document Identification : HPC1 Visual inspection acceptance Report	Type	Chrono	Issue	Page: 14 / 16
	MPP-REP	08-4101	(0)	

## 5. Appendix 2: Pictures of HPC1 after arrival

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








Revision	Author	Date	Checked by	Check date	Approved by	Approval Date	Reason for change
A1	A.Fossen	2009-06-15	R. Moyano	2009-06-26	E. Peiro	2009-06-26	



 <b>UAB</b> Universitat Autònoma de Barcelona		 <b>MELISSA</b> Memorandum of Understanding T20-MET2002/0161/rev.01	
Scale	1/1	Project Title	Higher Plants Compartment 1 – HPC1
		Drawing Title	P&ID
Document Code		MPP-4100-A-001-A1	

A1	RM	30/09/2010					
Edition	Author	Date	Checked by	Check date	Approved by	Approval Date	Reason for change
		<b>JLT Consultant</b> 47 ter route nationale 62580 VIMY - FRANCE					
<small>This document is confidential property of the MELISSA partners and shall not be used, duplicated, modified or transmitted without their authorization          Memorandum of Understanding TOS-MCT/2002/3161/Inv/CL</small>							
Scale	1/1	Project Title		M.P.P Facilities			
		Drawing Title		Cooling water loop diagram			
Document Code		<b>COOLING WATER LOOP</b>					

A1	RM	21/09/2010					
Edition	Author	Date	Checked by	Check date	Approved by	Approval Date	Reason for change
		<b>JLT Consultant</b> 47 ter route nationale 62580 VIMY - FRANCE					
<small>This document is confidential property of the MELISSA partners and shall not be used, duplicated, modified or transmitted without their authorization          Memorandum of Understanding TOS-MCT/2002/3161/Inv/CL</small>							
Scale	1/1	Project Title		M.P.P Facilities			
		Drawing Title		Hot water loop diagram			
		Document Code		<b>MPP-PID-10-9J24-A1</b>			

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---

## **1: General Product Description**

The HSE product line contains two different models. The 600 Watt type and the 1000 Watt type. The electronic driver is suitable for connection to a line voltage of 400 Volt nominal only. In order to achieve optimal illumination of the electronic lamps in these fixtures, they should be operated at high-frequency. This requires the exclusive use of the advanced lamps that have been specifically adapted for use in this driver.

The lamps can be identified through their specific product names (Philips GreenP.T600W EL400V and Philips GreenP.TD1000W EL400V). Over time, continued use of alternate lamps may cause damage to the driver, and is therefore not advisable. System efficiency will also be affected by the application of alternate lamps.

## **2: Technical Specifications of Electronics in HSE 600 and HSE 1000**

Refer to table 1 below for technical specifications.

### **General:**

Lamp capacity	: 600 and / or 1000 Watt HPS GP E lamp
Ballast loss	: 35 Watt in 1000 Watt / 20 Watt in 600 Watt
Weight total	: 600 Watt = 4.3 kg* / 1000 Watt = 4.7 *

\* weight of lamp and reflector included.

HSE systems have been tested in accordance with IEC 598 and meet all relevant requirements. The systems are rain watertight IP23 / class I. The electronics housing is dust- and watertight IP65.

All fixtures should be earthed prior to operation.

### **Filter Coils**

When using HSE (electronic) fixtures, the use of filter coils is no longer required.

## **3. Assembly Instructions**

### **The following parts will be supplied in a complete HSE fixture delivery:**

- HSE assembly unit
- HSE Cover to close off assembly unit
- Assembly kit consisting of: 2x swivel + ring M25  
2x bolt M6 + nut (assembly bracket)
- The HSE Driver, which is the complete electronics housing
- Reflector
- Lamp
- Brackets (if ordered) These can be either standard brackets or client specific brackets

### **3.1 HSE Assembly Unit**

Assembly units are supplied together with separate swivels M25 (plus nut), M6 bolt + nut and separate cover. The assembly unit can be mounted with the M6 bolt, to be placed in the support profile. Several different brackets and eye bolts are available for attachment to different types of profiles. In order to mount the fixture horizontally and in a straight line, it is advisable to attach

them in alternating fashion. This will prevent torsion of the support profile. When using the one-point suspension brackets (i.e. one per fixture), the attachment point (center of gravity) can easily be determined using the scale marked on the assembly unit.

The separate wire ends supplied are fitted with wire core end sleeves and are suitable for use in standard 2.5 to 4 mm<sup>2</sup> wago connectors or comparable products. The wires are not fit for use with 1.5 – 2.5 mm<sup>2</sup> wago connectors. The outer jacket of the feeder cable can be stripped for a maximum of 12 cm. The highest permissible cable diameter is 2x Vmvk 4x4mm<sup>2</sup>. After wiring, the cover of the assembly unit is fastened. Ensure that the length of the condensation loop in the wiring is kept to a minimum and that it does not rest upon the reflector.

### **3.2 Driver Unit HSE Complete**

The driver unit fits into the assembly unit with a special wedge construction. The electrical connection is realized by shoving the driver unit onto the assembly unit. Always disconnect power first. HSE 600W and 1000 Watt systems are fitted with identical assembly units and are therefore interchangeable. Once the driver unit has been positioned, the rotary knob at the back of the de assembly unit should be given a half turn (do not strain). This knob secures the assembly unit onto the driver unit, thus ensuring the correct contact pressure. All HSE fixtures should be locked prior to connection. Consequential damages caused by failure to comply with these directions do not qualify for compensation. (see **figure 1**)

The driver unit does not contain any service parts. The warranty will expire if the unit is opened. Exercise caution while handling the different “heavy” parts on the print, when working with the driver unit. Knocks to or improper use of the system may cause internal damage and should be avoided at all times.

### **4. Lamp**

The correct lamp should be used with the fixture.

HSE 600 Watt : Philips Greenpower-E 600W/400 Volt lamp.

HSE 1000 Watt : Philips Greenpower-E 1000W/400 Volt lamp.

Refer to the de type-sticker on the fixture for the correct lamp type. Use of incorrect lamps may cause damage to fixture components. During lamp replacement, always disconnect the fixture for a safe work situation. The lamp will light up only if the power has been disconnected prior to lamp replacement. Lamps that are replaced without disconnecting the power will not start up, regardless of the current mains voltage at the socket.

In 1000 Watt systems, the lamp has to be inserted into the special lamp bracket. Open the bracket by slightly lifting the clip making sure to guide the bracket during the opening process. Push the lamp with its flat honeycombed surfaces onto the ceramic parts. In order to do this, the H-ends of the lamp have to be pushed between the two lamp holder springs. Execute with precision. Place unfrayed loose wire ends straight into the appropriate slot groove on the ceramic block. (If performed improperly, the connection will either be poor, or cannot be made at all). Then close the bracket. Ensure that the bracket has been fully closed at both sides. Push down the bracket until an audible click is produced. Faulty or incomplete bracket closure will cause the contact to be burned in the socket. (**figure 2**)

To prevent lamp breakage, do not to operate the fixture in conditions of dense vapor or direct spraying. Also allow the installation to cool for at least 15 minutes to observe the cooling down

phase of the lamp, before bringing it back into use. The lamp can be cleaned with a damp cloth. This applies to both 600 and 1000W lamp types.

## **5: HSE Electronic Driver Unit / Service Status**

The electronic driver unit of the HSE fixtures has been specifically designed for horticultural purposes. It combines low energy consumption with a long life span. Any problems that may occur during use, will be displayed on the LED. For error codes, refer to the table below.

Status	LED	Description	Cause / action	Reset by:
Starter busy	Continuous blinking	Starter is busy	Lamp or ballast cools down until ballast is capacitated to restart	Not necessary
EOL timer expired	1 blink	Ballast switched off	Defective lamp / replace	Reset supply voltage
Lamp is cycling	2 blinks	Number of start up attempts exceeded	Lamp will not start up / replace	Reset supply voltage
Low supply voltage	3 blinks	Supply voltage is low	Lamp is switched off	Bring supply voltage within tolerance
High temperature	4 blinks	Ballast temperature high (max 115 degrees)	Lamp is switched off	Reset supply voltage
Lamp voltage outside parameters	5 blinks	EOL lamp *	Lamp is switched off / replace	Reset supply voltage

Error codes are displayed at an interval of approximately 4 seconds.

- = EOL  
( = end of lamp life)

### **5.1: Initial Start-Up with New Lamp**

After assembly (and / or lamp replacement), the electronics will signal the presence of a new lamp during the initial start up of the fixture. During the first few minutes of initial start up of the system, a number of lamps may not light up immediately. Since the lamp has to burn in, the electronics will automatically reduce the power fed to the lamps to a level that is below nominal specifications. They will continue to do so for the first 100 hours and this ensures that both lamp and electronics are protected from high currents generated in new lamps.

## **6: Maintenance / Safety:**

The HSE driver unit is basically maintenance free.

The electronics are equipped with 2 fuses. These have been fitted for safety reasons (fire safety) rather than for service purposes. The driver has no service parts and the electronics should not be subjected to any attempts at repair or otherwise. The combination of high-frequency and high

capacity is a potential hazard to (un)qualified personnel. Opening the housing will therefore cause the warranty to expire (see warranty conditions).

### **6.1: Assembly Instructions Reflector**

The following applies to the HSE 600 system; the reflector is attached to the housing through an adapter, which is located in the recess over the lamp holder. The adapter of the reflector is located beneath the spring that holds the reflector in place. The lamp should be screwed in place with some caution in order to ensure the correct lamp position. Lamps that have not been correctly placed will have an adverse effect on the light distribution. Screw in the lamp before placing the reflector at all times. (also refer to light plan notes)

The socket of the HSE 600 has been mounted to a plate, and is equipped with a special feature that allows you to slightly adjust the socket afterwards. This may be necessary to realize the perfect lamp position once the lamp has been placed.

In HSE 1000 systems, the reflector should be clicked into the frame. The footprint of the reflector fits underneath the snappers in the lamp bracket. This should be executed with some caution to prevent distortion of the reflector. (**figure 3**)

### **6.2 Operating the Lamp without the Reflector**

Lamps may not be operated without reflectors. The amount of heat generated is hazardous to any materials present above the fixture.

### **6.3 Maintenance of the Reflector**

Alcohol based cleaning product can be used for regular cleaning. Avoid the use of abrasive cleaners on the reflectors and always ensure to rinse thoroughly with demineralized water. In general, calcium pollution can be removed with cleaning vinegar. If, for some reason, cleaning activities have not been executed for an extended period of time, the pollution may have corroded the anodizing layer. In such cases it is important to take appropriate action to clean the polluted surface. The use of more aggressive cleaning agents should be avoided to prevent irreparable damage to the reflectors. In case of grave pollution, please consult a specialist cleaner. In case of even more severe pollution, the reflectors may have to be subjected to a re-anodizing treatment.

Reflectors constructed from scaled materials (such as the Alpha reflector) cannot be re-anodized. Also, cleaning activities should be executed carefully to prevent loss of shape of the reflector.

## **7: Liability**

In the interest of promotion of fire safety in greenhouses, a critical assessment should be made of the spacing of the fixtures and any screening fabrics present, possibly in consultation with your insurance adviser.

Hortilux Schröder b.v. is not liable for any damages caused by failure to comply with the installation instructions in accordance with NEN1 010 and EnergieNed. publications, or any provisions assigned to function as substitutes thereof, and the above assembly instructions.

The conditions of the Metaalunie (*Dutch Metal Association*) and Hortilux Schröder b.v. additional terms of guarantee are exclusively applicable to all deliveries.

**Technical specs.****Input (mains side)**

Parameter	Conditions	600W			1000W			unit
		Min	typ	Max	min	typ	max	
Mains voltage $V_{\text{mains}}$	Operational performance	368	400	424	368	400	424	V
	Operational safety <sup>1</sup>	360		440	360		440	V
Mains frequency $f_{\text{mains}}$	Operational performance	47.5		63	47.5		63	Hz
	Operational safety	45		66	45		66	Hz
Mains power $P_{\text{mains}}$	$P_{\text{la}} = 615\text{W}$		635	640				W
	$P_{\text{la}} = 1000\text{W}$					1032	1040	W
Mains current $I_{\text{mains}}$	At $P_{\text{la nom}}$ , $V_{\text{mains}} = 400\text{V}$		1.61	1.65		2.61	2.68	A
	At $P_{\text{la max}}$ , $V_{\text{mains}} = 368\text{V}$			1.9			3.0	A
Power factor	Within oper. performance mains and at $P_{\text{la nom}}$	0.95	0.98		0.95	0.98		
Distortion	Within oper. performance mains and at $P_{\text{la nom}}$	According to EN61000-3-2			According to EN61000-3-2			
THD			9	15		9	15	%
Inrush current $I_{\text{mains\_pk}}$	$V_{\text{mains}} = 424\text{V}$ , $Z_{\text{mains}} = 0.4\Omega + 0.8\text{mH}$		120			175		A



<sup>1</sup> The input of the driver is protected against transients and over voltage typically happens at 460V, Beyond that the driver will fail (not repairable by the end user).

**Output (lamp side)**

Parameter	Conditions	600W			1000W			unit
		Min	typ	Max	min	typ	Max	
Lamp power $P_{\text{la}}$	$V_{\text{mains}} \pm 3\%$ ( $P_{\text{la nom}}$ )	597	615	633	970	1000	1030	W
	$V_{\text{mains}} -8\%, +6\%$ ( $P_{\text{la nom}}$ )	585		645	950		1050	W
Ignition voltage $V_{\text{ign}}$	$C_{\text{load}} < 100\text{pF}$ <sup>2</sup>			2500			2500	$V_{\text{pk}}$



<sup>2</sup> Typically cable capacity 100pF/meter

**Dimensions**

- \* Center of gravity reflector
- \*\* Stop swivel M25
- \*\*\*Center of gravity fixture

**HSE1000™**

**HSE600™**


figure 1


figure 2


figure 3

# Growing your profits

Horticultural lighting

**PHILIPS**

## Horticultural lighting and Philips

Philips has been developing light sources for horticulture for many years already. For a deeper understanding of what is required in the greenhouse Philips has close contact with commercial growers. In order to meet these specific requirements we have our own laboratories and test stations and, to further advance our overall knowledge, we contribute to independent research and field testing. This approach has led to the development of no-fuss, highly efficient lamps which have been tailor-made for you.

## The role of light

The amount of natural light (global radiation) is in most cases measured in terms of energy (J or W) with a solar meter. Plants use a relatively small part of this radiation for growth and this we call growth light. The majority of the radiation is heat. When you use supplemental light to enhance plant growth, you need to ensure that the lamps are highly efficient at producing growth light, and are not, for example, mostly producing heat.

As evidence mounts that artificial light can increase productivity, more and more growers are turning to artificial light. It is already very commonly used in ornamental crops such as roses, chrysanthemums and lilies and is now increasingly used for vegetables such as tomatoes, cucumber, sweet pepper and lettuce. It is expected that in the future its use will be extended to cover a much wider range of crops.

## The essence of good lighting


Light is essential for plant growth. Natural sunlight is the cheapest source available but for horticulture it is not always available in sufficient quantities. Especially in regions between 40 and 80 degrees latitude the amount of daylight required for good plant growth is limited during the winter. Therefore, during this period, the use of artificial light has become very common in greenhouses to increase production and quality.

## Contents

The essence of good lighting	2-3
Artificial lighting	4-7
Translating knowledge into products	8-9
Product datasheets	10-19


## Growth light

Light is, for the human eye, the visible part of electromagnetic radiation. Most products for lighting are developed for human applications. For these purposes the intensity of visible light is expressed in lux. Lux is a photometric unit and is based on the average sensitivity of the human eye.

This sensitivity is maximal at green/yellow (555 nm) and is declining towards longer (red) and shorter wavelengths (blue).

Plants have a completely different sensitivity for light colours than the human eye. For plant growth it is important to define light as small light particles, also called photons or quantum. The energy content of photons is different, depending on wavelength (light colour). For one Watt of energy, almost twice as many red photons can be produced compared with blue. This means that although they still use the green and blue part for growth – or photosynthesis - they use the red part of the light much more efficiently. In fact we are dealing with a plant sensitivity curve for growth light.


So, contrary to common belief, plant growth is not determined by lux or energy, but by photons from the blue to red (400-700 nm) part of the spectrum. This is called growth light!

## Suitability for photosynthesis

Research at universities and applied research stations has demonstrated that the rate of photosynthesis is related to the amount of photons between 400 – 700 nm. This is called 'Photosynthetic Photon Flux' (PPF). It is the only reliable way of measuring if a light source is suitable for photosynthesis.

The higher the PPF value per Watt, the more efficient the light source for plant growth. This is why Philips specifies on all his light sources for horticultural use, the PPF value. This is expressed in micromole photons per second ( $\mu\text{mol/s}$ ).

The Philips MASTER GreenPower lamp is specially developed for maximal growth light and has the highest PPF per Watt available for horticulture.

## Artificial lighting

There are several ways in which artificial light can be used to improve growth and extend the growing season of commercial crops:

- 1 to supplement natural daylight and raise growth light levels in order to enhance photosynthesis and thereby improve growth and quality of plants in greenhouses (**supplemental growth light**).
- 2 to control the light period by extending the natural day length with artificial light (**photoperiodic lighting**).
- 3 to totally replace daylight with artificial light for ultimate climate control (**cultivation without daylight**).

Philips offers a wide range of lamps for all these horticultural applications.

### I Supplemental growth light in greenhouses

The amount of supplemental growth light required very much depends on plant type, desired plant growth and availability of natural daylight. For this reason Philips has designed a tool that calculates how much additional growth light is required in each individual situation. The tool is available at your local Philips office and allows you to calculate the optimal amount of light your greenhouse requires based on your wishes and/or possibilities.

Depending on plant type and desired plant growth for central European conditions, the following supplemental levels are suggested:

- 1 15 – 30  $\mu\text{mol}$  for improving quality, maintenance of the crop and limited production increase;
- 2 30 – 45  $\mu\text{mol}$  for seedlings, growth and production of pot plants;
- 3 40 – 100  $\mu\text{mol}$  for year-round cultivation, for example, of chrysanthemums and roses and multiple layer cultivation;
- 4 100 – 200  $\mu\text{mol}$  for production of plants with high light demand (fruit production of, for example, tomatoes and cucumbers);
- 5 100 – 800  $\mu\text{mol}$  for the production of plants under artificial light alone (for example growth chambers)

In the case of MASTER GreenPower 600W/400V: 1  $\mu\text{mol}$  growth light corresponds to 76 lux.

For supplemental growth light applications Philips recommends:

- **MASTER GreenPower**

Figure I Supplemental growth light


Figure I shows an example of how natural daylight is supplemented with MASTER GreenPower during winter. In this example plants are illuminated with 105  $\mu\text{mol}$  growth light (= circa 8000 lux) during 20 hours/day from November until February. In the remaining lighting period the operating hours is less.

**For the low dose of light applications  
Philips recommends:**

- Incandescent (Flower Power Pro or Superlux Agro Pro)
- Compact fluorescent (CFL)

## 2 Photoperiodic lighting

For many plants the moment of flowering is determined by the length of the light period. The use of artificial light for control of flowering is called photoperiodic lighting. With this method short-day and long-day plants can be cultivated all year round. For example, very good results have been achieved with photoperiodic lighting of chrysanthemum, euphorbia pulcherrima and kalanchoe as well as with gypsophila and carnations.

The most common ways of influencing the day length are:

- with growth light (also growth takes place during the day prolongation period) or
- with a low dose of light, primary for flower regulation (100 - 400 lux; 2 - 6  $\mu\text{mol}/\text{m}^2\cdot\text{s}$ ).

## 3 Cultivation without daylight

The total replacement of daylight by artificial light is primarily associated with climate controlled rooms. In these applications, it is essential that the spectral composition of the artificial light is balanced for optimal plant development.

For cultivation of plants without daylight Philips offers you three solutions:

- a. MASTER HPI-T Plus  
HPI-T Plus lamps have a spectrum closer to daylight, and plants develop very well under this light.
- b. A 1:1 mix of 400W MASTER GreenPower and 400W MASTER HPI-T Plus  
This combination of light sources has proven a good lighting solution for several species. Compared with HPI-T Plus lamps, MASTER GreenPower lamps offer a higher efficiency in growth light. The little blue in the spectrum of the MASTER GreenPower is compensated by HPI-T Plus.
- c. MASTER TL-D Super 80  
The use of MASTER TL-D Super 80 fluorescent lamps (for example colour 830 and 840) is also a very suitable option.  
TL-D Super 80 lamps offer the following advantages:
  - Low temperature of the tube. This allows short distance between lamps and plants (around 15 cm) and thus the possibility to grow plants in multiple layers.
  - TL-D Reflex is particularly suitable for this application. This tube has an internal reflector that creates an efficient lighting system without the need for an external reflector
  - Flexible light levels can be arranged. The light level can vary from very low (for tissue culture and seedlings) to high (approximately 800  $\mu\text{mol}$ ) with good light distribution.
  - Dimming is possible. The light output of a high frequency system (TL-D HF) is continuously dimmable between 100% down to approximately 8%. The efficiency of the fluorescent lamp is comparable with that of the HPI-T Plus.

MASTER TL-D for cultivation without daylight

**For cultivation without daylight  
Philips recommends:**

- MASTER HPI-T Plus
- Mix of MASTER GreenPower and MASTER HPI-T Plus
- MASTER TL-D Super 80

1:1 mix of MASTER GreenPower and MASTER HPI-T Plus

# Translating knowledge into products

Research to improve knowledge is essential. Philips supports several research projects that study the process of plant growth in relation to light. Leading universities in Europe and the USA, together with several growers and luminaire suppliers, have helped us to determine the most efficient lighting solution for ornamental crops and vegetables. Our international contacts ensure that our level of expertise is as broad as possible. A number of commercial growers with different crops allow us to test new concepts over lifetime. This approach means that new or adapted light sources are intensively evaluated by international experts as well as being tested in practice.

This approach has resulted in Philips designing horticultural lamps that offer you:

- Maximal efficiency in growth light
- Maximal output over life time
- Maximal benefit on plant growth
- Minimal early failure.

## Light measurements

Philips has its own independently and officially certified light measurement laboratory. This means measurements can be carried out on the performance of our lamps with the highest accuracy. Here we do all measurements on light sources specified in the IEC standards.

Light measurement laboratory Philips Turnhout.

## Stray light

With the increasing use of artificial light in greenhouses, stray light that reflects back out of the greenhouse, is becoming an environmental and social issue. In order to operate in a socially responsible manner, we recommend a range of technical measures, such as the installation of screens, to reduce this issue. Philips is also taking the issue of stray light into account as it develops new technologies for greenhouse lighting systems in the future.

85% screen (partially closed) installed to reduce light pollution

## Our customers

It is very important for Philips that you are satisfied with our products. We are developing lamps on the input we get from the market, a policy that has brought us very positive feedback from growers.


*Peter Klapwijk, Klapwijk GreenQ bv Tomato Grower - Monster, The Netherlands*  
"A good lamp is determined by yield times lifetime. That automatically leads you to Philips. Our contact is good and that's important, as technical development and optimal plant breeding should go hand in hand."

*T. Verheul, Rosa Plaza de Berckt Roses - Baarlo, The Netherlands*  
"I recently visited the Philips Turnhout factory where I received useful and clear information about the lamps and the use of these lamps in horticultural applications."

*Cor Boeters, Sunrise Pepper Growers - Wieringen, The Netherlands*  
"Reliability of the lamp and the manufacturer is our reason for choosing Philips. Our contact is good and people think along with test projects."

## Philips support

Your local Philips agent can demonstrate the benefits of using Philips MASTER GreenPower and other lamps for several crops.

## Website Horticultural Lighting


Detailed information and recommendations about using additional growth light in horticulture can be found on our horticultural website:

[www.horticultural.lighting.philips.com](http://www.horticultural.lighting.philips.com)


Product datasheets

Lifetime testing at Philips Turnhout

## MASTER GreenPower and MASTER Agro Leading technology

### ① PIA technology

The ceramic discharge tube with "Philips Integrated Antenna" technology contains no moving parts and therefore eliminates early lamp failures and guarantees a reliable and extended lamp life.

### ② ZrAl getter

The "Zirconium Aluminium" getter improves the vacuum in the outer bulb. This ensures a more consistent discharge temperature and lower sodium migration from the discharge tube resulting in excellent lumen maintenance over the total lifetime as well as fewer premature failures.

### ③ Simple and strong construction

The optimised design with only seven internal welds ensures a more robust construction to withstand vibration and adverse environmental conditions, enhancing reliability and prolonging the lifetime.

### ④ Lead free

The range is constructed with lead-free solder and is 96% "upward" recyclable (the whole lamp with the exception of the discharge tube) in an environmentally conscious way. This means that the recycled components are re-used in the production of new lamps.

MASTER GreenPower T  
400W / 600W EM for  
230V single phase  
systems

#### Lamp


- MASTER GreenPower lamps are High Pressure Sodium Lamps with a ceramic discharge tube, enclosed in a clear tubular outer bulb with optimised growth light ( $\mu\text{mol}$ ) output and maintenance.

#### Features

- GreenPower results in optimised growth light output.
- Ceramic discharge tube with PIA technology for long and reliable lifetime.
- ZrAl getter ensures excellent growth light maintenance over life and fewer premature failures.
- Simple and robust construction for enhanced reliability and longer life.
- Lead free solder.

#### Benefits

- High growth light maintenance safeguards a constant crop quality and quantity over life.
- Optimised performance for horticultural applications.

#### Comparison of MASTER GreenPower with MASTER SON-T PIA Plus

- MASTER GreenPower lamps are designed for optimal growth light output over lifetime based on an optimal spectral energy distribution for the light sensitivity curve of plants.
- MASTER SON-T PIA Plus lamps are designed for maximal lumen output over lifetime based on an optimal spectral energy distribution for the light sensitivity curve of the human eye.

#### Application

- Horticultural lighting, intended for the stimulation of  $\text{CO}_2$  uptake for improved photosynthesis and plant growth.

#### Gear

- The lamp requires a ballast and ignitor in accordance with the IEC HPS Plus standard.

#### Product information

	Growth light PPF* (initial)	Lifetime	Maintenance (nominal)
MASTER GreenPower 400W	725 $\mu\text{mol}/\text{sec}$	12.000 hrs	> 92% at 12.000 hrs
MASTER GreenPower 600W	1100 $\mu\text{mol}/\text{sec}$	10.000 hrs	> 92% at 10.000 hrs

\* PPF information: see page 3

Perfect light, perfect sense

MASTER GreenPower T  
600W EM for 400V  
phase/phase systems

#### Lamp

- MASTER GreenPower lamps are High Pressure Sodium Lamps with a ceramic discharge tube, enclosed in a clear tubular outer bulb with optimised growth light ( $\mu\text{mol}$ ) output and maintenance.

#### System

The system consists of a specially designed lamp, ballast and ignitor, suitable for phase/phase connection to the mains.

#### Features

- GreenPower results in optimised growth light output.
- Ceramic discharge tube with PIA technology for long and reliable lifetime.
- ZrAl getter ensures excellent growth light maintenance over life and fewer premature failures.
- Simple and robust construction for enhanced reliability and longer life.
- Lead free solder.

#### Benefits

- High growth light maintenance safeguards a constant crop quality and quantity over life.
- Optimised performance for horticultural applications.

#### Comparison of 400V system with 230V system

- Lower installation cost.
- Very low 3rd harmonic.
- Improved growth light, +4,5%.

#### Application

- Horticultural lighting, intended for the stimulation of  $\text{CO}_2$  uptake for improved photosynthesis and plant growth.

#### Gear

- The lamp requires a ballast and ignitor in accordance with the IEC HPS Plus standard.


#### Recommendation lamp replacement

- Luminaire must be disconnected from the mains before lamp replacement because outer ring of the lamp holder is connected to a phase and therefore live.

#### Product information

	Growth light PPF* (initial)	Lifetime	Maintenance (nominal)
MASTER GreenPower 600W/400V	1150 $\mu\text{mol}/\text{sec}$	10.000 hrs	> 92% at 10.000 hrs

\* PPF information: see page 3

Perfect light, perfect sense

MASTER Agro T  
400W for 230V single  
phase systems

Lamp

MASTER Agro lamps are High Pressure Sodium Lamps with a ceramic discharge tube, enclosed in a clear tubular outer bulb with an increased output of blue light designed for horticultural purposes.

Features

- Agro results in extra blue in the spectral energy distribution.
- Ceramic discharge tube with PIA technology for long and reliable lifetime.
- ZrAl getter ensures excellent growth light maintenance over life and fewer premature failures.
- Simple and robust construction enhancing reliability and longer life.
- Lead free solder.

Benefits

- High growth light maintenance safeguards a constant crop quality and quantity over life.
- Especially at lower lighting levels the spectrum of the MASTER Agro lamp will result in a more compact plant development for certain plants and can optimise plant development and quality with good leaf, bloom color and extra branching.

Comparison of MASTER Agro with MASTER GreenPower

- MASTER Agro lamps are designed for applications with low light levels. The extra blue in the spectral energy distribution results in more compact and sturdy plants.
- MASTER GreenPower lamps are designed for optimal growth light output over lifetime based on an optimal spectral energy distribution for the light sensitivity curve of plants.

Application

- Horticultural lighting, intended for the stimulation of CO<sub>2</sub> uptake for improved photosynthesis and plant growth.

Gear

- The lamp requires a ballast and ignitor in accordance with the IEC HPS Plus standard.

Product information

	Growth light PPF* (initial)	Lifetime	Maintenance (nominal)
MASTER SON-T PIA Agro 400W	660 µmol/sec	10.000 hrs	> 85% at 10.000 hrs

\* PPF information: see page 3

Perfect light, perfect sense

MASTER HPI-T Plus  
400W for 230V single  
phase systems

Lamp

MASTER HPI-T Plus lamps are Metal Halide Lamps with a quartz discharge tube, enclosed in a clear tubular outer bulb.

Features

- 3 band technology resulting in high luminous efficacy, both initially and over long lifetime.
- 3 band technology resulting in high color stability for stable plant growth over lifetime.
- Runs both on HPI gear and for higher growth light output on SON gear.

Benefits

- Minimal maintenance cost.
- Initial investment saving options when run on SON ballast resulting in higher growth light output.

Applications

- Horticultural lighting, intended for the stimulation of CO<sub>2</sub> uptake for improved photosynthesis and plant growth.
- Horticultural applications with low levels of daylight in combination with SON lamps.
- Horticultural applications without daylight like growing chambers.

Gear

- Can be run on HPI gear as well as SON gear, resulting in different light output and color temperature.


Recommendation lamp usage

- Luminaire with protective front glass is required.

Product information

	Growth light PPF* (initial)	Lifetime	Maintenance (nominal)
MASTER HPI-T Plus (on SON gear)	540 µmol/sec	7.000 hrs	> 85% at 7.000 hrs
MASTER HPI-T Plus (on HPI gear)	490 µmol/sec	7.000 hrs	> 85% at 7.000 hrs

\* PPF information: see page 3

MASTER HPI-T Plus 400W 230V

Perfect light, perfect sense

# MASTER TL-D Reflex Super 80

## Lamp

MASTER TL-D Reflex lamps are fluorescent lamps with an internal reflector to concentrate the lamplight in the direction in which it is needed.

## Features

- Internal reflector with an opening angle of 160° increases the light intensity by 60%.
- High growth light maintenance.
- MASTER TL-D Reflex lamps are 100% retrofit with all TL-D lamps with similar wattage.
- Recyclable; mercury, phosphor and glass can be re-used in production of new TL-D lamps.

## Benefits

- Lower initial investment because of built in reflector.
- Internal reflector reduces used space in multi layer cultivation.
- Light output is hardly affected in environments subject to dust accumulation.

## Comparison of MASTER TL-D Reflex with standard TL-D

- Up to 60% higher light output resulting in higher yield for existing installation or less battens for new installations.
- 60% longer service lifetime reduces the maintenance and relamping cost.

## Applications

- Cultivation of plants without daylight.
- Growth chamber or cabinet with possibility to grow plants in multiple layers.

## Gear

- Operates both on conventional but preferably on HF control gear.

## Product information

	Growth light PPF* (initial)	Lifetime	Maintenance (nominal)
MASTER TL-D Reflex Super 80 36W	47 µmol/sec	12.000 hrs	> 90% at 12.000 hrs
MASTER TL-D Reflex Super 80 58W	73 µmol/sec	12.000 hrs	> 90% at 12.000 hrs

\* PPF information: see page 3

Perfect light, perfect sense

# MASTER TL-D Secura Super 80

## Lamp

- MASTER TL-D Secura lamps are fluorescent lamps with an external transparent protective coating.

## Features

- External transparent protective coating.
- Easily identified by a blue ring at the end of the lamp.
- High growth light maintenance.
- MASTER TL-D Secura lamps are 100% retrofit with all TL-D lamps with similar wattage.
- Recyclable; mercury, phosphor and glass can be re-used in production of new TL-D lamps.

## Benefits

- Protective coating safely retains all glass and other lamp components if the lamp breaks.

## Comparison of MASTER TL-D Secura with standard TL-D

- Up to 30% higher light output resulting in higher yield for existing installation or less battens for new installations.
- 60% longer service lifetime reduces the maintenance and relamping cost.
- Protective coating ensures a safe environment without the possibility of glass fall.

## Applications

- Cultivation of plants without daylight.
- Growth chambers or cabinets where safety is essential.



## Gear

- Operates both on conventional but preferably on HF control gear.

## Product information

	Growth light PPF* (initial)	Lifetime	Maintenance (nominal)
MASTER TL-D Secura Super 80 36W	45 µmol/sec	12.000 hrs	> 90% at 12.000 hrs
MASTER TL-D Secura Super 80 58W	70 µmol/sec	12.000 hrs	> 90% at 12.000 hrs

\* PPF information: see page 3

Secura

Standard TL-D

Perfect light, perfect sense

## Flower Power Pro

100W

### Lamp

- Flower Power Pro is an incandescent lamp with internal mirror reflector for horticultural applications.

### Features

- Built-in reflector eliminates the need for an external reflector system.
- Corrosion free nickel-plated brass cap.

### Benefits

- Can directly replace existing 100W and 150W lamps.
- Internal reflector ensures that all light is directed to where it is needed.

### Comparison of Flower Power Pro with Superlux Agro Pro

- For less energy consumption;  
Flower Power Pro 100W with a lifetime of 1.000 hrs.
- For longer lifetime;  
Superlux Agro Pro 150W with a lifetime of 2.000 hrs.

### Applications

- Photoperiodic lighting in greenhouses.
- To prolong the short days during autumn, winter and spring.

### Application advice

- In an installation of 3 x 3.2 m the Flower Power Pro has an installed output of 10.4 w/m<sup>2</sup>.
- For stock plants; one lamp per 7.5 m<sup>2</sup>.
- For cut-flowers; one lamp per 9.6 m<sup>2</sup>.
- Minimal distance of 80 cm between possible obstacle and plants should be maintained due to shadow casting.
- Minimal distance between lamp and plant (Hnett) is 2 m, perpendicular to the plants.
- Nominal voltage variations can influence the lifetime of the lamps.

### Product information

	Light intensity at Hnett = 2 m and 1 lamp per 9,6 m <sup>2</sup>	Lifetime (average)
Flower Power Pro 100W	100 lux (min/max > 70%)	1.000 hrs

## Superlux Agro Pro

150W

### Lamp

- Superlux Agro Pro is an incandescent lamp with a mushroom shaped white coated bulb for horticultural applications.

### Features

- Internal white coating eliminates the need for an external reflector system.
- Corrosion free nickel-plated brass cap.
- Double lifetime compared with standard incandescent lamps.

### Benefits

- Can directly replace existing 150W lamps clear or diffuse.
- Internal white coating ensures that all light is directed to where it is needed.

### Comparison of Superlux Agro Pro with Flower Power Pro

- For longer lifetime;  
Superlux Agro Pro 150W with a lifetime of 2.000 hrs.
- For less energy consumption;  
Flower Power Pro 100W with a lifetime of 1.000 hrs.

### Applications



- Photoperiodic lighting in greenhouses.
- To prolong the short days during autumn, winter and spring.

### Application advice

- In an installation of 3 x 3.2 m the Superlux Agro Pro has an installed output of 15.6 w/m<sup>2</sup>.
- For stock plants; one lamp per 7.5 m<sup>2</sup>.
- For cut-flowers; one lamp per 9.6 m<sup>2</sup>.
- Minimal distance of 80 cm between possible obstacle and plants should be maintained due to shadow casting.
- Minimal distance between lamp and plant (Hnett) is 2 m, perpendicular to the plants.
- Nominal voltage variations can influence the lifetime of the lamps.

### Product information

	Light intensity at Hnett = 2 m and 1 lamp per 9,6 m <sup>2</sup>	Lifetime (average)
Superlux Agro Pro 150W	100 lux (min/max > 70%)	2.000 hrs

Philips Lighting

Lamps Europe - Professional Lamps

[www.horticultural.lighting.philips.com](http://www.horticultural.lighting.philips.com)

Printed in the Netherlands

Although this information has been prepared with extreme care, no rights can be derived from any information and illustrations in this brochure.

## **BREAKING NEWS ON SUPPLEMENTAL LIGHTING!**

Hortilux Schröder introduced the following step in supplemental lighting through their **HSE600** and **HSE1000**. The basic principles for developing the new fixtures, based on an electronic ballast device, was of course that it would work immediately and would be great and reliable. That is why Hortilux Schröder can guarantee more light for less energy and a new and more energy-efficient concept with a longer life span. The lamps shine better and are more efficiently and stable, which helps prolong the life span. An additional advantage is that the fixtures themselves have become lighter in weight and more compact, hence less shadows. The reflector concept has also been upgraded and improved, which helps to improve the light yield even more and creates an optimum uniformity.

### **HSE also in 1000 Watt!**

The electronic ballast device from Hortilux Schröder is available in two types: 400 Volt/600 Watt and 400 Volt/1000 Watt. The **HSE1000** is a real breakthrough – an innovation that allows you to use fewer fixtures. Fewer fixtures mean less shadow, less connection points and lower installation and maintenance costs. In this way, the introduction of Electronics offers a better use of capital as regards investments and the possibility to work with higher light levels.

## **New electronic supplemental lighting system:**

### **less energy, more light!**

Nom. ambient temperature                      -10 °C + 40 °C    -10 °C + 40 °C

## HS Electronics Technical Specifications

	<b>HS E600</b>	<b>HS E1000</b>
Type of lamp: high pressure sodium lamp	Master GreenPower-EI 600 Watt	Master GreenPower-EI 1000 Watt
Type of ballast device	GreenVision 600 Watt	GreenVision 1000 Watt
Type of fitting	E-40	Dubbel ended
Voltage (nom)	400 Volt	400 Volt
Voltage tolerance	-8/+6%	-8/+6%
Frequency	50 Hz	50 Hz
Electricity consumption (nom)	1,66 Amp	2,70 Amp
Absorbed voltage (nom)	645 Watts	1045 Watts
Power factor	>0.97	>0.97
Lamp's voltage (nom)	615 Watts	1000 Watts
Lamp's voltage tolerance	± 3%	± 3%
Weight	± 3,5 Kg	± 3,75 Kg

## HORTILUX SCHRÉDER'S DELIVERY PROGRAMME

Besides the **HS E 600** and the **HS E 1000**, of course the familiar Hortilux Schröder systems remain part of our programme. An overview:

System	Fixtures	Reflectors
Standard	HS 2000/400 Watt HS 2000/600 Watt HS Remote	Special Deep Deep Midi Medium Wide Super Wide
Electronics	<b>HS E 600</b> <b>HS E1000</b>	Alpha

## Light calculation

Customer	University of Guelph
City / Country	CAN
Pattern	1 m x 2 m H = 1 m + 1 m x 2 m H = 1 m
Calculator	P. van der Valk
Date	24-6-2005
Drawing number	NVT University of Guelph

### Luminaire information:

no	type	luminous flux [lumen]	x [m]	y [m]	z [m]	rx	ry	rz
1	981412 DEEP HS 600W grid: 1 x 3 armatures configuration: 1 m x 2 m H = 1 m	600W SONT 85000	0	0	1	0	0	90
2	921821 DEEP HS 400W grid: 1 x 2 armatures configuration: 1 m x 2 m H = 1 m	400W K P HPI/T M-H 37000	0	0	1	0	0	90


Maintenance factor = 1.00

### Luminaire information:

average (Eav):	23261 lx
minimum (Emin):	12932 lx
maximum (Emax):	29082 lx
Emin/Emax:	44,5 %
Emin/Eav:	55,6 %

### Uniformity in calculation field:

Emin 12932 lx Emax 29082 lx Eaverage 23261 lx  
Emin/Emax 44.5% Emin/Eaverage 55.6%


To this lighting plan (drawing + computer calculation the "Explanation to the Lightingplan",  
the supplementary guarantee provisions PL light systems Inc. apply.

**Calculation plane:**

14031	16663	17263	15442	12932
17306	21727	22989	21395	16774
20964	24968	26127	23813	19963
22202	25457	25423	24453	21345
23019	27366	28771	26212	22023
22735	27696	29082	27359	22203
23019	27366	28771	26212	22023
22202	25457	25423	24453	21345
20964	24968	26127	23813	19963
17306	21727	22989	21395	16774
14031	16663	17263	15442	12932

**Calculation field:**

**To this lighting plan (drawing + computer calculation the "Explanation to the Lightingplan", the supplementary guarantee provisions PL light systems Inc. apply.**

# P.L. LIGHT SYSTEMS

TECHNICAL SPECIFICATIONS  
AND REFERENCE GUIDE.

2006

**P.L. LIGHT SYSTEMS**  
**TECHNICAL SPECIFICATIONS**  
**AND REFERENCE GUIDE.**  
**DIRECTORY**

**1. General Information**

- [Definitions](#)
- [Fixture / Reflector Type, Dimensions, Weight](#)
- [Shipping Weights and Dimensions](#)
- [Greenhouse Shade Curtains](#)
- [Reflector Cleaning Instructions](#)

**2. Ballast Information**

- [Changing the Ballast Voltage](#)
- [Ballast Circuitry](#)
- Advance - [HPS MH](#)
- Venture - [HPS MH](#)
- [Ignitor Operating Life](#)
- [Ignitor Measurements](#)
- [Ignitor / Capacitor Cross Reference Guide](#)
- [Capacitor Info](#)

**3. Bulb Information**

- [GE](#)
- [Philips](#)

**4. Reflector Information**

- [Reflector Overview](#)
- [Deep](#)
- [Maxima](#)
- [Midi](#)
- [Superwide](#)
- [Wide](#)
- [Alpha \(HSE\)](#)

**5. Strut Information**

- [12awg, 14awg](#)
- [Outlet Boxes, End Cap, Closure Strip](#)

**6. Receptacles**

- [Straight Blade](#)
- [Twist Lock](#)

**7. [Request for Light Plan](#)**

**8. [Trouble Shooting Guide](#)**

# GENERAL INFORMATION

## DEFINITIONS

*VOLT* - The (electrical) potential difference between two points in a circuit.

*AMP* - The unit of measurement of electric current. It is proportional to the quantity of electrons flowing through a conductor past a given point in one second

*VA* - Power attained in an AC circuit as a product of effective voltage and current which reach their peak at different times.

*Watts* - The standard unit of measurement of electrical power. One watt is one ampere of current flowing at one volt. Watts are typically rated as AMPS x VOLTS or VOLT-AMP (V-A). However, this rating is only equivalent to watts when it applies to devices that absorb all the energy, such as electric heating coils or incandescent light bulbs. With inductive reactive loads, the actual watt rating is different than that of the VOLT-AMP rating.

*Power Factor* - In an alternating current system the voltage and current do not always reverse at the same instant in time. That is, they are not always "in phase." The current can be considered as being divided into two components: one in phase with the voltage and one out of phase with the voltage. The power factor of a circuit is the ratio of the in phase current to the total current. Usually expressed as a percentage. Power is the product of volts, ampere, and power factor.

*Total Harmonic Distortion* - The percentage of distortion to the fundamental frequency caused by harmonics within the electrical system.

*Line Loss* - Voltage drop over the length of electric line wire. Line loss robs your system of power when wire is too small for the load being run through the line or when voltage is too low for the distance the power must travel.

*Capacitor* - Capacitance refers to the ability of a device to store energy in the form of an electrostatic field. In its simplest form, a capacitor is a pair of parallel plates spaced apart with a dielectric material between them.

*Microfarad* - The term microfarad is used to note the storage capacity value of electrolytic capacitors and represents one millionth of a farad, or 0.000 001 farads. It is most often abbreviated  $\mu\text{F}$ , uF, UFD, MF, and MFD, although the latter two can be misleading, since it can be mistaken for mill farad. Since mill farad is almost never used in electronics engineering, this discrepancy is usually not a problem.

*Ballast* - A device used with an electric-discharge lamp to obtain the necessary circuit conditions (voltage, current and wave form) for starting and operating.


*Color Codes* - Conductors of a three phase system are usually identified by a color code, to allow for balanced loading. Colors used may adhere to old standards or to no standard at all, and may vary even within a single installation.

	L1	L2	L3	Neutral	Earth
North America Canada 600/ 347V	Black	Red	Blue	White	Green
North America USA (newer 277/480 installations)	Brown	Orange	Yellow	White	Green
UK until April 2006	Red	Yellow	Blue	Black	green/yellow striped (green on very old installs)
Europe (including UK) from April 2004	Brown	Black	Grey	Blue	green/yellow striped
Previous European (varies by country)	Brown or black	Brown or black	Brown or black	Blue	green/yellow striped
Europe, for bus bars	Yellow	Green	Purple	–	–

\* Note that in N.A. a green/yellow striped wire typically indicates an isolated ground.

*Circuit Numbering* - The circuit numbers of a three phase system are identified as follows. The following chart is true for three phase power (Black, Red, Blue) and 277/480 (Brown, Orange, Yellow) installations. The purpose of the standardized circuit number to wire color is for future troubleshooting. If an electrician knows that the wire running to a nonworking plug is blue, the electrician can instantly eliminate 2/3 of the remaining circuits.

Black/Brown	Red/Orange	Blue/Yellow
1/2	3/4	5/6
7/8	9/10	11/12
13/14	15/16	17/18
19/20	21/22	23/24
25/26	27/28	29/30
31/32	33/34	35/36
37/38	39/40	41/42

*Single phase loads* - Single-phase loads may be connected to a three-phase system, either by connecting across two live conductors (a phase-to-phase connection), or by connecting between a phase conductor and the system neutral which must be connected to the center of the Y(star)secondary winding of the supply transformer. Single-phase loads should be distributed evenly between the phases of the three-phase system for efficient use of the supply transformer and supply conductors.

The line-to-line voltage of a three-phase system is  $\sqrt{3}$  (1.732) times the line to neutral voltage. Where the line-to-neutral voltage is a standard utilization voltage, (for example in a 277 V/480 V system) individual single-phase utility customers or loads may each be connected to a different phase of the supply. Where the line-to-neutral voltage is not a common utilization voltage single-phase loads must be supplied by individual step-down transformers. In multiple-unit residential buildings in North America, lighting and convenience outlets can be connected line-to-neutral to give the 120 V utilization voltage, and high-power loads such as cooking equipment, space heating, water heaters, or air conditioning can be connected across two phases to give 208 V. This practice is common enough that 208 V single-phase equipment is readily available in North America. Attempts to use the more common 120/240 V equipment intended for three wire single phase distribution may result in poor performance since 240 V equipment will only produce 75% of its rating when operated at 208V.

*Connecting phase-phase* - Connecting between two phases provides  $\sqrt{3}$  or 173% of the single-phase voltage (208 VAC in US; 400 VAC in Europe) because the out-of-phase waveforms add to provide a higher peak voltage in the resulting waveform. Such connection is referred to as a line to line connection and is usually done with a two pole circuit breaker. This kind of connection is typically used for high power appliances, such as a 2kW, 208 volt baseboard heater.

*Three phase loads* - The most important class of three-phase load is the electric motor. A three phase induction motor has a simple design, inherently high starting torque, and high efficiency. Such motors are applied in industry for pumps, fans, blowers, compressors, conveyor drives, and many other kinds of motor-driven equipment. A three-phase motor will be more compact and less costly than a single-phase motor of the same voltage class and rating; and single-phase AC motors above 10 HP (7.5

Resistance heating loads such as electric boilers space heating may be connected to three-phase systems. Electric lighting may also be similarly connected. These types of loads do not require the revolving magnetic field characteristic of three-phase motors but take advantage of the higher voltage and power level usually associated with three-phase distribution.

**To be used as reference only.**

Please contact a Certified Electrician, Engineer or the Local Authority having Jurisdiction for requirements of installation practices in you area.

## FIXTURE, REFLECTOR, DIMENSIONS, WEIGHT

Fixture Type: <b>150W HPS</b>	Reflector Type	Length	Maximum Width	Depth:	Weight
	<i>Midi</i>	26.9" (67.25 cm)	11.4" (28.5 cm)	5.5" (13.7 cm)	14 lbs (6.3 kg)
	<i>Medium</i>	27.4" (68.5 cm)	13.4" (33.6 cm)	5.5" (13.7 cm)	14 lbs (6.3 kg)
	<i>Deep</i>	26.7" (66.75 cm)	9.2" (23.1 cm)	5.5" (13.7 cm)	14 lbs (6.3 kg)
	<i>Wide</i>	26.5" (66.25 cm)	16.9" (42.2 cm)	5.5" (13.7 cm)	14 lbs (6.3 kg)
	<i>Superwide</i>	25.8" (64.5 cm)	10.9" (27.2 cm)	5.5" (13.7 cm)	14 lbs (6.3 kg)

Fixture Type: <b>250W HPS</b>	Reflector Type	Length	Maximum Width	Depth:	Weight
	<i>Midi</i>	26.9" (67.25 cm)	11.4" (28.5 cm)	5.5" (13.7 cm)	18 lbs (8.2 kg)
	<i>Medium</i>	27.4" (68.5 cm)	13.4" (33.6 cm)	5.5" (13.7 cm)	18 lbs (8.2 kg)
	<i>Deep</i>	26.7" (66.75 cm)	9.2" (23.1 cm)	5.5" (13.7 cm)	18 lbs (8.2 kg)
	<i>Wide</i>	26.5" (66.25 cm)	16.9" (42.2 cm)	5.5" (13.7 cm)	18 lbs (8.2 kg)
	<i>Superwide</i>	25.8" (64.5 cm)	10.9" (27.2 cm)	5.5" (13.7 cm)	18 lbs (8.2 kg)

Fixture Type: <b>400W HPS</b>	Reflector Type	Length	Maximum Width	Depth:	Weight
	<i>Midi</i>	26.9" (67.25 cm)	11.4" (28.5 cm)	5.5" (13.7 cm)	23 lbs (10.5 kg)
	<i>Medium</i>	27.4" (68.5 cm)	13.4" (33.6 cm)	5.5" (13.7 cm)	23 lbs (10.5 kg)
	<i>Deep</i>	26.7" (66.75 cm)	9.2" (23.1 cm)	5.5" (13.7 cm)	23 lbs (10.5 kg)
	<i>Wide</i>	26.5" (66.25 cm)	16.9" (42.2 cm)	5.5" (13.7 cm)	23 lbs (10.5 kg)
	<i>Superwide</i>	25.8" (64.5 cm)	10.9" (27.2 cm)	5.5" (13.7 cm)	23 lbs (10.5 kg)

Fixture Type: <b>430W HPS</b>	Reflector Type	Length	Maximum Width	Depth:	Weight
	<i>Midi</i>	26.9" (67.25 cm)	11.4" (28.5 cm)	5.5" (13.7 cm)	26 lbs (11.8 kg)
	<i>Medium</i>	27.4" (68.5 cm)	13.4" (33.6 cm)	5.5" (13.7 cm)	26 lbs (11.8 kg)
	<i>Deep</i>	26.7" (66.75 cm)	9.2" (23.1 cm)	5.5" (13.7 cm)	26 lbs (11.8 kg)
	<i>Wide</i>	26.5" (66.25 cm)	16.9" (42.2 cm)	5.5" (13.7 cm)	26 lbs (11.8 kg)
	<i>Superwide</i>	25.8" (64.5 cm)	10.9" (27.2 cm)	5.5" (13.7 cm)	26 lbs (11.8 kg)

## FIXTURE, REFLECTOR, DIMENSIONS, WEIGHT

Fixture Type: <b>600W HPS</b>	Reflector Type	Length	Maximum Width	Depth:	Weight
	<i>Midi</i>	26.9" (67.25 cm)	11.4" (28.5 cm)	5.5" (13.7 cm)	34 lbs (15.5 kg)
	<i>Medium</i>	27.4" (68.5 cm)	13.4" (33.6 cm)	5.5" (13.7 cm)	34 lbs (15.5 kg)
	<i>Deep</i>	26.7" (66.75 cm)	9.2" (23.1 cm)	5.5" (13.7 cm)	34 lbs (15.5 kg)
	<i>Wide</i>	26.5" (66.25 cm)	16.9" (42.2 cm)	5.5" (13.7 cm)	34 lbs (15.5 kg)
	<i>Superwide</i>	25.8" (64.5 cm)	10.9" (27.2 cm)	5.5" (13.7 cm)	34 lbs (15.5 kg)

Fixture Type: <b>1000W HPS</b>	Reflector Type	Length	Maximum Width	Depth:	Weight
	<i>1000W</i>	31.6" (79.0 cm)	13.1" (32.7 cm)	5.5" (13.7 cm)	43 lbs (19.5 kg)
	<i>Maxima</i>	31.6" (79.0 cm)	13.1" (32.7 cm)	5.5" (13.7 cm)	43 lbs (19.5 kg)

Fixture Type: <b>250W MH</b>	Reflector Type	Length	Maximum Width	Depth:	Weight
	<i>Medium MH</i>	27.4" (68.5 cm)	13.4" (33.6 cm)	5.5" (13.7 cm)	16 lbs (7.3 kg)
	<i>Deep MH</i>	26.7" (66.75 cm)	9.2" (23.1 cm)	5.5" (13.7 cm)	16 lbs (7.3 kg)

Fixture Type: <b>400W MH</b>	Reflector Type	Length	Maximum Width	Depth:	Weight
	<i>Medium MH</i>	27.4" (68.5 cm)	13.4" (33.6 cm)	5.5" (13.7 cm)	20 lbs (9.0 kg)
	<i>Deep MH</i>	26.7" (66.75 cm)	9.2" (23.1 cm)	5.5" (13.7 cm)	20 lbs (9.0 kg)

Fixture Type: <b>1000W MH</b>	Reflector Type	Length	Maximum Width	Depth:	Weight
	<i>1000W MH</i>	31.6" (79.0 cm)	13.1" (32.7 cm)	5.5" (13.7 cm)	40 lbs (18.2 kg)

## SHIPPING WEIGHTS AND DIMENSIONS

COMMERCIAL FIXTURES			
TYPE	KG	LBS	BOX DIMS.
150W HPS	6.95	15.33	20"x 8"x 6"
250W HPS	8.033	17.7	20"x 8"x 6"
400W HPS	8.995	19.83	20"x 8"x 6"
430W HPS	9.745	21.48	20"x 8"x 6"
600W HPS	13.175	29.04	17"x 15"x 7"
750W HPS	13.32	29.36	17"x 15"x 7"
1000W HPS	16.07	35.42	17"x 15"x 7"
400W CONV.	12.66	27.91	17"x 15"x 7"
1000W CONV.	15.635	34.28	17"x 15"x 7"
250W MH	7.05	15.54	17"x 15"x 7"
400W MH	8.38	18.47	17"x 15"x 7"
1000W MH	12.395	27.32	17"x 15"x 7"

REFLECTOR HPS		
TYPE	X2 KG	BOX DIMS
MIDI	2	14"x 13"x 12"
MEDIUM	2	14"x 13"x 12"
WIDE	2	17"x 14"x 14"
SUPER WIDE	2	17"x 14"x 14"
DEEP	2	14"x 13"x 12"
1000W	2	17"x 14"x 14"
MAXIMA	2	19"x 14"x 14"

REFLECTOR MH		
TYPE	X2 KG	BOX DIMS.
MEDIUM	4	14"x 13"x 12"
DEEP	4	14"x 13"x 12"
1000W	4	17"x 14"x 14"

BULBS HPS & MH		
TYPE	PER CASE	KG
150W	12	3
250W	12	3
400W	12	3
430W	12	3
600W	12	3
750W	12	3
1000W	6	3

Customer satisfaction is very important to us. Please feel free to email your shipping related questions to us at <mailto:www.shipping@pllight.com>

## GREENHOUSE SHADING CURTAINS

In regards to all installations, P.L. Light Systems recommends that the installer contact both the local electrical authority having jurisdiction and the shading manufacturer for their specifications regarding issues pertaining to heat, rates of combustion and flame spread of their product(s). P.L. Light Systems also recommends that the Greenhouse Insurance Broker be contacted to limit potential liability.

For example the Ontario Electrical Code states:

### 30-1108 Spacing from Combustible Material

Conductors and lamp holders shall be maintained at a distance not less than 1 m from any combustible material, except for branch circuit conductors at the point of connection to bΔ 30-902 Spacing for Non-IC Type Luminaries.

ESA12363

## REFLECTOR CLEANING INSTRUCTIONS

### Instructions for cleaning reflectors

If the installation is in a clean environment or 1-3 years of age it is possible to clean the reflectors with a solution of water and vinegar at 1:100

If the installation is not a clean environment or beyond 3 years of age it may be possible to clean them with a solution of water and **RAYMAX No 5** at 1:120

Raymax advises cleaning your bulbs and reflectors annually for best performance.

### Work Instructions

1. Flush reflector with water inside and out to remove environmental contamination.
2. Fill basin with water and cleaning agent at ratios provided above. Submerge reflector and clean with a soft brush. *Do not submerge bulbs.*
3. Rinse reflector in a second basin of clean water to remove remaining cleaning solution.
4. Rinse reflector in third basin of distilled water to remove hard water residue. If residue is still present repeat steps 1 to 4.

For your safety, when cleaning reflectors always maintain the correct ratio and never add other chemicals to solution. Wearing hand and eye protection is advised.

## CHANGING THE VOLTAGE ON A BALLAST

1. Disconnect the power. Make sure that the unit is completely unplugged from the power outlet.
2. In a dry area, disassemble the housing of the fixture so that the wires are exposed. Read the label on the top of the ballast to ensure that the voltage to be changed to is present
3. Find the wires that run from your present cordset.
4. There will be a black wire off the cordset, attached to the present voltage on the ballast. A white wire off the cordset attached to the com wire on the ballast and a green ground wire off the cordset attached to the fixture housing.
5. Undo the 8/32 keps nut holding on the green wire and remove it. Cut off the amp splices on the identified wires.
6. Loosen off the cable strain relief holding the cord. Remove the old cordset and replace with new one meeting the requirements of the voltage and current of the fixtures new configuration.
7. Ensure that the new cordset is 4 to 6 inches past the strain relief and tighten down.
8. Find the green wire off the cordset and attach it to bolt in the fixture housing where the other green wire was removed. Reattach keps nut and tighten.
9. Find the original voltage wire from ballast and cap to preventing from shorting out. Ensure that no strands are exposed.
10. Find the new voltage wire from the ballast and strip off the end so that wires inside are showing about 1 cm or 1/2 inch.
11. Repeat step 10 for the COM wire you cut off cordset.
12. Take the black wire from the cordset and attach to new voltage wire from ballast. Take the white wire from the cordset and attach it to the COM wire on the ballast. Use rated amps or marrettes. **Do not use electrical tape.**

**Should be done by qualified personnel only.**

Changing the voltage on your fixture may void the warranty. Contact P.L. Light Systems if unsure.


# BALLAST INFORMATION

## BALLAST CIRCUITRY

The ballast in an HID lighting system has generally two purposes. 1) To provide the proper starting voltage to strike and maintain an arc; and 2) to provide the proper current to the lamp once the arc is established.

In addition to being designed to operate a particular type of HID lamp, a ballast design incorporates a basic circuitry to provide specific lamp/ballast operating characteristics. As an example, the effects of line voltage variations on resultant changes of lamp wattage are a function of the ballast circuit design. Requirements for a circuit which will provide a finer degree of lamp regulation generally result in higher ballast costs.

### Example of Wiring Diagrams.

**CWA**

**CWI**

Constant Wattage Autotransformer

Constant Wattage Isolated Transformer

## ADVANCE BALLASTS

<b>High Pressure Sodium Ballast Specifications</b>					
FIXTURE WATTAGE	VOLTAGE	WATTS	AMPS	LAMP TYPE	NEMA
<i>150 W CW</i>	120	190	1.7	150 W HPS	S55
	208	190	1.1	150 W HPS	S55
	240	190	0.8	150 W HPS	S55
<i>150 W CW</i>	120	190	2.8	150 W HPS	S55
	277	190	1.3	150 W HPS	S55
	347	190	0.9	150 W HPS	S55
<i>250 W CW</i>	120	295	2.7	250 W HPS	S50
	277	295	1.2	250 W HPS	S50
	347	295	0.9	250 W HPS	S50
<i>250 W CW</i>	120	295	2.5	250 W HPS	S50
	208	295	1.5	250 W HPS	S50
	240	295	1.3	250 W HPS	S50
	277	295	1.1	250 W HPS	S50
<i>400 W CW</i>	120	465	4.2	400 W HPS	S51
	208	465	2.4	400 W HPS	S51
	240	465	2.1	400 W HPS	S51
<i>400 W CW</i>	120	465	3.8	400 W HPS	S51
	277	465	1.7	400 W HPS	S51
	347	465	1.3	400 W HPS	S51
<i>430W CW</i>	120	490	5.2	430 W HPS	S51
	208	490	3	430 W HPS	S51
	240	490	2.6	430 W HPS	S51
	277	490	2.3	430 W HPS	S51
	347	490	1.8	430 W HPS	S51
<i>600 W CW</i>	480	490	1.3	430 W HPS	S51
	120	675	5.5	600 W HPS	S106
	208	675	3.3	600 W HPS	S106
	240	675	2.9	600 W HPS	S106
	277	675	2.5	600 W HPS	S106
<i>1000 W CW</i>	347	675	2	600 W HPS	S106
	480	675	1.4	600 W HPS	S106
	120	1100	9.5	1000W HPS	S52
	208	1100	5.5	1000W HPS	S52
	240	1100	4.8	1000W HPS	S52
<i>1000 W CW</i>	277	1100	4.2	1000W HPS	S52
	347	1100	2.5	1000W HPS	S52
	480	1100	1.78	1000W HPS	S52

These values are a general guideline to all ballasts produced by Advance Transformers. Values may vary according to availability at time of order and location of installation.

**ADVANCE BALLASTS**

Metal Halide Ballast Specifications					
FIXTURE WATTAGE	VOLTAGE	WATTS	AMPS	LAMP TYPE	NEMA
250 W CW	120	295	2.5	250 W MH	T15
	277	295	1.1	250 W MH	T15
	347	295	0.9	250 W MH	T15
250 W CW	120	295	2.5	250 W MH	T15
	208	295	1.4	250 W MH	T15
	240	295	1.3	250 W MH	T15
	277	295	1.1	250 W MH	T15
400 W CW	120	465	4.2	400 W MH	HPIT
	208	465	2.5	400 W MH	HPIT
	240	465	2.1	400 W MH	HPIT
400 W CW	120	460	4	400 W MH	HPIT
	277	460	1.7	400 W MH	HPIT
	347	460	1.4	400 W MH	HPIT
	480	460	1.4	400 W MH	HPIT
1000 W CW	120	1080	9	1000W MH	BT37
	208	1080	5.2	1000W MH	BT37
	240	1080	4.5	1000W MH	BT37
	277	1080	3.9	1000W MH	BT37
	347	1080	3.2	1000W MH	BT37
	480	1080	2.3	1000W MH	BT37

These values are a general guideline to all ballasts produced by Advance Transformers. Values may vary according to availability at time of order and location of installation.

**VENTURE BALLASTS**

<b>High Pressure Sodium Ballast Specifications</b>					
<b>FIXTURE WATTAGE</b>	<b>VOLTAGE</b>	<b>WATTS</b>	<b>AMPS</b>	<b>LAMP TYPE</b>	<b>NEMA</b>
<i>150 W CW</i>	120	195	3.2	150 W HPS	S55
	208	195	1.85	150 W HPS	S55
	240	195	1.6	150 W HPS	S55
<i>150 W CW</i>	120	195	3.2	150 W HPS	S55
	277	195	1.4	150 W HPS	S55
	347	195	1.1	150 W HPS	S55
<i>250 W CW</i>	120	290	2.7	250 W HPS	S50
	277	290	1.15	250 W HPS	S50
	347	290	0.95	250 W HPS	S50
<i>250 W CW</i>	120	290	2.7	250 W HPS	S50
	208	290	1.45	250 W HPS	S50
	240	290	1.25	250 W HPS	S50
	277	290	1.15	250 W HPS	S50
	347	290	1.15	250 W HPS	S50
<i>400 W CW</i>	120	465	3.8	400 W HPS	S51
	208	465	2.2	400 W HPS	S51
	240	465	1.9	400 W HPS	S51
<i>400 W CW</i>	120	465	3.8	400 W HPS	S51
	277	465	1.7	400 W HPS	S51
	347	465	1.3	400 W HPS	S51
<i>600 W CW</i>	120	675	6	600 W HPS	S106
	208	675	3.45	600 W HPS	S106
	240	675	3	600 W HPS	S106
	277	675	2.55	600 W HPS	S106
	347	675	2.05	600 W HPS	S106
<i>1000 W CW</i>	480	675	1.5	600 W HPS	S106
	120	1095	9.5	1000W HPS	S52
	208	1095	5.5	1000W HPS	S52
	240	1095	4.75	1000W HPS	S52
	277	1095	4.15	1000W HPS	S52
	347	1095	3.15	1000W HPS	S52
	480	1095	1.78	1000W HPS	S52

These values are a general guideline to all ballasts produced by Venture Transformers. Values may vary according to availability at time of order and location of installation.

**VENTURE BALLASTS**

Metal Halide Ballast Specifications					
FIXTURE WATTAGE	VOLTAGE	WATTS	AMPS	LAMP TYPE	NEMA
250 W CW	120	295	2.5	250 W MH	T15
	277	295	1.1	250 W MH	T15
	347	295	0.9	250 W MH	T15
250 W CW	120	295	2.5	250 W MH	T15
	208	295	1.45	250 W MH	T15
	240	295	1.25	250 W MH	T15
	277	295	1.1	250 W MH	T15
400 W CW	120	455	4.2	400 W MH	HPIT
	208	455	2.3	400 W MH	HPIT
	240	455	2	400 W MH	HPIT
400 W CW	120	455	4	400 W MH	HPIT
	277	455	1.75	400 W MH	HPIT
	347	455	1.45	400 W MH	HPIT
1000 W CW	120	1080	9	1000W MH	BT37
	208	1080	5.2	1000W MH	BT37
	240	1080	4.5	1000W MH	BT37
	277	1080	3.9	1000W MH	BT37
	347	1080	3.15	1000W MH	BT37
	480	1080	2.3	1000W MH	BT37

These values are a general guideline to all ballasts produced by Venture Transformers. Values may vary according to availability at time of order and location of installation.

## IGNITOR OPERATING LIFE

Under normal conditions, an ignitor actually operates for only a few cycles, once each day, when the lights are started. Under these conditions, the actual ignitor life expended is insignificant. Even if the lights were turned off momentarily, once a day, it requires only about 1 minute of pulsing by the ignitor to restrike the lamp. Assuming an operating temperature of 90C, an operating period of one minute a day would total only about 5 hours of operation per year. Since average ignitor life at 90C is a total of 800 hours, the use of 5 hours per year is only an insignificant portion of the total time.

However, ignitor life can be used up at a significant rate when an inoperative lamp remains in an energized socket for extended periods of time.

## IGNITOR MEASUREMENTS

### High Pressure Sodium Lamp Igniters

Igniters are used as a starting aid with all high pressure sodium and certain low wattage metal halide lamps. Measurement of the starting pulse voltage is beyond the capability of most instruments available in the field. In laboratory tests, an oscilloscope is used to measure the pulse height and width. In the field, some simple tests may be performed to determine if the ignitor is operable.

1. Replace the ignitor with one that is working. If the lamp then starts, the previous ignitor was either miswired or inoperative.
2. Remove the High Pressure Sodium lamp and replace with a known operable HPS lamp of the same Wattage. If the lamp lights, the ignitor is working properly and the bulb must be replaced.

**Warning!** This unit has the potential to destroy electrical measuring equipment and cause bodily harm. If unsure of any of the following directions, stop, and contact P. L. Light Systems before resuming work.

1. When testing for a faulty ignitor, the circuit **must** be turned off, and the ignitor isolated from the circuit.
2. Set your meter for capacitance.
3. Attach the black probe from the meter to the COM wire on the ignitor.
4. Attach the red probe from the meter to the LAMP wire on the ignitor.
5. The readings should be 3Mf +/- from the nominal values listed below.

<u>Advance Igniters</u>			<u>Venture Igniters</u>		
Model Number	Wattage	Nominal MF Value	Model Number	Wattage	Nominal MF Value
LI501-H4	400W	150	BVS-005	400W	205
LI501-JC	P.I. 1	150	BVS-005 P.I. 1	400W	180
LI501-J4	400 L.D. 2	150	BVS-038	600W	475
LI561-H5	600W	320	BVS-020	1000W	100
L1571-H5	1000W	85			

Note: 1. P.I designates Plug In Ignitor. 2. L.D. designates long distance ignitor.

## CAPACITOR INFO

All constant wattage Autotransformer, high power factor Reactor and Hi-Reactance ballast require a capacitor. With Core and Coil units, this capacitor is a separate component and must be properly connected to function properly.

There are two general types of capacitors currently used. Oil filled and dry type.

Oil filled capacitors today contain a non-PCB oil and are equipped with internal interrupters to prevent can rupture and resultant oil leakage in the event of a failure.

## IGNITOR / CAPACITOR CROSS REFERENCE GUIDE

### ADVANCE

<b>BALLAST</b>	<b>WATTS</b>	<b>VOLTAGE</b>	<b>HPS/MH</b>	<b>CAPS</b>	<b>UF</b>	<b>IGNITOR</b>
71A81E6	150	120/208/240	HPS	7C520P24	52MF	L1551J4
71A81A2	150	120/277/347	HPS	7C140M33R	14MF	L1551H4
71A8142	150	480	HPS	7C14M33R	14MF	L1551H4
71A82A1	250	120/277/347	HPS	7C35OP24	35MF	L1501H4
71A8291	250	120/208/240/277	HPS	7C35OP24	35MF	L1501H4
71A8241	250	480	HPS	7C35OP24	35MF	L1501H5
71A82E6	250	120/208/240	HPS	7C28OP30R	28MF	L1501J4
71A57A0	250	120/277/347	MH	7C150P40R	15MF	N/A
71A5790	250	120/208/240/277	MH	7C150P40R	15MF	N/A
71A84E6	400	120/208/240	HPS	7C48OP30	48MF	L1501J4
71A84A3	400	120/277/347	HPS	7C55OP24	55MF	L1501H4
71A8443	400	480	HPS	7C48OP30	48MF	L1501H4
71A60A1	400	120/277/347	MH	7C24OP40R	24MF	N/A
71A60E6	400	120/208/240	MH	MP2006100	20UF	N/A
71A6091	400	120/208/240/277	MH	7C24OP40R	24MF	N/A
71A6041	400	480	MH	7C24OP40R	24MF	N/A
71A9898	430	120/208/240	HPS	7C48OP30	48MF	L1501H4
71A85F5	430	277/347/480	HPS	7C48OP30	48MF	L1501H4
71A9942	600	120/208/240/277	HPS	7C66OS28	66MF	L1561H5
71A9943	600	277/347/480	HPS	7C66OS28	66MF	L1561H5
71A85E5	600	120/208/240	HPS	7C64OS28R	64MF	L1561H5
71A85F5	600	277/347/480	HPS	7C64OS28R	64MF	71561H5
71A86F5	750	277/347/480	HPS	7C75OS28R	75MF	L1561H5
71A86E5	750	120/208/240	HPS	7C75OS28R	75MF	L1561H5
71A87A3	1000	120/277/347	HPS	7C52OS28R	52MF	L1571H5
71A8793	1000	120/208/240/277	HPS	7C52OS28R	52MF	L1571H5
71A8743	1000	480	HPS	7C52OS28R	52MF	L1571H5
71A65A2	1000	120/277/347	MH	MD2409100	24UF	N/A
71A6592	1000	120/208/240/277	MH	MD2409100	24UF	N/A
71A6542	1000	480	MH	MD2409100	24UF	N/A

## IGNITOR / CAPACITOR CROSS REFERENCE GUIDE

### VENTURE



<b>BALLAST</b>	<b>WATTS</b>	<b>VOLTAGE</b>	<b>HPS/MH</b>	<b>CAPS</b>	<b>UF</b>	<b>IGNITOR</b>
V90D1435	150	120/208/240/277	HPS	ACG309	13.5	BVS-005-A
V90J1435	150	120/277/347	HPS	ACG309	14	BVS-006-A
V90J1711	250	120/277/347	HPS	ACG172	35	BVS-005-A
V90D1711	250	120/208/240/277	HPS	ACG172	35	BVS-005-A
V90J6212	250	120/277/347	MH	ACH277	15	N/A
V90D6212	250	120/208/240/277	MH	ACH277	15	N/A
V90J1911	400	120/277/347	HPS	ACG255	55	BVS-005-A
V90D1912	400	120/208/240/277	HPS	ACG255	55	BVS-005-A
V90D9410	400	120/208/240/277	CONV		33/25	BVS-045
V90J6411	400	120/277/347	MH		24	N/A
V90D6413	400	120/208/240/277	MH	ACH247	24	N/A
V90C2150	600	120/208/240	HPS	ACG302	45	BVS-038
V90AA2150	600	277/347/480	HPS	ACG302	45	BVS-038
V90J2310	1000	120/277/347	HPS	ACB264OV	26	BVS-020
V90D2311	1000	120/208/240/277	HPS	ACB264OV	26	BVS-020
V90D9511	1000	120/208/240/277	CONV		26	BVS-020
V90D6514	1000	120/208/240/277	MH	ACB220OV	24	N/A
V90DJ6514	1000	120/277/347	MH	ACB220OV	24	N/A

# 1 ZAWG STRUT INFORMATION

## P1000

Wt/100 Ft: 190 Lbs (283 kg/100 m)  
 Allowable Moment: 5,080 In-Lbs (570 N•m)  
 12 Gauge Nominal Thickness: .105" (2.7mm)

Channel Finishes:  
 PL, GR, HG, PG;  
 Standard Lengths:  
 10' & 20'

### COLUMN LOADING – P1000

Height In	Maximum Column Load Applied at C.G.				
	Unbraced Allowable Load at Slot Face Lbs	K = 0.65 Lbs	K = 0.80 Lbs	K = 1.0 Lbs	K = 1.2 Lbs
24	3,450	10,750	9,900	8,770	7,730
36	3,050	8,910	7,730	6,370	5,280
48	2,660	7,250	5,980	4,660	3,770
60	2,290	5,890	4,660	3,600	2,940
72	2,000	4,800	3,770	2,940	2,380
84	1,760	4,010	3,170	2,460	1,970
96	1,570	3,450	2,730	2,090	1,650
108	1,410	3,020	2,380	1,800	**
120	1,270	2,680	2,090	**	**

### BEAM LOADING – P1000

Span In	Max Allowable Uniform Load Lbs	Defl. at Uniform Load In	Uniform Loading at Deflection		
			Span/180 Lbs	Span/240 Lbs	Span/360 Lbs
24	1,690	0.06	1,690	1,690	1,690
36	1,130	0.13	1,130	1,130	900
48	850	0.22	850	760	510
60	680	0.35	650	490	320
72	560	0.50	450	340	220
84	480	0.68	330	250	170
96	420	0.89	250	190	130
108	380	1.13	200	150	100
120	340	1.40	160	120	80
144	280	2.01	110	80	60
168	240	2.74	80	60	40
192	210	3.57	60	50	NR
216	190	4.52	50	40	NR
240	170	5.58	40	NR	NR

### MATERIAL

Unistrut channels are accurately and carefully cold formed to size from low-carbon strip steel. All spot-welded combination members, except P1001T, are welded 3" (76 mm) maximum on center.

#### STEEL: PLAIN

12 Ga. (2.7 mm), 14 Ga. (1.9 mm) and 16 Ga. (1.5 mm)  
 ASTM A1011 GR33

#### STEEL: PRE-GALVANIZED

12 Ga. (2.7 mm), 14 Ga. (1.9 mm) and 16 Ga. (1.5 mm)  
 ASTM A653 GR 33

For other materials, see Special Metals or Fiberglass sections.

### FINISHES



All channels are available in:

- Perma Green II (GR)
- Pre-galvanized (PG), conforming to ASTM A653 G90
- Hot-dipped galvanized (HG), conforming to ASTM A123
- Plain (PL)

# 14AWG STRUT INFORMATION

## P1100

Wt/100 Ft: 142 Lbs (211 kg/100 m)  
 Allowable Moment: 4,170 In-Lbs (470 N•m)  
 14 Gauge Nominal Thickness: .075" (1.9mm)

Channel Finishes:  
 PL, GR, HG, PG;  
 Standard Lengths:  
 10' & 20'

### COLUMN LOADING – P1100

Unbraced Height In	Maximum Allowable Load at Slot Face Lbs	Maximum Column Load Applied at C.G.			
		K = 0.65 Lbs	K = 0.80 Lbs	K = 1.0 Lbs	K = 1.2 Lbs
24	2,770	8,120	7,450	6,540	5,660
36	2,410	6,650	5,660	4,480	3,520
48	2,040	5,240	4,140	3,040	2,390
60	1,690	4,050	3,040	2,270	1,830
72	1,440	3,140	2,390	1,830	1,480
84	1,260	2,560	1,980	1,530	1,240
96	1,120	2,170	1,700	1,310	1,060
108	1,000	1,880	1,480	1,140	**
120	910	1,670	1,310	**	*

### BEAM LOADING – P1100

Span In	Max Allowable Uniform Load Lbs	Defl. at Uniform Load In	Uniform Loading at Deflection		
			Span/180 Lbs	Span/240 Lbs	Span/360 Lbs
24	1,390	0.06	1,390	1,390	1,390
36	930	0.13	930	930	720
48	700	0.23	700	610	410
60	560	0.36	520	390	260
72	460	0.51	360	270	180
84	400	0.70	270	200	130
96	350	0.91	200	150	100
108	310	1.15	160	120	80
120	280	1.42	130	100	70
144	230	2.05	90	70	50
168	200	2.79	70	50	30
192	170	3.65	50	40	30
216	150	4.62	40	30	NR
240	140	5.70	30	NR	NR

### MATERIAL

Unistrut channels are accurately and carefully cold formed to size from low-carbon strip steel. All spot-welded combination members, except P1001T, are welded 3" (76 mm) maximum on center.

#### STEEL: PLAIN

12 Ga. (2.7 mm), 14 Ga. (1.9 mm) and 16 Ga. (1.5 mm)  
 ASTM A1011 GR33

#### STEEL: PRE-GALVANIZED

12 Ga. (2.7 mm), 14 Ga. (1.9 mm) and 16 Ga. (1.5mm)  
 ASTM A653 GR 33

For other materials, see Special Metals or Fiberglass sections.

### FINISHES


All channels are available in:

- Perma Green II (GR)
- Pre-galvanized (PG), conforming to ASTM A653 G90
- Hot-dipped galvanized (HG), conforming to ASTM A123
- Plain (PL)



## OUTLET BOXES, END CAPS, CLOSURE STRIPS

### Receptacle Box - mounts on strut

- Reduces installation time
- Cleaner installation ( no feeders showing)
- Approved by UL and CSA


### Closure Strip – Modifies Strut into an approved Raceway


Finish: Green, pre-galvanized, plain.  
Standard length: 10' (3m).

### End Cap for Strut Raceway System

### 1/2\" data-bbox="142 634 750 653"/>

Call us for more products, details and how this system can work for you!

# REFLECTOR INFORMATION

## Reflector Overview

### The Most Important Part Of A Greenhouse Fixture.

With more reflectors than any other supplemental lighting supplier, P.L. Light Systems sets itself apart from its competitors through their capability to direct light to achieve the best possible light intensity and uniformity for your plants with the fewest number of fixtures.

#### Reflectors:

#### 60° 1000W REFLECTOR

- Reflection angle of 0° - 60°
- Designed for a more rectangular illumination field
- Best performances at a mounting height at 8' or higher depending on light level required
- Deep penetration of light into the crop

#### 60° DEEP REFLECTOR

- Reflection angle of 0 - 60°
- Designed for a square illumination field
- Deep penetration of light into the crop
- Suitable for a height of 5' depending on the light level required

#### 65° MIDI REFLECTOR

- Reflection angle of 0 - 65°
- Specially designed for the modern greenhouse to allow you to work with one row per bay in 26' and 32' bays
- For suspension heights over 6' depending on the required light level
- Remarkable performance in 21' bays as well

#### 70° MEDIUM REFLECTOR

- Reflection angle of 0° -70°
- Designed for a more rectangular illumination field
- Suitable for mounting heights over 6' with a maximum distance of 21' between rows, depending on the light level required

#### 70° WIDE REFLECTOR


- Reflection angle of 0 - 70°
- Designed for a more rectangular illumination field
- Suitable for mounting heights between 4' and 5' depending on the light level required
- May be equipped with one of two mirror diaphragms to optimize light distribution

#### 75° SUPERWIDE REFLECTOR

- Reflection angle of 0 - 75°
- Designed for a more rectangular illumination field
- Suitable for very low mounting heights between 3' and 4' depending on the light level that is required
- Provided with a standard screen to optimize light distribution
- Commonly used for photo-period lighting

#### MAXIMA REFLECTOR

- Designed for a more square illumination field
- Suitable for low mounting heights
- Suitable for higher light levels

#### Warranty Period

FIXTURE COMPONENTS ONLY:  
5 Years from Date of Purchase

LAMPS: 1 Year from Date of Purchase

P.L. Light Systems Canada Inc. will determine warranty eligibility of returned defective items and provide the buyer with an RMA number.

P.L. LIGHT SYSTEMS






For technical specifications or other information on P.L. Light Systems products, please visit

[www.pllight.com](http://www.pllight.com) or call us at

905.563.4133 or 1.800.263.0213

# RECEPTACLES

## STRAIGHT BLADE

DESCRIPTION		15 AMPERE	20 AMPERE	30 AMPERE	50 AMPERE	60 AMPERE	
		RECEPTACLE	RECEPTACLE	RECEPTACLE	RECEPTACLE	RECEPTACLE	
2-POLE 3-WIRE GROUNDING	125 V	5	5-15R 	5-20R 	5-30R 	5-50R 	
	125 V	5A		5-20RA ALTERNATE 			
	*250 V	6	6-15R 	6-20R 	6-30R 	6-50R 	
	*250 V	6A		6-20RA ALTERNATE 			
	277 V AC	7	7-15R 	7-20R 	7-30R 	7-50R 	
	347 V AC	24	24-15R 	24-20R 	24-30R 	24-50R 	
3-POLE 4-WIRE GROUNDING	125/250 V	14	14-15R 	14-20R 	14-30R 	14-50R 	14-60R 
	3 ∅ 250 V	15	15-15R 	15-20R 	15-30R 	15-50R 	15-60R 

\*For configurations 6-15R, 6-20R, 6-20RA, 6-30R, and 6-50R, Y denotes identified terminal when used on circuits derived from 3-phase, 4-wire 416 V circuits.

Note: Except as noted above, in Diagrams 1 and 2.

(a) G represents the terminal for bonding to ground;

(b) W represents the identified terminal; and

(c) X, Y, and Z represent the terminals for ungrounded conductors.

***When ordering plugs for your fixtures please denote the following criteria.***

- *Straight blade or twist lock*
- *Voltage of circuit fixture to be installed on*
- *Amperage of circuit fixture to be installed on*
- *Length of cord required for easy installation*

# RECEPTACLES

## TWIST LOCK

DESCRIPTION		15 AMPERE	20 AMPERE	30 AMPERE	50 AMPERE	60 AMPERE
		RECEPTACLE	RECEPTACLE	RECEPTACLE	RECEPTACLE	RECEPTACLE
2-POLE 3-WIRE GROUNDING	125 V	L5 	L5-20R 	L5-30R 	L5-50R 	L5-60R 
	250 V	L6 	L6-20R 	L6-30R 	L6-50R 	L6-60R 
	277 V AC	L7 	L7-20R 	L7-30R 	L7-50R 	L7-60R 
	480 V AC	L8 	L8-30R 	L8-50R 	L8-60R 	
	600 V AC	L9 	L9-30R 	L9-50R 	L9-60R 	
3-POLE 4-WIRE GROUNDING	125/250 V	L14 	L14-30R 	L14-50R 	L14-60R 	
	3 Ø 250 V	L15 	L15-30R 	L15-50R 	L15-60R 	
	3 Ø 480 V	L16 	L16-30R 	L16-50R 	L16-60R 	
	3 Ø 600 V	L17 	L17-50R 	L17-60R 		
4-POLE 5-WIRE GROUNDING	3 Ø 208 Y / 120 V	L21 	L21-30R 	L21-50R 	L21-60R 	
	3 Ø 480 Y / 277 V	L22 	L22-30R 	L22-50R 	L22-60R 	
	3 Ø 600 Y / 347 V	L23 	L23-30R 	L23-50R 	L23-60R 	

Note: Except as noted above, in Diagrams 1 and 2.  
 (a) G represents the terminal for bonding to ground;  
 (b) W represents the identified terminal; and  
 (c) X, Y, and Z represent the terminals for ungrounded conductors.


**When ordering plugs for your fixtures please denote the following criteria.**

- Straight blade or twist lock
- Voltage of circuit fixture to be installed on
- Amperage of circuit fixture to be installed on
- Length of cord required for easy installation

**Info grower:**

Company name \_\_\_\_\_  
 Contact person \_\_\_\_\_  
 Correspondence address \_\_\_\_\_  
 Postcode / City / Country \_\_\_\_\_  
 Mobile tel. \_\_\_\_\_  
 Telephone / Fax \_\_\_\_\_  
 E-mail \_\_\_\_\_  
 Address Nursery \_\_\_\_\_  
 Postcode / City / Country \_\_\_\_\_  
 Telephone / Fax \_\_\_\_\_  
 Dealer / contactperson \_\_\_\_\_

P.L. Light Systems  
 4800 Hlan Dr.  
 Beamsville, ON  
 Canada L0R 1B1  
 Tel 001 9055634133  
 Fax 001 9055630445  
**WWW.PLLIGHT.COM**

Please sketch in any spaces, e.g. the location of the shed, basin, CHP installation and paths in the plan. In the case of more complex spaces, please include the plan.

Heating pipes, irrigation pipes etc. affect the uniformity and lightlevel. Please sketch in the exact locations and their numbers on **page 2**.

- (A) Nbr. of bays \_\_\_\_\_
- (B) Bay size \_\_\_\_\_
- (C) Total width \_\_\_\_\_
- (D) Nbr. of sections \_\_\_\_\_
- (E) Section size \_\_\_\_\_
- (F) Total length \_\_\_\_\_

	Yes	No		Yes	No
Screen	<input type="radio"/>	<input type="radio"/>	Newly built greenhouse	<input type="radio"/>	<input type="radio"/>
Heating pipes	<input type="radio"/>	<input type="radio"/>	CHP installation	<input type="radio"/>	<input type="radio"/>
Irrigation pipes	<input type="radio"/>	<input type="radio"/>	Energy from the mains	<input type="radio"/>	<input type="radio"/>
Other obstacles	<input type="radio"/>	<input type="radio"/>	Fixtures already present	<input type="radio"/>	<input type="radio"/>
If so, what kind: _____			If so, what type: _____		


The lightplan is designed on the basis of the available information. We are not liable for any lacking data. Please enter the locations and the number of obstacles on **page 2**.

To the lightplan (drawing + computer calculation) the "Explanation to the Lightingplan", the Metaalunie conditions and supplementary guarantee provisions of P.L. Light Systems apply.

Company name: \_\_\_\_\_


P.L. Light Systems  
 4800 Hlnan Dr.  
 Beamsville, ON  
 Canada L0R 1B1  
 Tel 001 9055634133  
 Fax 001 9055630445  
**WWW.PLLIGHT.COM**

- (G) Gutterheight \_\_\_\_\_
- (H) Height of concrete pile \_\_\_\_\_
- (I) Gutter to truss \_\_\_\_\_
- (J) Truss height \_\_\_\_\_
- (K) Underside truss to soil \_\_\_\_\_

Crop  
 Top of crop - soil \_\_\_\_\_ ( mtr )  
 Lightlevel \_\_\_\_\_ ( lux )

**Example**

Heating pipes, Irrigation pipes etc. affect the uniformity and lightlevel. Please sketch in the exact locations and their numbers.

- Type of fixture:
- |                                       |  |
|---------------------------------------|--|
| PL2000 1000Watt <input type="radio"/> | HS2000 600Watt <input type="radio"/>   |
| 600Watt <input type="radio"/>         | 400Watt <input type="radio"/>          |
| 400Watt <input type="radio"/>         | HSE600 600Watt <input type="radio"/>   |
|                                       | HSE1000 1000Watt <input type="radio"/> |

The lightplan is designed on the basis of the available information. We are not liable for any lacking data. Please enter the locations and the number of obstacles above.

To the lightplan (drawing + computer calculation) the "Explanation to the Lightingplan", the Metaalunie Conditions and supplementary guarantee provisions of P.L. Light Systems apply.

Submit your **P.L. Light Plan** today!

## **TROUBLE SHOOTING GUIDE**

The following questions may provide one or more solutions to fixture problems.

Q: WHY WILL MY FIXTURE NOT LIGHT?

- A: 1. Check to see if there is power to the ballast  
2. The components may be faulty.  
3. The bulb may be faulty  
4. You may have the wrong ballast/ bulb combination check to see if the bulb and ballast are the same type HPS or MH.  
5. Make sure that the fixture has been properly wired it may have defective or improper wiring.

Q: WHY DOES MY FIXTURE HAVE 50% LUMINATION?

- A: 1. The capacitor may be faulty.  
2. The voltage to the fixture may not coincide with voltage used.  
3. The ballast type may not be the same as the bulb type.  
4. Your bulb could just be dirty.

Q: HOW DO I KNOW IF THE BULB IS BLOWN?

- A: Take the bulb out and replace it with a new bulb, if it lights then it was the bulb if it does not light then there is different problem.

Q: HOW DO I KNOW IF THE IGNITOR IS BLOWN?

- A: If it is possible take out the bulb and replace it with a metal halide bulb, if it fires then it is the ignitor.

Q: THE REPLACEMENT PARTS I RECEIVED HAS A PLASTIC CASING BUT THE ONE I TOOK OUT IS METAL IS THAT OK?

- A: The replacement parts are ok. Just make sure that the part numbers are the same.

Q: WHY IS MY BULBS STARTING SHOW SIGNS OF SCORHING EARLY?

- A: The normal end of life is indicated by the scorching or blackening at the ends of the arc tube. If it turns starts to go black earlier that could be a good indication that you may have a faulty capacitor.

Q: WHAT DOES EVERY FIXTURE COME WITH?

- A: Every fixture comes with an endplate, bulb, reflector, and hanging brackets with nuts and bolts. Cord is optional.

## **P.L. LIGHT SYSTEMS**

*All your Hard Work. Your Livelihood. Your Pride. The perfect crop. A better season. Higher yields. All the time you spent searching for the ideal greenhouse, irrigation system, fertilizer, pesticide; now your plants are under performing because their not receiving enough light... not very amusing is it? **A simple quiz. You knew the answer.***

P.L. Light Systems have been servicing the North American and South American Greenhouse Supplementary Lighting markets since 1981. Our reputation goes into the products our customers rely on every day for their livelihood and the health and well being of their crops. That why every light fixture and every reflector is tested to ISO 9000/2001 standards before leaving our facilities.

Through longstanding integral relationships with our European Commercial Divisions; Hortilux [www.hortilux.com](http://www.hortilux.com) and our R&D facility in The Netherlands, Raymax [www.raymax.nl](http://www.raymax.nl) our reflectors are tested for consistency in reflectivity, low lumen losses, greater lux output and pattern for the ultimate performance on the market. Our guarantee to our customer is that every reflector you purchase meets the specifications of P.L. Light Systems has been imported from the Netherlands to maintain a consistent level of quality and performance. We know that the best performing reflector on the market needs the best performing components to maximize yield, that's why the components we surround our reflector with are provided through committed relationships with Philips Lighting [www.philipslighting.com](http://www.philipslighting.com) and Advance Transformers. [www.advancetransformer.com](http://www.advancetransformer.com)

The success of supplemental lighting depends on other environmental parameters such as temperature, relative humidity, carbon dioxide concentrations, crop nutrition and plant architecture. All these parameters need to be favorable before supplementary lighting can be optimized to growing better plants. For this reason P.L. Light Systems works with an intelligent distributor network that really understands growing and can help with customer questions and concerns. Call P.L. Light Systems and find out what we can do for you and your business.

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
# GRUNDFOS COMFORT

## Circulator pump

GRUNDFOS COMFORT circulator pumps are available in 2 pump housing versions and lengths incorporating isolating and non-return valves or prepared for subsequent fitting of such valves.

The water-conduction part of the pump is hermetically separated from the stator with a stainless steel spherical separator.

The motor can be separated from the pump housing, enabling easy maintenance and replacement.

## Applications

- Domestic hot water systems in single and two-family houses
- Small heating systems
- Cooling and air-conditioning systems

## Options

### 24-hour timer

The timer is built into the pump. The table below shows the functions of the timer.

Timer set to...	The pump...
OFF	is switched off.
TIMER	starts and stops automatically at set intervals of minimum 20 minutes.
ON	runs continuously.

### Adjustable thermostat

The built-in thermostat of the pump types BT, BUT, BXT, BXUT can be set to stop the pump at a preset liquid temperature.

Setting range: 35-65°C.

The thermostat function can be interrupted by turning the thermostat to the position .

Factory setting: 35°C.

### Various fittings

- Extension sets
- Non-return valve and shut-off valve
- Union sets
- Venting flange

TM01 9237 1500

## Pumped liquids

- Thin, clean, non-aggressive and non-explosive liquids without solid particles or fibres.
- Cooling liquids, not containing mineral oil.
- Domestic hot water.
- Softened water.

The kinematic viscosity of water is 1 mm<sup>2</sup>/s (1 cSt) at 20°C. If the circulator pump is used for a liquid with a higher viscosity, the hydraulic performance of the pump will be reduced.

Example: 50% glycol at 20°C means a viscosity of approx. 10 mm<sup>2</sup>/s and a reduction of pump performance by approx. 15%.

When selecting a pump, the viscosity of the pumped liquid must be taken into consideration.

## Ambient and liquid temperatures

Liquid temperature: +2°C to +95°C

It is recommended to keep the operating temperature as low as possible (e.g. 65°C) to avoid lime precipitation.

The ambient temperature should always be lower than the liquid temperature, as otherwise condensation may form in the stator housing.

## Maximum system pressure

PN 10: 1.0 MPa (10 bar).

## Inlet pressure

To avoid cavitation noise and damage to the pump bearing, at high temperatures the following minimum pressures are required at the pump suction port.

Liquid temperature	85°C	95°C
Inlet pressure	0.5 m head	2.8 m head
	0.049 bar	0.27 bar



**GRUNDFOS**

### UP 15-14

80

1 x 230 V, 50 Hz

<b>P<sub>1</sub> [W]</b>	<b>I<sub>n</sub> [A]</b>
25	0.11

The motor incorporates thermal overload protection.



Connections: Various fittings are available  
 System pressure: Max. 10 bar  
 Liquid temperature: +2°C to +95°C (TF 95)

Pump type	Dimensions [mm]									Weights [kg]		Ship. vol. [m³]
	L1	L2	H1	H2	H3	H4	B1	B2	G	Net	Gross	
UP 15-14 B	80		25	13.5	133		79.5	84	Rp ½	1.00	1.12	0.0026
UP 15-14 BU	80	90	25	13.5		205	79.5	84	Rp ½	1.15	1.31	0.0034
UP 15-14 BT	80		25	13.5	155		79.5	84	Rp ½	1.05	1.24	0.0034
UP 15-14 BUT	80	90	25	13.5		205	79.5	84	Rp ½	1.16	1.32	0.0034

### UP 20-14

110

1 x 230 V, 50 Hz

<b>P<sub>1</sub> [W]</b>	<b>I<sub>n</sub> [A]</b>
25	0.11

The motor incorporates thermal overload protection.

Connections: Various fittings are available  
 System pressure: Max. 10 bar  
 Liquid temperature: +2°C to +95°C (TF 95)

Pump type	Dimensions [mm]								Weights [kg]		Ship. vol. [m³]
	L1	H1	H2	H3	H4	B1	B2	G	Net	Gross	
UP 20-14 BX	110	25	21	133		79.5	84	G 1¼	1.20	1.35	0.0026
UP 20-14 BXU	110	25	21		205	79.5	84	G 1¼	1.35	1.51	0.0034
UP 20-14 BXT	110	25	21	155		79.5	84	G 1¼	1.25	1.44	0.0034
UP 20-14 BXUT	110	25	21		205	79.5	84	G 1¼	1.36	1.52	0.0034

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
**IC200UDR001**

**14 Point Micro PLC, (8) 24VDC In, (6) Relay Out, 120/240VAC Power Supply**

VersaMax Micro PLC IC200UDR001 accepts eight DC inputs and provides six normally-open 2 Amp relay outputs that can control 5-30VDC or 5-250VAC output devices.

### Features

- 100VAC to 240VAC nominal input power.
- Eight configurable DC inputs that can be used as positive or negative logic standard inputs or High-speed Counter inputs.
- Six Form A relay outputs (SPST-single pole single throw).
- High-speed Counter inputs can be configured as 4 Type A Counters or 1 Type A and 1 Type B Counter.
- +24VDC output power available for field devices, up to 200mA maximum.
- Expandable to 70 points using up to four 14-point Expansion Units.
- Two removable screw-down “barrier-style” terminal strips with protective covers.
- RS-232 serial port that supports SNP/SNPX, RTU slave protocols, and Serial I/O.
- Run/Stop mode switch that can be configured as a run/stop switch, a memory protect switch, and also used for clearing faults when a fatal fault exists.
- Two analog potentiometers.
- Full-featured programming Instruction Set with floating point math The application program can be either Ladder Diagram (LD) or Instruction List (IL) format.
- 9K words of program memory, 256 words of registers.
- Flash memory (ROM) for non-volatile program storage and for system firmware.
- Configurable to read configuration at powerup from either RAM or flash memory (ROM). Can also be configured to read application program from flash at powerup.
- Capacitor backs up RAM for at least 3 days.

**IC200UDR001**  
**14 Point Micro PLC, (8) 24VDC In, (6) Relay Out, 120/240VAC Power Supply**

**Micro PLC IC200UDR001 Specifications**

Weight	380 grams (0.84lb)
Module Dimensions	Height: 90mm (3.6 inches) Depth: 76mm (3.0 inches) Width: 95mm (3.8 inches)
Typical Scan Rate	1.1 ms/K for Boolean logic (release 1.10 or later)
Inputs	Eight 24 VDC positive/negative logic input circuits
Outputs	Six normally open 2 Amp relay circuits
High-speed Counters	Up to four Type A or one Type A and one Type B
Maximum number of slave devices per RS-485 network	8 (can be increased with a repeater). Requires IC200ACC415.
Output Power Supplies	24VDC for input circuits & user devices, 200mA max. +5VDC on pin 7 of Serial Port, 100mA max.
Realtime clock accuracy (for Timer contacts and Timer function blocks)	+/-0.5%

**Immunity and Emission Specifications**

Description	Standards	Specifications
Electrostatic Discharge	EN 61000-4-2	± 4.0 kV (Contact) ± 8.0 kV (Air)
RF Susceptibility	EN 61000-4-3	10 V/m (unmodulated), 80-1000 MHz, 80% AM, 1 kHz sine wave
RF Susceptibility from Digital Radio Telephones	ENV 50204	10 V/m (unmodulated), 900±5Mhz, 100% AM 200 Hz square wave, 50% duty cycle)
Fast Transient	EN 61000-4-4	± 2.2 kV (PS), ± 1.1 kV (I/O), ± 2.2 kV (RS232), ± 2.2 kV (RS485)
Voltage Surge	EN 61000-4-5	± 2.2 kV, common mode (PS) ± 1.1 kV, differential (PS)
Conducted RF	EN 61000-4-6	10 V <sub>rms</sub> , 0.15-80 MHz, 80% AM, 1 kHz sine wave (PS, I/O)
Voltage Dip Voltage Interrupt Voltage Variation	EN 61000-4-11	30% Nom., 10 ms >95% Nom., 10ms 10% Nom. 10 sec.
Radiated Emissions	EN 55011*	30 dBμV/m, 30 – 230 MHz (measured @ 30m) 37 dBμV/m, 230 – 1000 MHz (measured @ 30m)
Conducted Emissions	EN 55011*	79/66 dBμV, 0.15 – 0.5 MHz 73/60 dBμV, 0.5 – 30 MHz

\* EN 55011 limits are equivalent to limits specified in EN 55022, CISPR 11, CISPR 22, and 47 CFR 15.

**IC200UDR001**

**14 Point Micro PLC, (8) 24VDC In, (6) Relay Out, 120/240VAC Power Supply**

### AC Power

#### AC Power Specifications

Range Frequency	100 -15% to 240 +10% VAC 50 -5% to 60 +5% Hz
Hold-up	10mS at 85 to 100VAC. 20mS at 100 to 264VAC
Inrush Time	2mS for 40A
Inrush Current	18 Amp maximum at 120 VAC 30 Amp maximum at 200 VAC 40 Amp maximum at 265 VAC
Input Current	0.06 Amp typical at 200 VAC 0.10 Amp typical at 100 VAC
Input Power Supply Rating	13 VA
Isolation	1500VAC RMS field-side to logic (power supply input).

### Inputs

Whether used as a standard or HSC input, each input can have positive or negative logic characteristics. Current into an input point results in a logic 1 in the input status table (%I). When used as standard inputs, the input characteristics are compatible with a wide range of input devices, such as pushbuttons, limit switches, and electronic proximity switches.

#### DC Input Specifications

Number of Inputs	8
Rated Input Voltage	24 volts DC
Input Voltage Range	0 to 30 volts DC
Input Current	7.5mA typical
Input Resistance	2.8 Kohms
Input Threshold Voltage	ON 15VDC minimum OFF 5VDC maximum
Input Threshold Current	ON 4.5mA maximum OFF 1.5mA minimum
Response Time	0.5 to 20ms configurable as regular input; 100 $\mu$ s as HSC input
Isolation Voltage	1500V RMS field side to logic side 500V RMS between groups

The Micro PLC provides a +24 VDC supply that can be used for input devices and to power the DC input circuits at about 7.5 mA per input. The combination of input circuit current and external device current must not exceed 200 mA.

**IC200UDR001**  
**14 Point Micro PLC, (8) 24VDC In, (6) Relay Out, 120/240VAC Power Supply**

**Outputs**

The six normally-open relay outputs can control a wide range of devices such as motor starters, solenoids, and indicators. Power for the internal relay coils is provided by the internal supply. An external source of AC or DC power must be supplied to operate field devices.

Connections and specifications for HSC outputs are the same as for standard relay outputs.

**Relay Output Specifications**

Operating Voltage	5 to 30 VDC or 5 to 250 VAC		
Isolation	1500 V RMS between field side and logic side 500 V RMS between groups		
Leakage Current	15 mA at 240 VAC maximum		
Maximum UL Pilot Duty Rating	2 amps at 24 VDC and 240 VAC		
Maximum Resistive Load Rating	2 amps at 24 VDC and 240 VAC		
Minimum Load	10 mA		
Maximum Inrush	5 amps per half cycle		
On Response Time	15 ms maximum		
Off Response Time	15 ms maximum		
Fuse	None		
Contact Life: Mechanical	20 x 10 <sup>6</sup> mechanical operations		
Contact Life: Electrical Voltage 240VAC, 120VAC, 24VDC	Current: Resistive 2A	Current: Lamp and Solenoid 0.6A	Typical Operations 200,000

**IC200UDR001**

**14 Point Micro PLC, (8) 24VDC In, (6) Relay Out, 120/240VAC Power Supply**

### **High Speed Counters**

VersaMax Micro PLC UDR001 can be configured to provide built-in high-speed counter operation.

When configured for High-speed Counter operation, inputs can be set up as:

- Up to four Type A Counters or
- One Type A and one Type B Counter.

Each counter provides direct processing of rapid pulse signals up to 10kHz for industrial control applications such as meter proving, turbine flowmeter, velocity measurement, material handling, motion control, and process control.

Each counter can be enabled independently. Type A counters can be configured for up or down counting (default is up) and for positive or negative edge detection (default is positive). The type B counter provides an A Quad B counting function.

The relay outputs can be configured as up to four counter outputs. They cannot be used as Pulse Train or PWM outputs.


### **High-speed Counter Specifications**

Available High Speed Counters	Four Type A or One Type A and One Type B
Maximum Counter Frequency	10kHz
Input Voltage	ON     15V OFF    5V
Count Pulse Width	20% to 80% duty cycle
Count Registers	16 bits
<b>Outputs</b>	
Available Pulse Outputs	None
Load Voltage	Refer to relay output specifications

**IC200UDR001**  
**14 Point Micro PLC, (8) 24VDC In, (6) Relay Out, 120/240VAC Power Supply**

**Field Wiring and Circuit Diagrams**

\* When I1-18 are used as High-speed Counter inputs, the input switches should be solid state to prevent switch bouncing, which could cause unintended high speed counter counts or strobe input signals.

**Typical 24 VDC Positive/Negative Logic Input Circuit**

\* Positive connection shown; reverse polarity of 24VDC power supply connections for negative connection.

The diagram shows a terminal strip with a 24VDC source connected to a COM terminal and an I/O terminal. The I/O terminal is connected to a high-frequency filter, a 2.8kΩ resistor, and a diode. The other side of the diode is connected to a 5V supply and an LED. The LED is connected to the CPU. The COM terminal is also connected to other circuits.

**Typical Relay Output Circuit**

The diagram shows a terminal strip with a 24V source connected to a COM terminal and an I/O terminal. The I/O terminal is connected to a 24V source through a 2.8kΩ resistor and a diode. The diode is connected to a 5V supply and an LED. The LED is connected to the CPU. The COM terminal is also connected to other circuits. The I/O terminal is connected to a relay coil, which is connected to a 24V source. The relay contact is connected to a terminal strip, which is connected to a fuse, a load (LOAD), and a common terminal. The common terminal is connected to a power source.

**High Speed Counter Circuit - Positive Logic Connection**

The diagram shows a terminal strip with a 24VDC source connected to a COM terminal and a Count or Strobe terminal. The Count or Strobe terminal is connected to a high-frequency filter, a 2.8kΩ resistor, and a diode. The other side of the diode is connected to a 5V supply and an LED. The LED is connected to the CPU. The COM terminal is also connected to other circuits.

**High Speed Counter Circuit - Negative Logic Connection**

The diagram shows a terminal strip with a 24VDC source connected to a COM terminal and a Count or Strobe terminal. The Count or Strobe terminal is connected to a high-frequency filter, a 2.8kΩ resistor, and a diode. The other side of the diode is connected to a 5V supply and an LED. The LED is connected to the CPU. The COM terminal is also connected to other circuits.

**DC Load**

The diagram shows a relay output connected to a DC supply. A 1A, 100V diode is connected in parallel with the load to suppress inductive spikes.

**AC Load**

The diagram shows a relay output connected to an AC supply. A .022μf 100Ω capacitor is connected in parallel with the load to suppress inductive spikes.

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**96.1**

## ***APPENDIX 7***

This document is confidential property of the MELISSA partners and shall not be used, duplicated, modified or transmitted without their authorization

Memorandum of Understanding ESTEC 4000 100 293/10/NL/PA

Wiring Panel Module 'A' Column 1					Wiring Panel Module 'A' Column 2					Wiring Panel Module 'A' Column 3					
1	RED ZS_4100_01 CATHODEs of LEDs	Door A Contact Power	+12VDC	WHITE +12VDC	1	RED AT_4108_01	EC Probe PWR	+24VDC	YELLOW See EC diagram	1					1
2	WHITE ZS_4100_01 NC WHITE ZS_4100_02 NC GREEN	Door A Contacts	Argus Input VIA Diode	YELLOW TC1-1-IN3	2	WHITE AT_4108_01	EC Probe Signal	Argus Via Signal Conditioner	GREEN See EC diagram TC1-4-IN4	2					2
3	ZS_4100_01 COM GREEN ZS_4100_02 NO	Door A Contacts	JUMPER	N/C	3	BLACK / SHIELD AT_4108_01	EC Probe Shield	+24VDC COM	Panel 'A' Column 2 Row 10	3					3
4	BLACK ZS_4100_02 NC	Door A Contact COMMON	+12VDC COM	BROWN TC1-1-IN3 COM	4					4					4
5					5	RED AT_4107_01	pH Probe PWR	+24VDC	WHITE See pH diagram	5					5
6	RED PT_4102_01	Airlock A Pressure Power	+24VDC	WHITE TC1-9-24V	6	WHITE AT_4107_01	pH Probe Signal	Argus Signal Conditioner	BROWN See pH Diagram TC1-4-IN1	6					6
7	BLACK PT_4102_01	Airlock A Pressure Signal	Argus Input 4-20mA VIA 250ohm	BROWN TC1-9-IN1	7	BLACK / SHIELD AT_4107_01	pH Probe Shield	N/C	N/C	7					7
8					8					8					8
9	BLACK PS_4102_01	Airlock A Pressure Switch SIGNAL	Argus Input	GREEN TC1-9-IN4	9	RED TT_4109_01	Reservoir Temperature	Argus Input	WHITE TC1-4-IN7	9					9
10	WHITE	N/C	N/C	N/C	10	BLACK TT_4109_01	Reservoir Temperature	Argus Input	BROWN TC1-4-IN7 COM	10					10
11	RED PS_4102_01	Airlock A Pressure Switch COMMON	Argus Input COMMON	YELLOW TC1-9-IN4 COM	11					11					11
12					12	BLACK LSL_4107_01	Acid Tank Low Level	Argus Input	GREY TC1-4-IN2	12	SV_4109_01	Nutrient Tank Chilled Water Valve	N/C	TC1-2-OP4-B	12
13					13	BROWN LSL_4107_01	Acid Tank Low Level COMMON	Argus Input COMMON	BROWN TC1-4-IN2 COM	13	SV_4109_01	Nutrient Tank Chilled Water Valve COMMON	N/C	TC1-2-OP4-B COM	13
14					14	RED LSL_4107_02	Base Tank Low Level	Argus Input	RED TC1-4-IN3	14					14
15					15	ORANGE LSL_4107_02	Base Tank Low Level COMMON	Argus Input COMMON	PINK TC1-4-IN3 COM	15	RED SV_4108_01	Stock 'A' Inject Valve	Argus Output	GREY TC1-2-OP3-B	15
16					16	YELLOW LSL_4108_01	Stock 'A' Tank Low Level	Argus Input	YELLOW TC1-4-IN5	16	BLACK SV_4108_01	Stock 'A' Inject Valve COMMON	Argus Output	BROWN TC1-2-OP3-B COM	16
17					17	GREEN LSL_4108_01	Stock 'A' Tank Low Level COMMON	Argus Input COMMON	GREEN TC1-4-IN5 COM	17	RED SV_4108_02	Stock 'B' Inject Valve	Argus Output	RED TC1-2-OP4-A	17
18	RED TT_4105_01	Loft A Temperature	Argus Input	GREY TC1-3-IN3	18	BLUE LSL_4108_02	Stock 'B' Tank Low Level	Argus Input	BLUE TC1-4-IN6	18	BLACK SV_4108_02	Stock 'B' Inject Valve COMMON	Argus Output	PINK TC1-2-OP4-A COM	18

Wiring Panel Module 'A' Column 1					Wiring Panel Module 'A' Column 2					Wiring Panel Module 'A' Column 3					
19	BLACK TT_4105_01	Loft A Temperature COMMON	Argus Input COMMON	BROWN TC1-3-IN3 COM	19	WHITE LSL_4108_02	Stock 'B' Tank Low Level COMMON	Argus Input COMMON	WHITE TC1-4-IN6 COM	19	RED SV_4107_01	Acid Inject Valve	Argus Output	YELLOW TC1-2-OP2- B	19
20	RED FSL_4105_01	Loft A Air Flow POWER	+24VDC	BROWN	20	RED LSL_4110_01	Reservoir Low Level	Argus Input	GREY TC1-5-IN1	20	BLACK SV_4107_01	Acid Inject Valve COMMON	Argus Output	GREEN TC1-2-OP2- B COM	20
21	GREEN FSL_4105_01	Loft A Air Flow GOOD SIGNAL	Argus Input	YELLOW TC1-2-IN7	21	BLACK LSL_4110_01	Reservoir Low Level COMMON	Argus Input COMMON	BROWN TC1-5-IN1 COM	21	RED SV_4107_02	Base Inject Valve	Argus Output	BLUE TC1-2-OP3- A	21
22					22	WHITE LSH_4110_01	Reservoir High Level	Argus Input	RED TC1-5-IN2	22	BLACK SV_4107_02	Base Inject Valve COMMON	Argus Output	WHITE TC1-2-OP3- A COM	22
23	WHITE FSL_4105_01	Loft A Air Flow FAN ENABLE	Argus Output	GREEN TC1-1-OP4-B	23	GREEN LSH_4110_01	Reservoir High Level COMMON	Argus Input COMMON	PINK TC1-5-IN2 COM	23					23
24	BLACK FSL_4105_01	Loft A Air Flow COMMON	24VDC Power COMMON	WHITE	24	RED LSH_4110_02	Condensate Tank Low Level	Argus Input	BLUE TC1-5-IN4	24					24
25	A1 BLACK TT_4112_04	Module A Temperature 02	Argus Input	RED TC1-6-IN1	25	BLACK LSH_4110_02	Condensate Tank Low Level COMMON	Argus Input COMMON	WHITE TC1-5-IN4 COM	25					25
26	A1 BROWN TT_4112_04	Module A Temperature 02 COMMON	Argus Input COMMON	PINK TC1-6-IN1 COM	26	WHITE LSL_4110_02	Condensate Tank High Level	Argus Input	YELLOW TC1-5-IN3	26					26
27	A1 RED TT_4112_05	Module A Temperature 03	Argus Input	YELLOW TC1-6-IN4	27	GREEN LSL_4110_02	Condensate Tank High Level COMMON	Argus Input COMMON	GREEN TC1-5-IN3 COM	27					27
28	A1 ORANGE TT_4112_05	Module A Temperature 03 COMMON	Argus Input COMMON	GREEN TC1-6-IN4 COM	28					28	RED GP_4110_01	Condensate PUMP	Argus Output	WHITE TC1-3OP1-A	28
29	A1 YELLOW TT_4112_06	Module A Temperature 04	Argus Input	BLUE N/C	29					29	BLACK GP_4110_01	Condensate PUMP	Argus Output	BROWN TC1-3OP1-A COM	29
30	A1 GREEN TT_4112_06	Module A Temperature 04 COMMON	Argus Input COMMON	WHITE N/C	30					30					30
31	A1 BLUE				31					31					31
32	A1 WHITE				32					32					32
33	A2 BLACK RT_4104_01	Module A PAR	Argus Input	WHITE TC1-2-IN4	33	A3 BLACK				33	A5 BLACK				33
34	A2 BROWN w/SHIELD RT_4104_01	Module A PAR	Argus Input COMMON	BROWN +SHIELD TC1-2-IN4 COM	34	A3 BROWN				34	A5 BROWN				34
35	A2 RED				35	A3 RED				35	A5 RED				35
36	A2 ORANGE	Camera A POWER			36	A3 ORANGE				36	A5 ORANGE				36

Wiring Panel Module 'A' Column 1					Wiring Panel Module 'A' Column 2					Wiring Panel Module 'A' Column 3					
37	A2 YELLOW	Camera A SIGNAL			37	A3 YELLOW				37	A5 YELLOW				37
38	A2 GREEN	Camera A GROUND			38	A3 GREEN				38	A5 GREEN				38
39	A2 BLUE	Camera A COMMON			39	A3 BLUE				39	A5 BLUE				39
40	A2 WHITE				40	A3 WHITE				40	A5 WHITE				40
41	RED TT_4112_01	Module A T/RH 01 POWER	+24VDC	WHITE TC1-8-24V	41	A4 BLACK				41					41
42	GREEN TT_4112_01	Module A Temperature 01 SIGNAL	Argus Input 4-20mA VIA 250ohm	GREEN TC1-8-IN1	42	A4 BROWN				42					42
43	BLACK TT_4112_01	Module A Relative Humidity 01 SIGNAL	Argus Input 4-20mA VIA 250ohm	BROWN TC1-7-IN4	43	A4 RED				43					43
44					44	A4 ORANGE				44					44
45					45	A4 YELLOW				45					45
46					46	A4 GREEN				46					46
47					47	A4 BLUE				47					47
48					48	A4 WHITE				48					48
49					49					49					49
50					50					50					50

Wiring Panel Module 'C' Column 1				Wiring Panel Module 'C' Column 2				Wiring Panel Module 'C' Column 3					
RED ZS_4101_01 CATHODES of LEDS	Door C Contact Power	+12VDC	WHITE +12VDC	1				1	B0 BLACK PS_4104_01	Module B Pressure Switch SIGNAL	Argus Input	YELLOW TC1-9-IN5	
WHITE ZS_4101_01 NC	Door C Contacts	Argus Input VIA Diode	YELLOW TC1-1-IN7	2				2	B0 BROWN PS_4104_01	Module B Pressure Switch COMMON	Argus Input COMMON	GREEN TC1-9-IN5 COM	
GREEN ZS_4101_02 NC	Door C Contacts	N/C	N/C	3				3	B0 RED PS_4104_01		+12VDC	WHITE	
BLACK ZS_4101_02 NC	Door C Contact COMMON		BROWN TC1-1-IN7 COM	4				4	B0 ORANGE PS_4104_01	Module B Vent Detect	Argus Input	BROWN TC1-9-IN7	
RED PT_4103_01	Airlock C Pressure Power	+24VDC	WHITE TC1-9-24V	5				5	B0 YELLOW				
BLACK PT_4103_01	Airlock C Pressure Signal	Argus Input 4-20mA VIA 250ohm	BROWN TC1-9-IN3	6	RED PT_4114_01	Chamber Pressure Power	+24VDC	WHITE TC1-9-24V	6	B0 GREEN			
				7	BLACK PT_4114_01	Chamber Pressure Signal	Argus Input 4-20mA VIA 250ohm	BROWN TC1-9-IN2	7	B0 BLUE			
				8					8	B0 WHITE			
RED PS_4103_01	Airlock C Pressure Switch SIGNAL	Argus Input	YELLOW TC1-9-IN5	9					9				
BLACK PS_4103_01	Airlock C Pressure Switch COMMON	Argus Input COMMON	GREEN TC1-9-IN5 COM	10	BLACK TT_4112_15	Hot Coil Temperature 1	Argus Input	GREY TC1-7-IN3	10	MVFD_4111_01	Main Blower Speed	Argus Output 4-20 mA	TC1-5-OP1 4-20mA
				11	BROWN TT_4112_15	Hot Coil Temperature 1 COMMON	Argus Input COMMON	BROWN TC1-7-IN3 COM	11	MVFD_4111_01	Main Blower Signal Return	Argus Output Return	TC1-5-OP1 COM
				12	RED TT_4112_19	Hot Coil Temperature 2	Argus Input	PINK TC1-7-IN7	12	MVFD_4111_01	Main Blower Start / Stop	Argus Output Contact	
				13	ORANGE TT_4112_19	Hot Coil Temperature 2 COMMON	Argus Input COMMON	RED TC1-7-IN7 COM	13	MVFD_4111_01	Main Blower Start / Stop	Argus Output Return	
				14	YELLOW TT_4112_16	Cold Coil Temperature 1	Argus Input	YELLOW TC1-7-IN2	14				
WHITE FSL_4105_03	Loft C Air Flow FAN ENABLE	Argus Output	GREEN TC1-2-OP1- B	15	GREEN TT_4112_16	Cold Coil Temperature 1 COMMON	Argus Input COMMON	GREEN TC1-7-IN2 COM	15				
BLACK FSL_4105_03	Loft C Air Flow COMMON	24VDC Power COMMON	WHITE	16	BLUE TT_4112_20	Cold Coil Temperature 2	Argus Input	BLUE TC1-8-IN7	16				
RED TT_4105_03	Loft C Temperature	Argus Input	GREY TC1-3-IN5	17	WHITE TT_4112_20	Cold Coil Temperature 2 COMMON	Argus Input COMMON	WHITE TC1-8-IN7 COM	17				
BLACK TT_4105_03	Loft C Temperature COMMON	Argus Input COMMON	BROWN TC1-3-IN5 COM	18					18				

Wiring Panel Module 'C' Column 1				Wiring Panel Module 'C' Column 2				Wiring Panel Module 'C' Column 3					
RED FSL_4105_03	Loft C Air Flow POWER	+24VDC	BROWN	19	RED TT_4105_02	Loft B Temperature	Argus Input	BROWN TC1-3-IN4	19				
BLACK FSL_4105_03	Loft C Air Flow GOOD SIGNAL	Argus Input	YELLOW TC1-3-IN1	20	BLACK TT_4105_02	Loft B Temperature COMMON	Argus Input COMMON	WHITE TC1-3-IN4 COM	20	RED S3CV_4112_02	Hot Water Valve Power	+24VDC	WHITE
				21					21	WHITE S3CV_4112_02	Hot Water Valve SIGNAL	Argus Output 0-10 VDC	TC1-4-OP2 0-10V
C1 BLACK/SIELD TT_4112_10	Module C Temperature 02	Argus Input	RED TC1-6-IN3	22	WHITE FSL_4105_02	Loft B Air Flow FAN ENABLE	Argus Output	GREEN TC1-2-OP1- A	22	BLACK S3CV_4112_02	Hot Water Valve COMMON	Argus Output 0-10 VDC RETURN	TC1-4-OP2 0-10V COM
C1 BROWN TT_4112_10	Module C Temperature 02 COMMON	Argus Input COMMON	PINK TC1-6-IN3 COM	23	BLACK FSL_4105_02	Loft B Air Flow COMMON	24VDC Power COMMON	WHITE TC1-2-OP1- A COM	23	RED S3CV_4112_01	Chilled Water Valve Power	+24VDC	WHITE
C1 RED TT_4112_11	Module C Temperature 03	Argus Input	YELLOW TC1-6-IN6	24					24	WHITE S3CV_4112_01	Chilled Water Valve SIGNAL	Argus Output 0-10 VDC	TC1-4-OP1 0-10V
C1 ORANGE TT_4112_11	Module C Temperature 03 COMMON	Argus Input COMMON	GREEN TC1-6-IN6 COM	25	RED FSL_4105_02	Loft B Air Flow POWER	+24VDC	BROWN	25	BLACK S3CV_4112_01	Chilled Water Valve COMMON	Argus Output 0-10 VDC RETURN	TC1-4-OP1 0-10V COM
C1 YELLOW TT_4112_12	Module C Temperature 04	Argus Input	BLUE	26	BLACK FSL_4105_02	Loft B Air Flow GOOD SIGNAL	Argus Input	YELLOW TC1-3-IN1	26				
C1 GREEN TT_4112_12	Module C Temperature 04 COMMON	Argus Input COMMON	WHITE	27	B1 BLACK TT_4112_07	Module B Temperature 02	Argus Input	GREY TC1-6-IN2	27	RED FT_4111_01	Air Flow POWER	+24VDC	WHITE
C1 BLUE TT_4112_12				28	B1 BROWN TT_4112_07	Module B Temperature 02 COMMON	Argus Input COMMON	BROWN TC1-6-IN2 COM	28	WHITE FT_4111_01	Air Flow SIGNAL	Argus Input Via 250ohm Resistor	GREEN TC1-5-IN5
C1 WHITE				29	B1 RED TT_4112_08	Module B Temperature 03	Argus Input	RED TC1-6-IN5	29	GREEN FT_4111_01	Air Flow COMMON	24VDC Power COMMON	BROWN
C2 BLACK RT_4104_03	Module C PAR	Argus Input	WHITE TC1-2-IN4 COM	30	B1 ORANGE TT_4112_08	Module B Temperature 03 COMMON	Argus Input COMMON	PINK TC1-6-IN5 COM	30	BLACK FT_4111_01			
C2 BROWN w/SIELD RT_4104_03	Module C PAR	Argus Input COMMON	BROWN w/SIELD TC1-32-IN4	31	B1 YELLOW TT_4112_09	Module B Temperature 04	Argus Input	YELLOW	31	RED TT_4111_13	Chilled Water Entry Temperature	Argus Input	BLUE TC1-6-IN7
C2 RED				32	B1 GREEN TT_4112_09	Module B Temperature 04 COMMON	Argus Input COMMON	GREEN	32	BLACK TT_4111_13	Chilled Water Entry Temperature COMMON	Argus Input COMMON	WHITE TC1-6-IN7 COM
C2 ORANGE	Camera C POWER			33	B1 BLUE				33	RED TT_4111_14	Hot Water Entry Temperature	Argus Input	GREY TC1-7-IN1
C2 YELLOW	Camera C SIGNAL			34	B1 WHITE				34	BLACK TT_4111_14	Hot Water Entry Temperature COMMON	Argus Input COMMON	BROWN TC1-7-IN1 COM
C2 GREEN	Camera C GROUND			35	B2 BLACK RT_4104_02	Module B PAR	Argus Input	WHITE TC1-2-IN5	35	WHITE TT_4111_22	Hot Water Loop Temperature	Argus Input	YELLOW TC1-5-IN6
C2 BLUE	Camera C COMMON			36	B2 BROWN /SIELD RT_4104_02	Module B PAR	Argus Input COMMON	BROWN TC1-2-IN5 COM	36	BROWN TT_4111_22	Hot Water Loop Temperature COMMON	Argus Input COMMON	GREEN TC1-5-IN6 COM

Wiring Panel Module 'C' Column 1				Wiring Panel Module 'C' Column 2				Wiring Panel Module 'C' Column 3					
C2 WHITE				37	B2 RED				37	RED TT_4111_18	Hot Water Exit Temperature	Argus Input	RED TC1-3-IN6
				38	B2 ORANGE				38	BLACK TT_4111_18	Hot Water Exit Temperature COMMON	Argus Input COMMON	PINK TC1-3-IN6 COM
				39	B2 YELLOW				39	BROWN TT_4111_21	Chilled Water Loop Temperature	Argus Input	BROWN TC1-5-IN7
				40	B2 GREEN				40	WHITE TT_4111_21	Chilled Water Loop Temperature COMMON	Argus Input COMMON	WHITE TC1-5-IN7 COM
				41	B2 BLUE				41	RED TT_4111_17	Chilled Water Exit Temperature	Argus Input	YELLOW TC1-3-IN7
RED TT_4112_03	Module C T/RH 01 POWER	+24VDC	WHITE TC1-7-24V	42	B2 WHITE				42	BLACK TT_4111_17	Chilled Water Exit Temperature COMMON	Argus Input COMMON	GREEN TC1-3-IN7 COM
GREEN TT_4112_03	Module C Temperature 01 SIGNAL	Argus Input 4-20mA VIA 250ohm	BROWN TC1-8-IN3	43	C6 BLACK	Nutrient Flow Meter COMMON	+24VDC COMMON	BROWN	43				
BLACK TT_4112_03	Module C Relative Humidity 01 SIGNAL	Argus Input 4-20mA VIA 250ohm	GREEN TC1-7-IN6	44	C6 BROWN	Nutrient Flow Meter POWER	+24VDC	WHITE	44				
				45	C6 RED	Nutrient Flow Meter SIGNAL	Argus Input	YELLOW TC1-3-IN6	45				
RED TT_4112_02	Module B T/RH 01 POWER	+24VDC	WHITE TC1-7-24V	46	C6 ORANGE				46				
BLACK TT_4112_02	Module B Temperature 01 SIGNAL	Argus Input 4-20mA VIA 250ohm	BROWN TC1-8-IN5	47	C6 YELLOW				47				
WHITE TT_4112_02	Module B Relative Humidity 01 SIGNAL	Argus Input 4-20mA VIA 250ohm	GREEN TC1-7-IN2	48	C6 GREEN				48				
SHIELD				49	C6 BLUE				49				
				50	C6 WHITE				50				

# MELISSA

**TECHNICAL NOTE**

**96.1**

## ***APPENDIX 8***


This document is confidential property of the MELISSA partners and shall not be used, duplicated, modified or transmitted without their authorization

Memorandum of Understanding ESTEC 4000 100 293/10/NL/PA


# UAB Argus Panel Wiring

2009-02-06


TC1-1


OP1A


OP1B


OP2A


OP2B


OP3A


OP3B

OP4A



OP4B


# UAB Argus Panel Wiring

2009-02-06

## TC1-2

24VDC  
 Common COM  
 with  
 TC1-3





# UAB Argus Panel Wiring

2009-02-06

## TC1-3

CHA +	
CHA -	
COM	
CHB +	
CHB -	
COM	
IN1	Loft 'B' Fan Status
COM	
IN2	Loft 'C' Fan Status
COM	
IN3	Loft 'A' Temperature
COM	
IN4	Loft 'B' Temperature
COM	
IN5	Loft 'C' Temperature
COM	
IN6	Hot Water Exit Temperature
COM	
IN7	Cold Water Exit Temperature
COM	
24VDC	
Common	COM






with  
TC1-2


# UAB Argus Panel Wiring

2009-02-06



## TC1-4

## OP3A

## OP3B


# UAB Argus Panel Wiring

2009-02-06

## TC1-5

# UAB Argus Panel Wiring

2009-02-06

TC1-6		
<hr/>		
CHA +		
CHA -		
COM		
<hr/>		
CHB +		
CHB -		
COM		
<hr/>		
IN1	A2	
COM	Temperature	
<hr/>		
IN2	B2	
COM	Temperature	
<hr/>		
IN3	C2	
COM	Temperature	
<hr/>		
IN4	A3	
COM	Temperature	
<hr/>		
IN5	B3	
COM	Temperature	
<hr/>		
IN6	C3	
COM	Temperature	
<hr/>		
IN7	Chilled Water Source	
COM	Temperature	
<hr/>		
24VDC		
COM		
<hr/>		

OP1A

OP1B

OP2A

OP2B

OP3A


OP3B

OP4A

OP4B


# UAB Argus Panel Wiring

2009-02-06

# UAB Argus Panel Wiring

2009-02-06

Tie MFC 24VDC common to COM

OP4A


OP4B

# UAB Argus Panel Wiring

2009-02-06

TC1-9

OP1A


OP4A

OP4B

# UAB Argus Panel Auxiliary Power Wiring

2009-02-06

# MELISSA

**TECHNICAL NOTE**


**96.1**

## ***APPENDIX 9***

This document is confidential property of the MELISSA partners and shall not be used, duplicated, modified or transmitted without their authorization

Memorandum of Understanding ESTEC 4000 100 293/10/NL/PA

# HPC1 Airlock and CO2 Gas Control

# MELISSA

**TECHNICAL NOTE**

**96.1**

## ***APPENDIX 10***

This document is confidential property of the MELISSA partners and shall not be used, duplicated, modified or transmitted without their authorization

Memorandum of Understanding ESTEC 4000 100 293/10/NL/PA

**Angstrom Engineering Inc.  
91 Trillium Drive  
Kitchener, Ontario  
N2E 1W8**

**Attention: Dave Pitts**

**INVESTIGATION OF CORROSION OF A 316SS GROWTH  
CHAMBER**

**PO# Master Card  
Laboratory Ref. #G808171  
October 21, 2008**

**BODYCOTE TESTING GROUP**

Per:  **Bob Elliott, MSc  
Corrosion Scientist, Metallurgical Engineering**

Per:  **Carl Fleck, CET  
Metallurgical Department Manager**

Samples will be retained for a period of three months prior to disposal.

## **1.0 INTRODUCTION**

Angstrom Engineering Inc. requested Bodycote Testing Group to examine a grow chamber located at the University of Guelph to determine the cause of the corrosion. The grow chamber was made from 316L SS and was in service for only six months. The chamber components had been glass beaded and were washed with detergent (Simple Green Concentrate or Alkonox) before use. The chamber normally operated at about 75% relative humidity; however, the air recirculation fan failed at some point and the humidity increased to about 100%. The corrosion appeared during the last twenty days. It was initially reported that the chamber had been disinfected by circulation of a 0.5% bleach solution to remove algae prior to the onset of the corrosion; this was subsequently clarified to the bleach solution having been used following appearance of the corrosion in an attempt to remove it. It was also reported that the chamber was not rinsed after application of the bleach solution. A nutrient solution is dripped onto the trays containing the plants; percolates through the rock wool used to support the plants and is collected and recirculated.

## **2.0 INVESTIGATION AND RESULTS**

R. W. Elliott visited the facility on 17 Sep 08. Figures 1 and 2 are photographs illustrating the condition of the air recirculation chamber floor and wall. Samples were taken of the rust and white coloured deposits on the floor and the rust coloured deposits on the wall. Figure 3 illustrates the nodular corrosion products that were liquid filled.

The open circuit potential (OCP) of the frame was measured in a clean and in a corroded area using a Saturated Calomel Reference Electrode (SCE) with a cotton ball soaked with 5% sodium chloride solution. The approximate location is indicated in Figure 1. The open circuit potential of the clean area was -0.250 to -0.300 Volts versus and that of the rusted area was -0.175 to -0.200 V vrs SCE. The stainless steel was in the passive condition.

Figure 4 shows the interior of the grow chamber and a tray link that was provided as a sample. The tray link shows the same type of liquid filled nodules as the wall of the air recirculation chamber.

Both of the air locks at the ends of the grow chamber had been left closed during the last grow cycle. One had been cleaned and the other which had not been cleaned had numerous white coloured deposits that could not be scrapped off for analysis. Figure 5 shows a photograph of these including an area that had been cleaned using an abrasive pad. Pitting had started to initiate in the cleaned area.

The panel from the air recirculation chamber also showed similar liquid filled nodules (Figure 6) as those on the tray link and the wall of the air recirculation chamber. The areas opposite welds were more affected than the bulk metal.

The samples of rust coloured deposits from the floor and walls of the air recirculation chamber and on the tray link and the white deposit from the floor of the air recirculation chamber were analyzed using a JEOL JSM 5600 scanning electron microscope (SEM) equipped with an energy dispersive x-ray spectrometer (EDS) MII#B05028. Figures 7 and 8 show EDS scans taken of the deposits. The elements detected are listed in Table 1.

Table 1  
Elements Detected By EDS

Sample	Major	Minor	Trace
Floor – rust coloured	Chlorine [Cl], iron [Fe], chromium [Cr]	Nickel [Ni], sodium [Na]	Sulphur [S], carbon [C], oxygen [O], silicon [Si], calcium [Ca], phosphorus [P]
Floor – white	S, Na, K	Magnesium [Mg], P, Cl, Ca O	Si
Wall – rust coloured	Cl, Fe, Cr	Ni, S	Na, Ca, Si
Tray clip - nodules	Fe, Cl, Cr	Si, Ni	C, O, Na, P, K, Ca, Aluminum [Al]

A sample of the deposit from the wall of the recirculation chamber was added to a test solution specific for sulphate reducing bacteria. The test was negative. The tray clip was examined for the presence of iron oxidizing bacteria. The test was negative.

### 3.0 DISCUSSION

The appearance of the corrosion products on the walls and panel of the air recirculation chamber and on the tray clip suggested microbiologically induced corrosion (MIC). The nodules on the tray link were examined for the presence of iron oxidizing bacteria. These bacteria can cause pitting corrosion of stainless steel and are known to be in certain municipal water supplies in the general area of Guelph, Kitchener, Waterloo, and Cambridge. Tests for sulphate reducing bacteria and for iron oxidizing bacteria were both negative.

Stainless steels owe their corrosion resistance to a passive film of chromium oxide that forms on the surface. The chloride ion has the ability to penetrate this passive film and create a pitting site. The shallow pits that had started to form under the nodules on the tray link were round and slightly undercut, characteristic of chloride pits.

The nutrient solution composition provided is listed in Table 2.

Table 2  
 Nutrient Solution Composition

Component	Molecular Weight	Feed Strength (mM)
$\text{Ca}(\text{NO}_3)_2 \cdot 4\text{H}_2\text{O}$	236.16	3.62
$\text{MgSO}_4 \cdot 7\text{H}_2\text{O}$	246.48	1
$\text{KNO}_3$	101.1	5
$\text{NH}_4\text{H}_2\text{PO}_4$	115.08	1.5
$(\text{NH}_4)_2\text{SO}_4$	132	1
$\text{FeCl}_3$ (DTPA)	162.2	0.025
$\text{H}_3\text{BO}_4$	61.83	0.02
$\text{MnSO}_4 \cdot \text{H}_2\text{O}$	169.01	0.005
$\text{ZnSO}_4 \cdot 7\text{H}_2\text{O}$	289.54	0.0035
$\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$	249.68	0.0008
$\text{H}_2\text{MoO}_4$ (85% $\text{MoO}_3$ )	161.97	0.0005

The nutrient solution was considered as a possible source of some of the elements detected by the EDS scans. The concentrations in the nutrient solution of elements detected by EDS were potassium (1.9 ppm), calcium (0.6 ppm), phosphorus (0.4 ppm), sulphur (0.4 ppm), magnesium (0.1 ppm) and chlorine (0.02 ppm). These concentrations are quite low. For the nutrient solution

to be a source of these elements, it would have to be entrained in the recirculating air and deposited on the walls of the chamber where it could concentrate by evaporation. The very low chlorine concentration in the nutrient solution and the prominence of the chlorine peak in the EDS scans suggest that the chlorine is from a source other than the nutrient solution.

Bleach (sodium hypochlorite, NaClO) was used at a 0.5% (500 ppm) concentration reportedly only after the appearance of the corrosion. It was reported that it was not rinsed off after application. Sodium, as well as chlorine was detected in all of the EDS scans. Reference [1] offers the following advice for the use of bleach to ensure the best sanitizing conditions while minimizing stainless steel pitting. The recommended sodium hypochlorite concentration is 200 ppm, the contact time no longer than 30 minutes and the temperature no higher than 49°C. Some of the corrosion observed during the site visit and the sodium and chlorine detected by the EDS may be the result of the circulation of the 500 ppm bleach solution with no subsequent rinsing.

It has been reported that the grow chamber has been disassembled and passivated at a commercial established specialising in this operation. It was further reported that there were some areas in which a rust coloured ring remained around pitting sites.

It has also been reported that pickling paste (gel containing nitric and hydrofluoric acids) was used to clean rusted areas but that the rust reappeared. It is important to ensure that pickling pastes are thoroughly rinsed from the surface of the stainless steel and that it be allowed to dry to form the protective passive film.

The interior of the chamber can be cleaned using an abrasive pad such as Scotchbrite. The surface should be rinsed and left to dry in air for at least 24 hours to re-establish the passive film. Alternatively, the passive film can be regenerated by application of a solution of 1% citric acid and 1% sodium nitrate at room temperature with a contact time of one hour followed by rinsing and drying. This cleaning should be performed whenever similar rusting appears.

It is suggested that the air flow velocity in the chamber be reviewed to determine if it is entraining nutrient solution when in operation.

It is also suggested that the bead blasting operation be reviewed to ensure that only clean glass beads have been used and that the resulting surface is free from laps (folds in the metal) that could trap solutions.

The MSDS for the two detergents, Simple Green Concentrate and Alkonox, used during construction of the grow chamber were reviewed and neither reported chlorine containing compounds among the ingredients. However, it may be advisable to analyse these materials to determine if they are a source of the chloride ion.

#### **4.0 CONCLUSION**


The corrosion of the interior surfaces of the grow chamber appear to result from deposition of chloride ions on the surfaces. The source is not clear as the 0.5% bleach solution was reported not to have been applied until after the appearance of the rust. This solution was also reported not to have been rinsed off after application; therefore, the chlorine detected in the corrosion products may have been deposited after the corrosion occurred. The appearance of the pits under the rust nodules support chloride corrosion and the appearance and nature of the nodules support microbiologically induced corrosion; however tests for sulphate reducing bacteria and for iron oxidizing bacteria were both negative.

#### **5.0 REFERENCES**

- [1] ASM Handbook Volume 13C, "Corrosion: Environments and Industries", ASM International, Materials Park, Ohio, 2006, page 808



**Figure 1:** The lower air recirculation chamber floor. Samples were taken of the rust coloured deposits on the floor and wall and of the white coloured deposits on the floor. The open circuit potential was measured in the corroded area on top of the frame and the clean area on the front of the frame.


**Figure 2:** The lower air recirculation chamber wall.


**Figure 3: Close-up photographs of the corrosion on the air recirculation chamber wall. Note the nodules (red arrows in the upper photograph). These were filled with liquid, which smudged when sampled (lower photograph).**


**Figure 4:** The upper photograph shows the interior of the growth chamber with a tray link indicated by the red arrow. The lower photograph is a tray link showing the same type of corrosion as apparent on the wall of the air recirculation chamber. This tray link was taken for further examination.

**Figure 5: Two photographs of the wall in the air lock that had not been cleaned. The upper photograph illustrates a large number of white coloured deposits. These were too adherent to sample for further examination. The lower photograph shows an area of the wall that had been cleaned. There were faint signs of pitting under the white deposits.**

**Figure 6: Two photographs of the panel from the air recirculation chamber. Similar nodules to those found on the walls of the chamber and on the tray link were apparent and were more numerous opposite welds.**


Figure 7: EDS scans taken of the rust deposits from the floor (top) and the white deposits from the floor (bottom).


Figure 8: EDS scans taken of the deposit sample taken from the wall (top) and a crusted deposit on the tray clip (bottom).

# MELISSA

**TECHNICAL NOTE**

**96.1**

## ***APPENDIX 11***



This document is confidential property of the MELISSA partners and shall not be used, duplicated, modified or transmitted without their authorization

Memorandum of Understanding ESTEC 4000 100 293/10/NL/PA

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


**96.1**

## ***APPENDIX 12***

This document is confidential property of the MELISSA partners and shall not be used, duplicated, modified or transmitted without their authorization

Memorandum of Understanding ESTEC 4000 100 293/10/NL/PA

# 600

SERIES

# NDIR/O<sub>2</sub>

## Infrared/Oxygen Multi-Component Analyzer

### APPLICATIONS

- Stack Gases (CEM)
- Combustion Efficiency
- Turbine/Generator Feedback Control
- Process Chemical Gas Analysis
- Personnel Safety
- Fuel Cell Analysis
- Vehicle Emissions

### OPTIONS

- Internal Sample Pump
- Sample Flow Control
- Multiple Inputs
- 19 Inch Rack Mount Slides

### FEATURES

- Measures From Low ppm up to 100% Full Scale
- Measures Oxygen From 0-1% up to 0-100%
- Multi-Component—Up to Three IR Channels or Two IR Channels Plus Oxygen
- Up to Four User-Definable Ranges Per Channel (Within Factory Pre-Set Ranges)
- CE Mark and ETL Listed—Conforms to UL STD 61010-1, Certified to CAN/CSA C22.2 STD No. 610610.1
- Auto Calibration and Ranging
- Fast Response Time
- Temperature and Pressure Compensation
- Comprehensive Diagnostics/Data Archiving
- Output Options: Voltage, Current, RS-232, TCP/IP
- Remote Monitoring and Control

1312 West Grove Avenue  
Orange, CA 92865-4134  
Phone: 714-974-5560 Fax: 714-921-2531  
[www.gasanalyzers.com](http://www.gasanalyzers.com)

# 600 Series NDIR/O<sub>2</sub> Analyzer

# NDIR/O<sub>2</sub>

## DESCRIPTION

The California Analytical Models 601-2-3 NDIR/O<sub>2</sub> digital analyzer product line is designed around a state of the art 16 bit microprocessor, with 16 digital inputs, 16 digital outputs, 16 analog inputs and 4 analog outputs. The analyzer can be manually operated from the keypad or remotely via TCP/IP, RS-232C communications and discrete inputs. The analyzer display includes screen presentation of all analyzer alarms. Four levels of password protection are provided. For precision measurements, the analyzer's accuracy is increased by entering calibration curve fit polynomials. Automatic calibration may be activated locally or remotely and includes auto calibration via preset times.

## METHOD OF OPERATION—NDIR

The California Analytical Instruments' NDIR analyzer is based on the infrared absorption characteristics of gases. Using a single infrared beam to measure gas concentrations, this analyzer produces highly stable and reliable results. A single infrared light beam is modulated by a chopper system and passed through a sample cell of predetermined length containing the gas sample to be analyzed. As the beam passes through the cell, the sample gas absorbs some of its energy. The attenuated beam (transmittance) emerges from the cell and is introduced to the front chamber of a two-chamber infrared microflow detector. The detector is filled with the gas component of interest and consequently the beam experiences further energy absorption. This absorption process increases the pressure in both of the chambers. The differential pressure between the front and rear chambers of the detector causes a slight gas flow between the two chambers. This flow is detected by a mass-flow sensor and is converted into an output signal.

## METHOD OF OPERATION - Oxygen

The California Analytical Instruments oxygen analyzer section utilizes either the paramagnetic or fuel cell method to determine the percent level of oxygen contained in the sample gas. The oxygen level is displayed on the LCD panel in percent concentration.

## SPECIFICATIONS

**IR Analysis Method:** Non-Dispersive Infrared (NDIR)

**NDIR Components:** CO / CO<sub>2</sub> / CH<sub>4</sub> / C<sub>3</sub>H<sub>8</sub> / SO<sub>2</sub>

**Detector Type:** Microflow

**NDIR Ranges:** From 0-50\* ppm up to 0-100%

\* CH<sub>4</sub>, C<sub>3</sub>H<sub>8</sub> 0-250 ppm minimum, SO<sub>2</sub> 1000 ppm minimum

**Range Ratio:** 10:1

**Response Time (IR):** 90% of Full Scale in < 1 Second\*\*

\*\*Depending on Cell Length, Flow Rate, and Time Constant

**IR Sample Cell:** Stainless Steel w/ Replaceable gold cell liner

**Resolution:** Typically 0.1% of Full Scale

**Repeatability:** Better than 1.0% of Full Scale

**Linearity:** Better than 1.0% of Full Scale

**Noise:** Less than 1% of Full Scale

**Zero & Span Drift:** Less than 1% of Full Scale per 24 Hours

**Zero & Span Adjustment:** Via front panel, TCP/IP or RS-232

**Sample Flow Rate:** 0.5 to 2.0 LPM

**Oxygen Analysis Method:** Paramagnetic or Fuel Cell

**O<sub>2</sub> Ranges:** 0-1% (Paramagnetic Only) up to 0-100% O<sub>2</sub> Full Scale, Four Definable Ranges

**O<sub>2</sub> Response Time:** T90—2 Seconds Paramagnetic, 16 Seconds Fuel Cell

### 600 Series Features:

**Outputs available:** TCP/IP, RS232 0-1, 0-5, 0-10 VDC, 4-20mA (selectable)

**Discrete Control:** Remote/Local Control, Range Change, Range Sense Mode

**Discrete Alarms/Control:** 15 definable optically isolated solid state relays (60 VDC max @ 600 mA max)

**Digital Diagnostics:** Pressure, Pressure Control Voltages, Temperatures, Flow Parameters

**Keypad Displays:** Factory Settings, TCP/IP address, Passwords(4), Scalable Analog Output Voltages, Full Scale Range Select, Auto Cal Times

**Special Features:** Auto Ranging, Auto Calibration (adjustable through internal clock), Data Archiving

**Display:** 3" x 5" Back lit LCD

**Sample Temperature:** Up to 50° C, Non-condensing

**Ambient Temperature:** 5° to 40° C

**Ambient Humidity:** Less than 90% RH (Non-condensing)

**Fittings:** ¼ inch Tube

**Power Requirements:** 115 VAC/60 Hz or 230 VAC/50 Hz

**Dimensions:** 5¼"Hx19"Wx23"D

**Weight:** 30-50 lbs. (Depending on configuration)

Specifications subject to change without notice.

1312 West Grove Avenue  
Orange, CA 92865-4134  
Phone: 714-974-5560 Fax: 714-921-2531  
www.gasanalyzers.com

# 600

SERIES

# NDIR

Infrared/Oxygen Multi-Component Analyzer

# USER'S

# MANUAL

**The Model 600 NDIR Series Instruments starting with Serial Number U06081 have several new Hardware and Software features. For a complete explanation, see Appendix 2 starting on page 89.**

**NOTE:** Instrument Serial Number 7013 is an exception and does not contain the new features described above

**Safety Alert**  
Caution or Warning

**Temperature Hazard**  
Caution or Warning

**Electrical Shock Hazard**  
Caution or Warning

## Safety Information in this Manual

Note, caution and warning symbols appear on the instrument and throughout this manual to draw your attention to important operational and safety information.

A “**NOTE**” marks a short message to alert you to an important detail.

A “**CAUTION**” safety alert appears with information that is important for protecting your equipment and performance.

A “**WARNING**” safety alert appears with information that is important for protecting you, others and equipment from damage. Pay very close attention to all warnings that apply to your application.

The  symbol (an exclamation point in a triangle) precedes a general CAUTION or WARNING statement.

The  symbol (wavy vertical lines with an underline in a triangle) precedes an elevated temperature hazard CAUTION or WARNING statement.

The  symbol (a lightning bolt in a triangle) precedes an electric shock hazard CAUTION or WARNING statement.

Some or all of the above symbols may appear in this manual or on the equipment. This manual should be consulted whenever one of these symbols is encountered on the equipment.

**ALWAYS REMOVE POWER BEFORE CONNECTING OR DISCONNECTING SIGNAL CABLES OR WHEN SERVICING THE EQUIPMENT.**

**The 600 series NDIR instruments meet or exceed the following directives and standards.**

Application of Council Directive(s):

*Electrical Safety:*

Low Voltage Directive 73/23/EEC

*Electromagnetic Compatibility:*

EMC Directive 89/336/EEC

Standard(s) to which Conformity is Declared:

*Electrical Safety:*

*Standard for Electrical Equipment for Measurement, Control, and Laboratory Use [EN 61010-1:2001 (2nd Edition)*

*Electromagnetic Compatibility:*

*EN 61326:1997 Electrical equipment for measurement, control and laboratory use - EMC requirements (Amendment A1: 1998 to EN 61326:1997; Amendment A2:2001 to EN 61326:1997)*

---

**600 NDIR Quick Start Guide**

- 1) Plug in the analyzer and turn the power on.
- 2) Connect the appropriate gas lines and vents to the analyzer.
- 3) Allow the analyzer to stabilize for at least one (1) hour.
- 4) During the analyzer's stabilization period, setup the analyzer to the desired configuration.
  - a) Set the analyzer's output as desired.
    - i) From the Main Menu press F5 (Setup) then F7 (System Settings.)
    - ii) Press F3 (output assignments) to setup the output channels.
      - (1) Set to real time.
      - (2) AUX 1
      - (3) AUX 2
      - (4) AUX 3
      - (5) Press back to return to the system settings menu.
    - iii) Press F4 (output range) to setup the output ranges
      - (1) Set the Min and Max to zero (0). This will cause the outputs to default to the current ranges.
      - (2) Set the outputs for milliamps (mA) or voltage (V) as desired.
        - (a) The mA setting will give 4-20 mA or 2-10V (if 500 ohm resistors are installed).
        - (b) The V setting will give 0-20mA or 0-10V (if 500 ohm resistors are installed).

**Disconnect power before proceeding**

- 5) Connect all appropriate analog outputs.
  - a) Pin 1 is the output common (ground).
  - b) Pins 2, 3, 4, and 5 are output channels 1, 2, 3, and 4 as setup in step 4.a.ii.
  - c) If the output is set for voltages but there is no voltage output, you will need to install a 500 ohm resistor between pin 1 and the channel you are trying to measure. You will have to measure the voltage drop across this resistor.
- 6) Reconnect the power and turn on the analyzer. Press F1 from the analyzer's main Menu to get to the measurement screen.
- 7) Supply sample gas to the analyzer.
- 8) The measurement screen should indicate the concentration of the sample gas. If the concentration is incorrect, the analyzer will need to be calibrated using zero and span calibration gases. Please reference

the appropriate section of the manual for a description on how to zero and span the analyzer.

---

## TABLE OF CONTENTS

Low Voltage Directive 73/23/EEC .....	4
EMC Directive 89/336/EEC .....	4
<b>1. Introduction.....</b>	<b>16</b>
1.1. Overview.....	16
1.2. Unpacking Instructions .....	16
1.3. Reporting Damage .....	16
1.4. Contact Information .....	16
1.5. Warranty Certificate.....	17
1.6. Proper Operation .....	18
1.7. Possible Explosion Hazard.....	18
1.8. Electrical Shock Hazard .....	18
1.9. Plug Removal .....	18
<b>2. Features.....</b>	<b>19</b>
2.1. Description .....	19
2.2. Features-General .....	19
2.3. Features-General .....	19
2.4. Infrared Gas Analyzers.....	20
2.5. Interference Gases.....	21
2.6. Electronics .....	22
2.7. Model 600 NDIR Specifications.....	23
2.8. Model 600 NDIR Features.....	23
<b>3. Installation.....</b>	<b>24</b>
3.1. General.....	24
3.2. Site and Mounting.....	24
3.3. Electrical .....	24
3.4. Analog Output Connections (Appendix) .....	25
3.5. Gases .....	25
3.6. Gas Handling Equipment.....	26
3.7. Gas Connections .....	26
3.8. Sampling Requirements .....	26
3.8.1. Filtration .....	26
3.8.2. Condensation.....	26

---

3.8.3.	Presence of Corrosive Gases.....	26
3.8.4.	Gas Temperature.....	26
3.8.5.	Sample Gas Bypass Outlet (Vent).....	26
<b>4.</b>	<b>Basic Operation.....</b>	<b>27</b>
4.1.	Menu Trees.....	29
4.2.	Keyboard.....	32
4.2.1.	Operation with the Cursor Keys and the Enter Key.....	32
4.2.2.	Operation with the Function Keys.....	32
4.2.3.	Read/Change Parameters.....	32
<b>5.</b>	<b>Operating Structure.....</b>	<b>33</b>
5.1.	The Main Menu.....	34
<b>6.</b>	<b>Menu Structure.....</b>	<b>35</b>
6.1.	User Functions (Level 1).....	35
6.2.	Advanced User Functions (Level 2).....	35
6.3.	Maintenance Functions (Level 3).....	35
6.4.	System User Functions (Level 4).....	35
<b>7.</b>	<b>Main Menu Function Descriptions.....</b>	<b>36</b>
7.1.	F1 Measurements.....	36
7.1.1.	F1 Measurement.....	36
7.1.2.	Range Select.....	37
7.2.	F2 Purge Analyzer.....	38
7.3.	F3 Diagnostics.....	38
7.4.	F4 Calibrations.....	39
7.4.1.	F1 Automatic Calibration.....	39
7.4.2.	F2 Manual Calibration.....	40
7.4.3.	Display Deviations – from Calibration menu F3.....	41
7.4.4.	Absolute Zero Gas Deviation.....	41
7.4.5.	Relative Zero Gas Deviation.....	41
7.4.6.	Absolute Span Gas Deviation.....	42
7.4.7.	Relative Span Gas Deviation.....	42
7.4.8.	F4 Check Calibration.....	42
7.4.9.	F5 Reset Calibration Values.....	42
7.4.10.	F6 Range Select.....	42

---

7.5.	F5 Setup.....	43
7.5.1.	F1 Span Gas Concentration .....	43
7.5.2.	Calibration Settings.....	44
7.5.3.	F1 Times.....	44
7.5.4.	F2 Measuring Deviations .....	44
7.5.5.	F3 Deviations .....	44
7.5.6.	F3 Range Limits.....	44
7.5.7.	F1 Range 1-4 (Change Upper Range Limits) .....	45
7.5.8.	F2 Change Auto Range Limits.....	45
7.5.9.	F4 Alarms .....	46
7.5.10.	F5 Password.....	46
7.5.11.	F1 Enter Password .....	47
7.5.12.	F2 Change Password .....	47
7.5.13.	F3 Reset Passwords.....	47
7.5.14.	F6 Linearization .....	48
7.6.	F7 System Settings .....	49
7.6.1.	F1 Real Time Clock .....	49
7.6.2.	System Setup F2 Displays TCP/IP Address .....	51
7.6.3.	Systems Setup F3 Displays Output Signal Assignments .....	51
7.6.4.	System Setup F4 Displays Output Ranges .....	51
7.6.5.	F5 Turns Status Line On or Off.....	52
7.6.6.	F8 Measure Settings.....	52
7.6.7.	F1 Set Lowpass filter .....	52
7.6.8.	F2 Purge Time .....	53
7.6.9.	F3 Set Temperature Compensation .....	53
7.6.10.	F4 Pressure Compensation .....	53
7.6.11.	F3 Low Pass Filter Time Constant.....	54
7.6.12.	F10 Displays the Current Analyzer and Software Versions .....	54
7.7.	F7 Remote / Manual Control .....	54
7.8.	F8 Standby .....	55
<b>8.</b>	<b>Analyzer Components.....</b>	<b>56</b>
8.1.	Rear Panel.....	56
8.2.	Rear Panel Connectors .....	57
8.3.	Digital Outputs .....	58

8.3.1.	RS-232 (Standard 9 Pin DIN Connector).....	58
8.3.2.	TCP/IP (8 Pin RJ-47 Connector) .....	58
8.4.	Internal Component Locations.....	59
8.5.	Main Electronics Board (Potentiometers).....	60
8.6.	NDIR Detector Assembly.....	61
<b>9.</b>	<b>Operation.....</b>	<b>62</b>
9.1.	External Wiring .....	62
9.2.	External Piping .....	62
9.3.	Operation & Calibration .....	62
9.3.1.	Power On:.....	62
9.3.2.	Zero/San Adjustment: .....	62
9.3.3.	Start-Up and Routine Maintenance: .....	62
<b>10.</b>	<b>MAINTENANCE .....</b>	<b>63</b>
10.1.	Zero and Span Calibration.....	63
10.2.	Cleaning of the Optical Bench Measuring Cell (Infrared Analyzers Only) .....	63
10.3.	Optical Bench Configuration.....	63
10.4.	Removal of Pipe Cell .....	65
10.5.	Removal of Block Cell .....	67
10.6.	Disassembly of Combination Pipe & Block Type Cells.....	69
<b>11.</b>	<b>ADJUSTMENTS CHECKS AND REPAIRS .....</b>	<b>71</b>
11.1.	Adjustment of Detector Voltage (NDIR'S only).....	71
11.2.	Coarse Zero Adjustment.....	71
11.3.	Span Gain Adjustment.....	71
11.4.	Check and Repair Detector .....	72
11.4.1.	Problem: .....	72
11.4.2.	Symptom:.....	72
11.4.3.	Check and/or replace:.....	72
11.5.	Check and Repair Infrared Light Source Unit.....	73
11.5.1.	Problem: .....	73
11.5.2.	Symptom:.....	73
11.5.3.	Check and/or replace:.....	73
11.6.	Check and Replace Chopper .....	73
11.6.1.	Problem: .....	73
11.6.2.	Symptom:.....	73

11.6.3.	Check and/or replace:.....	73
11.7.	Check and Repair Measuring Cell.....	74
11.7.1.	Problem: .....	74
11.7.2.	Symptom:.....	74
11.7.3.	Check and or replace:.....	74
11.8.	Check and Repair Tubing Trouble.....	74
11.8.1.	Problem: .....	74
11.8.2.	Symptom.....	74
11.8.3.	Check and/or replace:.....	74
11.9.	Check and Repair Amplifier Circuit.....	75
<b>12.</b>	<b>Communication Master Computer / Analyzer (AK Protocol).....</b>	<b>76</b>
12.1.	Serial Interface and AK-Commands .....	76
12.2.	Interface Specifications .....	76
12.3.	Protocol Description .....	77
12.3.1.	Instruction command .....	77
12.3.2.	Acknowledgement command.....	77
12.3.3.	Data Description .....	78
12.3.4.	Error Handling .....	78
12.4.	Scan Commands.....	79
12.4.1.	<b>AKON:</b> Measured concentration value.....	79
12.4.2.	<b>AEMB:</b> Get measuring range .....	79
12.4.3.	<b>AMBE:</b> Measuring range limit.....	79
12.4.4.	<b>AKAK:</b> Calibration gas concentrations.....	79
12.4.5.	<b>AMBU:</b> Upper and lower range switchover values for auto range.....	79
12.4.6.	<b>ASTZ:</b> Normal device status.....	80
12.4.7.	<b>ASTF:</b> Error status.....	80
12.4.8.	<b>AKEN:</b> Device identification .....	80
12.4.9.	<b>ARMU:</b> Raw value .....	80
12.4.10.	<b>ATEM:</b> Temperatures .....	81
12.4.11.	<b>ADRU:</b> Pressures/ Valve voltage .....	81
12.4.12.	<b>ADUF:</b> Flows .....	81
12.4.13.	<b>AGRD:</b> Polynom coefficients .....	81
12.4.14.	<b>AANG:</b> Deviation from zero point after autocalibration .....	81
12.4.15.	<b>AAEG:</b> Deviation from end point after autocalibration.....	81

12.4.16.	<b>AFDA:</b> Purge and Autocalibration times.....	81
12.4.17.	<b>APAR:</b> Request Autocalibration tolerance values .....	82
12.4.18.	<b>AKAL:</b> Deviations from calibration .....	82
12.4.19.	<b>ASYZ:</b> Respond System Time.....	82
12.4.20.	<b>AT90:</b> Respond Lowpass filter time.....	82
12.4.21.	<b>ADAL:</b> Diagnostic alarm limits.....	82
12.4.22.	<b>ATCP:</b> Query TCP/IP settings .....	82
12.4.23.	<b>AVER:</b> Query Software version .....	83
12.4.24.	<b>AH2O:</b> Query H <sub>2</sub> O correction parameter .....	83
12.4.25.	<b>ACO2:</b> Query CO <sub>2</sub> correction parameter .....	83
12.4.26.	<b>AUPD:</b> Query UDP data streaming parameter .....	83
12.5.	Control commands .....	84
12.5.1.	<b>SRES:</b> Reset .....	84
12.5.2.	<b>SPAU:</b> Pause .....	84
12.5.3.	<b>STBY:</b> Standby.....	84
12.5.4.	<b>SNGA:</b> Open valve for zero gas calibration.....	84
12.5.5.	<b>SEGA:</b> Open valve for end gas calibration.....	84
12.5.6.	<b>SSPL:</b> Purge Analyzer with zero gas .....	84
12.5.7.	<b>SATK:</b> Start automatic calibration .....	84
12.5.8.	<b>SEMB:</b> Set measuring range.....	84
12.5.9.	<b>SARE:</b> Auto range on .....	85
12.5.10.	<b>SARA:</b> Auto range off.....	85
12.5.11.	<b>SREM:</b> Remote mode for AK-commands.....	85
12.5.12.	<b>SMAN:</b> Manual control to control device manually .....	85
12.5.13.	<b>SMGA:</b> Start measuring .....	85
12.5.14.	<b>SNKA:</b> Saves measured value as new offset.....	85
12.5.15.	<b>SEKA:</b> Saves measured value as new span value .....	85
12.5.16.	<b>SUDP:</b> Start /Stop UDP data streaming .....	85
12.6.	Settings .....	86
12.6.1.	<b>EKAK:</b> The four span gas concentration values are set .....	86
12.6.2.	<b>EMBE:</b> The four measuring range end values are set .....	86
12.6.3.	<b>EMBU:</b> The upper and the lower range switchover for auto range are set..	86
12.6.4.	<b>EKEN:</b> Set new device identification and information .....	86
12.6.5.	<b>EGRD:</b> Set polynom coefficients .....	86

---

12.6.6.	<b>EFDA:</b> Set autocalibration and purge times .....	86
12.6.7.	<b>EPAR:</b> Set autocalibration tolerance values.....	86
12.6.8.	<b>ESYZ:</b> Set System Time .....	87
12.6.9.	<b>ET90:</b> Set Lowpass Filter Time .....	87
12.6.10.	<b>EDAL:</b> Diagnostic alarm limits .....	87
12.6.11.	<b>ETCP:</b> Set TCP/IP Parameters.....	87
12.6.12.	<b>EH2O</b> Set H <sub>2</sub> O correction parameters .....	87
12.6.13.	<b>ECO2</b> Set CO <sub>2</sub> correction parameters .....	87
12.6.14.	<b>EUDP</b> Set UDP Data streaming parameters .....	88
12.6.15.	Format of the streaming Data via UDP: .....	88
12.7.	Abbreviations used .....	88
<b>13.</b>	<b>Appendix 1 Electrical Block Diagram .....</b>	<b>89</b>
<b>14.</b>	<b>Appendix 2- Starting With Serial Number U06081 .....</b>	<b>90</b>
	<b>Typical Relay .....</b>	<b>102</b>
	<b>Press F1 to observe the remaining seven outputs.....</b>	<b>103</b>

## Table of Figures

Figure 2-1 Single NDIR Analyzer.....	20
Figure 2-2 Detector Output Signal .....	21
Figure 2-3 Absorption Characteristic of Detector.....	21
Figure 2-4 Block Diagram .....	22
Figure 3-1 AC Power Switch, Connector, and Fuse.....	24
Figure 3-2 EMI Suppressor.....	Error! Bookmark not defined.
Figure 4-1 LCD Display.....	27
Figure 4-2 Keyboard .....	32
Figure 5-1 Main Menu on Power Up Screen.....	34
Figure 5-2 Main User Menu (Level 4) .....	34
Figure 7-1 Main Menu Screen .....	36
Figure 7-2 Measurements Screen .....	36
Figure 7-3 Set to Auto Range (601).....	37
Figure 7-4 Analyzer set to Range 1 (603).....	37
Figure 7-5 Main Menu (User Level 4) .....	38
Figure 7-6 Purge Screen .....	38
Figure 7-7 Diagnostics Screen I .....	38
Figure 7-8 Calibration Screen.....	39
Figure 7-9 Auto Calibration Screen .....	39
Figure 7-10 Manual calibration .....	40
Figure 7-11 Manual zero calibration .....	40
Figure 7-12 Manual span calibration .....	40
Figure 7-13 Display deviations .....	41
Figure 7-14 Zero gas deviations .....	41
Figure 7-15 Span gas deviations .....	41
Figure 7-16 Setup menu screen I .....	43
Figure 7-17 setup menu screen II.....	43
Figure 7-18 Change Span Gas Settings .....	44
Figure 7-19 Change Auto Calibration Settings .....	44
Figure 7-20 Change Range Limits .....	45
Figure 7-21 Change Upper Range Limits .....	45
Figure 7-22 Change Auto Range Limits.....	45
Figure 7-23 Alarm screen I .....	46
Figure 7-24 Alarm screen II .....	46
Figure 7-25 Enter / Change Password.....	46
Figure 7-26 Access Level Screen.....	47
Figure 7-27 Linearization Screen.....	48
Figure 7-28 Coefficients Range Select .....	48
Figure 7-29 Change Coefficients .....	48
Figure 7-30 Linearized and raw data .....	48

Figure 7-31 System Setup Screen.....	49
Figure 7-32 Clock and Timing Setup Screen.....	49
Figure 7-33 Clock set screen.....	49
Figure 7-34 Set Auto Cal Timing.....	50
Figure 7-35 Set Auto Cal Ranges.....	50
Figure 7-36 F4 Toggles Auto Cal ON of OFF.....	50
Figure 7-37 TCP/IP Address.....	51
Figure 7-38 Output Assignments.....	51
Figure 7-39 Output Ranges.....	51
Figure 7-40 Status line.....	52
Figure 7-41 Status line on/off.....	52
Figure 7-42 Measure setup.....	52
Figure 7-43 Lowpass filter.....	52
Figure 7-44 Purge Time.....	53
Figure 7-45 T-compensation.....	53
Figure 7-46 P-Compensation.....	53
Figure 7-47 Low pass filter time constant.....	54
Figure 7-48 Analyzer Information.....	54
Figure 7-49 Software Version.....	54
Figure 7-50 F6 Remote/Manual.....	54
Figure 7-51 F7 Standby.....	55
Figure 7-52 Standby Mode.....	55
Figure 8-1 Internal components.....	59
Figure 8-2 Main Electronic Board Potentiometers.....	60
Figure 8-3 NDIR detector assembly.....	61
Figure 10-1 Optical Bench with Pipe Sample Cell.....	64
Figure 10-2 Optical Bench with Block Sample Cell.....	66
Figure 10-3 Optical Bench with Pipe and Block Type Sample Cells.....	68
Figure 11-1 Amplifier circuit ac wave form.....	75
Figure 12-1 Serial Interface.....	76

## Tables

Table 10-1 Optical Bench Configuration.....	63
Table 12-1 Structure of an instruction command.....	77
Table 12-2 Structure of an acknowledgement command.....	77
Table 12-3 Acknowledgment response in case of error.....	78

## 1. Introduction

### 1.1. *Overview*

Congratulations and thank you! You have just purchased one of the most reliable gas analyzers in the world. Before using the analyzer, please familiarize yourself with its operation by reading this manual. If you have any questions, please do not hesitate to call California Analytical Instruments for assistance. We want you to be a member of our thousands of satisfied customers.

### 1.2. *Unpacking Instructions*

Open the shipping container and carefully remove the analyzer from the packing materials. Inspect the instrument for any sign of damage. Remove the Top Cover retaining screws. Visually check for loose parts or connectors that are not properly seated. Verify all circuit boards and circuit board connections are secure. If all internal components look normal, re-install the cover.

### 1.3. *Reporting Damage*

Should there be any apparent damage either to the inside or outside of the instrument due to shipping or handling, immediately notify the shipper. The shipping container or packing materials should be retained for inspection by the shipper.

### 1.4. *Contact Information*

California Analytical Instruments, Inc.

1312 West Grove Avenue

Orange, CA 92865

714 974-5560

Fax 714 921-2531

Website: [www.gasanalyzers.com](http://www.gasanalyzers.com)

**1.5. Warranty Certificate**

Subject to the exceptions and upon the conditions stated below, California Analytical Instruments (CAI) warrants that the products sold under this sales order shall be free from defects in workmanship and materials for one year after delivery of the product to the original Buyer by CAI and if any such product should prove to be defective within such one year period, CAI agrees, at its option, either (i) to correct by repair or, at CAI's election, by replacement with equivalent product any such defective product, provided that investigation and factory inspection discloses that such defect developed under normal and proper uses, or (ii) to refund the purchase price. The exceptions and conditions mentioned above are as follows:

- a. components or accessories manufactured by CAI that by their nature are not intended to or will not function for one year are warranted only to give reasonable service for a reasonable time. What constitutes reasonable time and reasonable services shall be determined solely by CAI. A complete list of such components and accessories is maintained at the factory;
- b. CAI makes no warranty with respect to components or accessories not manufactured by it; in the event of defect in any such component or accessory CAI will give reasonable assistance to Buyer in obtaining from the respective manufacturer whatever adjustment is authorized by the manufacturer's warranty;
- c. any product claimed to be defective must be returned to the factory transportation charges prepaid and CAI will return the repaired or replaced product freight collect;
- d. if the product claimed to be defective requires on-site repair, such warranty labor will be provided at no charge; however, transportation and living expenses will be charged to Buyer;
- e. if the product is a consumable or the like, it is warranted only to conform to the quantity and content and for the period (but not in excess of one year) stated on the label at the time of delivery or 90 days;
- f. CAI may from time to time provide a special printed warranty with respect to a certain product, and where applicable, such warranty shall be deemed incorporated herein by reference;
- g. CAI shall be released from all obligations under all warranties, either expressed or implied, if any product covered hereby is repaired or modified by persons other than its own authorized service personnel unless such repair by others is made with the written consent of CAI.

IT IS EXPRESSLY AGREED THAT THE ABOVE WARRANTY SHALL BE IN LIEU OF ALL WARRANTIES OF FITNESS AND OF THE WARRANTY OF MERCHANTABILITY AND THAT CAI SHALL HAVE NO LIABILITY FOR SPECIAL OR CONSEQUENTIAL DAMAGES OF ANY KIND OR FROM ANY CAUSE WHATSOEVER ARISING OUT OF THE MANUFACTURE USE, SALE, HANDLING, REPAIR, MAINTENANCE OR REPLACEMENT OF ANY OF THE PRODUCTS SOLD UNDER THIS SALES ORDER. SOME STATES DO NOT ALLOW THE EXCLUSION OR LIMITATION OF INCIDENTAL OR CONSEQUENTIAL DAMAGES, SO THAT THE ABOVE LIMITATIONS OR EXCLUSIONS MAY NOT APPLY. THIS WARRANTY GIVES YOU SPECIFIC LEGAL RIGHTS, AND YOU MAY HAVE OTHER RIGHTS, WHICH VARY FROM STATE TO STATE.

Representations and warranties made by any person, including dealers and representatives of CAI which are inconsistent or in conflict with the terms of this warranty, shall not be binding upon CAI unless produced in writing and approved by an expressly authorized officer of CAI.

**1.6. Proper Operation**

Personnel should be trained in the proper operation of this equipment before attempting to operate the equipment.

**1.7. Possible Explosion Hazard**

Do not apply power to the analyzer or attempt to energize the analyzer until determining the analyzer environment to be non-hazardous.

Use this analyzer in a **NON-HAZARDOUS** environment.

This analyzer has not been designed for use with a hazardous sample.

Tampering or use of substitute components may cause a safety hazard. Use only factory authorized replacement parts.

**1.8. Electrical Shock Hazard**

Disconnect power before removing cover. Servicing requires access to live electrical components that can cause death or serious injury. Refer servicing to qualified service personnel. For safety and proper performance, connect this instrument to a properly grounded three-wire receptacle.

**Caution****1.9. Plug Removal**

***Do not operate this analyzer UNTIL REMOVING the red plastic ¼ inch plugs from the sample inlet and exhaust fittings on the rear panel.***

**Use of this equipment in a manner not approved by California Analytical Instruments is not recommended and may cause harm to the equipment or operating personnel.**

## 2. Features

### 2.1. Description

The Model 600 series of NDIR analyzers incorporate a single-beam photometric system and a detector with a microflow sensor assuring high reliability, sensitivity, accuracy, and stability. The microflow detector is a sealed unit filled with the same gas as the component of interest (CO, CO<sub>2</sub>, and CH<sub>4</sub>). The length of the sample cell determines the most sensitive range for each component.

### 2.2. Features-General

The Model 600 series of NDIR analyzers have a 3 by 5 inch liquid crystal display and a 20 key data/operation input keyboard. The 16-bit microprocessor control board consists of the MSR-Card with 16 digital inputs, 16 digital outputs, 16 analog inputs and 4 analog outputs. The analyzer can be manually operated from the keypad or remotely via TCP/IP or RS-232C communications. After turning on the analyzer, it needs at least 30 seconds for initialization. During this time, the screen is illuminated.

**IMPORTANT TIP: When the analyzer is powered up, it defaults to access level 1 (User). To operate ALL parameters, check the access level. See Section 5.**

### 2.3. Features-General

High stability is provided by an improved photometric system, which assures less influence due to contamination of the measuring cell and higher long-term stability than conventional dual-beam analyzers.

- A dual-chamber type detector effectively minimizes influence due to concomitant gas components.
- A microflow sensor within the detector features high reliability, long service life, very low noise, and excellent resistance to vibration.
- The easily serviced single-beam photometric system does not require delicate adjustment of the optical-balance.
- Simple construction assures reliable performance.
- Modular component design simplifies maintenance. Independent elements are easily removed for maintenance.
- Low Power Consumption - The instruments are of energy-saving design with power consumption as low as 30 VA.

### 2.4. Infrared Gas Analyzers


The infrared gas analyzer measures gas concentration based on the principle that each type of gas component shows a unique absorption line spectrum in the infrared region.

The instrument consists of an infrared light source, a chopper, a measuring cell, and a detector filled with a gas mixture containing the gas component to be measured. The operating principle of the instrument is described with reference to Figure 2-1. The infrared light source emits infrared light in all directions. The light emitted forward is transmitted and reflected into the detectors.

The infrared light emitted backward is reflected by a reflecting surface and is added to the infrared light emitted forward. Arranged between the infrared light source and measuring cell is a chopper blade which rotates to modulate the infrared light beam at regular frequency. The modulated infrared light beam thus formed passes through the measuring cell filled with a sample gas where the light energy is partially absorbed or attenuated before it reaches the front chamber of the detector. Both the front and rear chambers of the detector are filled with the gas component to be measured.

The infrared light energy is partially absorbed in the front chamber and residual light is absorbed in the rear chamber, thereby increasing pressure in both chambers. Since the detector is designed to produce a pressure difference between the front and rear chambers, a slight gas flow is produced through a path connecting these chambers with each other.

This slight flow is converted into an AC electrical signal by a microflow sensor arranged in the path connecting the chambers with each other. The AC signal is amplified and rectified to a DC voltage supplied to the output terminals and indicator (Figure 2-4). Figure 2-2 shows the detector output signal with the greatest amplitude when zero gas is flowing in the measuring cell. Amplitude is reduced as the concentration of measured gas component increases (Figure 2-3)

**Figure 2-1 Single NDIR Analyzer**


Figure 2-2 Detector Output Signal


Figure 2-3 Absorption Characteristic of Detector

### 2.5. Interference Gases

Whenever a sample gas contains a gas component that has an absorption spectrum that overlaps the spectrum of the gas to be measured, that gas is commonly referred to as an interference gas. The Microflow Detector was specifically designed to minimize the effect of interference gases. When these gases are present, the pressure rises in the front and rear chamber of the detector cancel each other minimizing any response to the interference gases.

### 2.6. Electronics

The sinusoidal output signal of the detector's microflow sensor is transmitted to the AC amplifiers on the main circuit board. The signal frequency is related to the rate of the beam-interrupting chopper blade. The signal amplitude is related to the measured gas concentrations in the sample cell.

This signal is amplified by successive AC amplifiers and then demodulated and filtered. The resulting DC signal is further amplified and fed into two output buffer amplifiers. The DC signal output of the printed circuit board is the input to a microprocessor. Here it is digitized and linearized for digital display

The digitized information is then fed to a D/A Converter so it can be isolated and converted to a 0-10 VDC or 4-20 mA output. This output (along with optional alarm contacts) is sent to the 28-pin output connector located on the rear panel of the analyzer for customer connection.


Figure 2-4 Block Diagram

**2.7. Model 600 NDIR Specifications**

<b>IR ANALYSIS METHOD</b>	Non-Dispersive Infrared (NDIR)
<b>NDIR COMPONENTS</b>	CO /CO <sub>2</sub> / CH <sub>4</sub> / C <sub>3</sub> H <sub>8</sub> /SO <sub>2</sub>
<b>DETECTOR TYPE</b>	Microflow
<b>RANGE RATIO</b>	50:1 (Highest Range/50 = Lowest Range)
<b>RESPONSE TIME (IR)</b>	T90 < 2 Seconds to 60 Seconds Adjustable (Depending on configuration)
<b>IR SAMPLE CELL</b>	Stainless Steel with Replaceable Gold Cell Liner
<b>RESOLUTION</b>	Displays Five Significant Digits
<b>REPEATABILITY</b>	Better than 1.0% of Full Scale
<b>LINEARITY</b>	Better than 0.5% of Full Scale of Factory Calibrated Ranges
<b>NOISE</b>	Less than 1% of Full Scale of Factory Calibrated Ranges
<b>ZERO &amp; SPAN DRIFT</b>	Less than 1% of Full Scale per 24 Hours
<b>ZERO &amp; SPAN ADJUSTMENT</b>	Via front panel, TCP/IP or RS-232
<b>SAMPLE FLOW RATE</b>	0.25 to 2.0 Liters/minute (LPM) (Consult Factory for other flow rates)
<b>OXYGEN ANALYSIS METHOD</b>	Paramagnetic
<b>O<sub>2</sub> RANGES</b>	0 — 1% up to 0 – 100% O <sub>2</sub> Full Scale, Four Definable Ranges
<b>O<sub>2</sub> RESPONSE TIME</b>	90 < 2 Seconds

**2.8. Model 600 NDIR Features**

<b>OUTPUTS AVAILABLE</b>	TCP/IP, RS232, Four Scalable Analog 0-10 V / 4-20 mA (Allows Offset and Expandable Range DC Analog Outputs)	
<b>DISCRETE CONTROL</b>	Remote/Local Control, Range Change, Range Sense Mode (AI TTL Logic)	
<b>DISCRETE ALARMS (Local &amp; Remote Adjustable)</b>	General Fault/ TTL Logic 0-5 VDC(Ground True) Calibration Failure/ TTL Logic 0-5 VDC (Ground True) High Concentration (2 each)/ TTL 0-5 VDC Logic (Ground True)	
<b>KEYPAD DISPLAYS</b>	Factory Settings TCP/IP Address Passwords (4)	Scalable Analog Output Voltages Full Scale Range Select Auto Cal Times
<b>SPECIAL FEATURES</b>	Auto Ranging Auto Calibration (adjustable through internal clock)	
<b>DISPLAY</b>	3" x 5" Back lit LCD	
<b>SAMPLE TEMPERATURE</b>	Up to 50°C Non-condensing	
<b>AMBIENT TEMPERATURE</b>	5 to 40°C	
<b>AMBIENT HUMIDITY</b>	Less than 90% RH Non-condensing	
<b>WARM-UP TIME</b>	1 Hour (Typical)	
<b>FITTINGS</b>	1/4 Inch Tube	
<b>POWER REQUIREMENTS</b>	115 V 60Hz (Option: 230V 50Hz), ±10%, 600W	
<b>DIMENSIONS</b>	5¼ H x 19 W x 23 D (Inches)	
<b>WEIGHT</b>	30-45 Pounds (Depending on configuration)	

**SPECIFICATIONS ARE SUBJECT TO CHANGE WITHOUT NOTICE**

### 3. Installation

#### 3.1. General

The design of this instrument is for industrial applications. These installation instructions are for a typical site. Direct any questions regarding specific installation situations to Technical Service of California Analytical Instruments, Inc.

#### 3.2. Site and Mounting

**NOTE: The carefully observe the following precautions:**

Select a site free from direct sunlight, radiation from a high temperature surface, or abrupt temperature variations.

This analyzer is not suitable for installation outdoors.

Select a site where the air is clean. Avoid exposing the instrument to corrosive or combustible gases.

The instrument must not be subject to severe vibration. If severe vibration is present, use isolation mounts.

The instrument is designed for rack-mounting. Optional rack mount slides are available.

Do not install near equipment emitting electromagnetic interference (EMI).

**NOTE: A rear supporting brace or equivalent is required if the optional rack mount slides were not purchased.**

**The power on/off switch is accessible from the rear of the instrument only. DO NOT mount such that the power on/off switch is inaccessible.**

#### 3.3. Electrical

All wiring is connected at the rear of the instrument. The AC power is connected to the power/fuse/switch as shown below


Figure 3-1 AC Power Switch, Connector, and Fuse.

**NOTE: A defective ground may affect the operation of the instrument. The output voltages are connected per Table 8.1.1. Shielded wiring is recommended for output signals.**

**Replace fuses with recommended fuse size indicated on rear panel of instrument. Replacement with any other size fuse may cause damage to the instrument and possible injury to operating personnel.**

**3.4. *Analog Output Connections (Appendix)***

See Appendix for connector pinouts located on the analyzer rear panel. Remote range identification and range selection are obtained via the rear panel connections. When a range is selected, the corresponding control line is pulled low to zero VDC. Ranges not selected will remain at approximately 5 VDC. When remote range control is selected on the front panel switch, a contact closure is provided at the rear panel connector. Remote range selection is made by connection of the control line for the desired range to the analyzers zero VDC line provided in the connector. Five VDC is also provided.

**3.5. *Gases***

- 1) Nitrogen or (zero air) in pressurized cylinder.
- 2) Standard span gas(es) near full-scale concentration with a nitrogen balance, in a pressurized, certified cylinder.

### 3.6. Gas Handling Equipment

- 1) Pressure regulators for zero gas (Air or N<sub>2</sub>), and span gas cylinders.
- 2) Corrosive resistant gas tubing.

### 3.7. Gas Connections

The tubing from the sampling system to the gas analyzer should be made from corrosive-resistant material such as Teflon or stainless steel. Even when the gases being sampled are corrosive themselves, rubber, or soft vinyl tubing should not be used since readings may be inaccurate due to gas absorption into the piping material. To obtain fast response, the tube should be as short as possible. Optimum tube internal diameter is 0.16 inch (4 mm). Couplings to the instrument are ¼-Inch tube.

**Note: Be sure tubing and joints are clean. Dust entering the instrument may cause it to malfunction.**

### 3.8. Sampling Requirements

#### 3.8.1. Filtration

Dust must be eliminated completely. Use filters as necessary. The final filter must be capable of removing particles larger than 4 microns.

#### 3.8.2. Condensation

Dew point of the sample gases must be lower than the instrument temperature to prevent accidental condensation within the instrument. Bypass the sample through a dehumidifier to reduce the dew point to about 2 to 4° C or less. If the sample contains an acid mist, use an acid mist filter, cooler or similar device to remove all traces of the mist.

#### 3.8.3. Presence of Corrosive Gases

High concentrations of corrosive gases such as Cl<sub>2</sub>, SO<sub>2</sub>, F<sub>2</sub>, HCl, in the sample gas shorten the useful service life of the instrument.

#### 3.8.4. Gas Temperature

Do not exceed the maximum rating of the instrument 104° F (50° C) when measuring high temperature gases.

#### 3.8.5. Sample Gas Bypass Outlet (Vent)

A sample gas bypass outlet connector is located on the rear panel (¼ Inch Tube). Pressure at this outlet should be kept at atmospheric level. **ANY** backpressure will cause an error in reading.


## 4. Basic Operation

The operation of the digital microprocessor conforms to the guidelines of the AK committee, originally developed in the German automotive industry. Via the serial port of the MSR-Card, the analyzer can be remote-controlled by a master computer. The serial communication fully corresponds to the specifications of the AK protocol. TCP/IP communication is also available.

### Display

The analyzer's LCD display includes 16 lines with 30 characters each. The display also has background lighting that can be switched on and off via the Display key on the keyboard. The following example shows the measurement screen that is formatted into 4 information areas.


**Figure 4-1 LCD Display**


- 1) THE TOP INFORMATION AREA CONTAINS:
- 2) The AK Protocol Information. This capability is for advanced uses and may be toggled on and off in the setup screen, F5. Next to the symbol for the active operating mode, the device status is indicated. The status field is also displayed on all other screens.
- 3) SARE Auto range enabled
- 4) SMGA Measuring gas is flowing
- 5) SMAN Device is in manual operation status
- 6) SWET Device
- 7) Shown on the right is the Password Entry level with 1 to 4 horizontal lines.
- 8) THE LARGE INFORMATION AREA CONTAINS the data portion of the screen.
- 9) THE THIRD INFORMATION AREA CONTAINS the help information for the parameter selected, ranges, etc.
- 10) THE LOWER INFORMATION AREA CONTAINS The time and date and any error condition.
- 11) The symbol in the bottom right corner indicates the keyboard mode. In the example shown, the keyboard is in the function key mode. For input fields, the mode is usually switched to numerical input. Then, an


N appears in the lower right of the screen. This symbol is displayed on all screens.

4.1. Menu Trees


Model 601 Menu Tree

## 4.2. Keyboard

**Figure 4-2 Keyboard**

### 4.2.1. Operation with the Cursor Keys and the Enter Key

When operating the unit with the cursor keys, you select the various functions with the up/down cursor keys and start them with the Enter key. This method is particularly suitable for less proficient users since the system displays a short on-line help for nearly every function selected. The actual cursor position is shown as a black horizontal bar.

**TIP:** *If you are not yet familiar with the screens and their fields, just press any cursor key after a screen appears. This moves the cursor from field to field and displays the corresponding online help.*

### 4.2.2. Operation with the Function Keys

When using the function keys (F1 through F10), functions are directly accessed by pressing their corresponding function keys. This method is suitable for the advanced user since it is faster than the operation with the cursor keys.

### 4.2.3. Read/Change Parameters

To read and/or change parameters, you must switch to the parameter input mode by pressing the Enter key after calling the corresponding parameter screen. The input cursor (horizontal bar under the first character) then appears in the active edit field (black background). The cursor can be positioned with the right and left cursor keys, and the value displayed (number or letter) can be changed with the up and down cursor keys or entered directly. Every input has to be concluded by pressing the Enter key again, which causes the cursor to disappear.

## 5. Operating Structure

The analyzer's operation can be divided into 4 operating levels. The current level is always displayed as a stack of 1 to 4 horizontal bars in the top right corner of the screen. In the access level menu, you can choose between the following operating levels:

F1	User	(operating level 1)
F2	Advanced user	(operating level 2)
F3	Maintenance	(operating level 3)
F4	System user	(operating level 4)

A password can be assigned to each operating level. Only the system user, who normally has the highest operating priority, can assign the password. At the factory, the default passwords for the CAI analyzers are set as follows:


User:	111
Advanced user:	222
Maintenance:	333
System:	444

The default setting can be changed only by the system user. This manual is written to include all information for the advanced system user.

***TIP: Because of the user settings, some of the parameters shown in this manual may not appear on your analyzer. Check the access level.***

### 5.1. The Main Menu

Upon power up, the CAI logo is first displayed and then the main menu appears as below:


**NOTE: Access Level Indication**

Figure 5-1 Main Menu on Power Up Screen

**NOTE: F6 is not available because, on initial start up, the analyzer reverts to ONLY Level 1 access. See Section 7.5.5 for Password information.**

All functions can be selected with the cursor keys and activated by pressing the Enter key, or directly with the function keys F1 through F7. A ">" to the right of a function means that one or more sub-menus are available. If this sign is missing, the function starts immediately after the activation.

**NOTE: Access level is 4.**

**NOTE: Access Level Indication**

Figure 5-2 Main User Menu (Level 4)

## 6. Menu Structure

There are four operating levels based on the level of your password. This section shows the access rights of the single levels.

### 6.1. User Functions (Level 1)

Main Menu		Setup Menu		Password Menu	
F1	Measurements	F5	Password	F1	Enter password
F2	Purge Analyzer	F10	Version		
F3	Diagnostics				
F4	Calibrations				
F5	Setup				
F7	Standby				

### 6.2. Advanced User Functions (Level 2)

Main Menu		Setup Menu		Password Menu	
F1	Measurements	F3	Range Limits	F1	Enter password
F2	Purge Analyzer	F5	Password		
F3	Diagnostics	F10	Version		
F4	Calibrations				
F5	Setup				
F7	Standby				

### 6.3. Maintenance Functions (Level 3)

Main Menu		Setup Menu		Password Menu		System Settings Menu	
F1	Measurements	F1	Span Gas Concentration	F1	Enter password	F1	Real Time Clock
F2	Purge Analyzer	F3	Range limits	F2	Reset password	F5	Status Line on/off
F3	Diagnostics	F5	Password	F7	Auto Startup		
F4	Calibrations	F7	System Settings				
F5	Setup	F8	Measure Settings				
F7	Standby	F10	Version				

### 6.4. System User Functions (Level 4)

All Function described in this manual may be accessed from Level 4.

### 7.1.2. Range Select

With the arrow keys, the ranges 1 to 4 can be selected and locked in which will disable the auto range capability. Continue pressing the arrow keys will recycle the analyzer back to auto range. The range and/or auto range is displayed on the measurement screen. If the limits are exceeded while not in the auto range mode, a warning "Over Range" appears on the screen.


Figure 7-3 Set to Auto Range (601)


Figure 7-4 Analyzer set to Range 1 (603)

### 7.2. F2 Purge Analyzer


Figure 7-5 Main Menu (User Level 4)


Figure 7-6 Purge Screen

F2 from the Main Menu activates the Purge (analyzer) function if equipped.

### 7.3. F3 Diagnostics

F3 from the Main Menu activates the Diagnostics function. F3 brings up the two diagnostics screens. The Diagnostics screens may be brought up from **EITHER** the Main Menu or the Measurements screen.


Figure 7-7 Diagnostics Screen I

#### 7.4. F4 Calibrations

F4 from the Main Menu activates the Calibrations screen. Calibrations may be automatic or manual. Deviations can also be displayed. Calibration values can be reset to default values and the range to be calibrated can be changed.


Figure 7-8 Calibration Screen

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From the Calibrations screen, F1 starts automatic calibration. If auto range is selected, the actual range in use will be calibrated. Auto calibration works as follows: First, zero gas is purged a certain time, called purge-time. Then the measurement begins. The measured value must be a minimum-time, called measuring-time, and within an upper and a lower limit to be saved as new offset value. The maximum length of measuring time is 9 seconds. If the measured value was constant during calibration time, it is checked to determine if this value deviates from the preceding value. If the deviations are too large, a warning "Deviation error!" appears and the user can choose if the new value is saved or not. At last, the zero gas is flowed a further time, verifying time, so it can be checked if the signal is still constant. All of these times can be changed. After zero gas calibration, the same happens with span gas. During auto calibration, "Calibration in progress" is displayed. It also shows which gas is flowing and which time runs. When auto calibration has finished it is displayed. If the span value of the selected range is 0 (see section 5.6.1), then it will not be calibrated. If one range is calibrated and the span value for the lower ranges is zero, calibration parameters will be copied to this range. To calibrate all ranges with the same span gas, you must enter the gas concentration in the Span Gas Calibration screen for ALL RANGES. You must also calibrate each range. Offsets and scalors are NOT copied to other ranges.



Figure 7-9 Auto Calibration Screen

### 7.4.2. F2 Manual Calibration

From the Calibration screen, F2 starts manual calibration. If auto range is selected, calibration is not possible, and the appropriate range can be selected.


**Figure 7-10 Manual calibration**


In the manual calibrations menu, two options are possible:

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F2 Flow span gas


**Figure 7-11 Manual zero calibration**

**Figure 7-12 Manual span calibration**

When zero or span gas is flowing, the measured value can be saved by pressing F1. If the screen is left by pressing the buttons "Main" or "Back", the measured value is not saved. Solenoids are closed by pressing F2.


**7.4.3. Display Deviations – from Calibration menu F3**


**Figure 7-13 Display deviations**

After every calibration, the deviations are calculated for zero and for span gas.

- F1 Zero gas deviations
  - F2 Span gas deviations
  - F3 Deviations of zero gas during verifying
  - F4 Deviations of span gas during verifying
- F1 and F2 deviations are displayed in percent.

**Figure 7-14 Zero gas deviations**

**Figure 7-15 Span gas deviations**

During calibration, there is verification for zero and span gas. With option F3 and F4, you can view the deviations during the verification time. Absolute deviation is the absolute average difference from the saved value in ppm. Relative deviation is the absolute average difference in percent, related to the range limit.

**7.4.4. Absolute Zero Gas Deviation**

Absolute zero gas deviation is zero gas content calculated by the factory polynomial related to the range limit of the calibrated range.

**7.4.5. Relative Zero Gas Deviation**

Relative zero gas deviation is the actual deviation minus the deviation of the previous calibration related to the range limit of the calibrated range.

**7.4.6. Absolute Span Gas Deviation**

Absolute span gas deviation is span gas bottle value minus span gas value calculated by the factory-polynomial related to the range limit of the calibrated range.

**7.4.7. Relative Span Gas Deviation**

Relative span gas deviation is the actual deviation minus the deviation of the previous calibration related to the range limit of the calibrated range.

**7.4.8. F4 Check Calibration**

There is a default calibration. Pressing F4 activates an automatic zero and span check for verification.

**7.4.9. F5 Reset Calibration Values**

There is a default calibration. Pressing F5, a new screen appears and asks if the user is sure to reset calibration values to the default calibration values. F1 confirms and the calibration values are reset to default calibration values. F2 leaves this menu without resetting to default values. This function will overwrite all calibrations with factory values. In addition, the linearization polynomial will be overwritten with the factory values.

**7.4.10. F6 Range Select**

This allows a range change to be activated from the calibration menu.

### 7.5. F5 Setup

From the Main Menu, F5 brings up the setup menu. Span gas concentrations, calibration settings, range limits, alarms, password, linearization, system and measure settings can be changed. The Setup menu begins as shown below. A description of each parameter is shown in the information box. NOTE: Use the down arrow key to obtain the additional setup parameters.


Figure 7-16 Setup menu screen I



Figure 7-17 setup menu screen II

#### 7.5.1. F1 Span Gas Concentration


**Note: If you do not have a specific gas for a specific range, the calibration will use the previous ranges calibration. You must have a least one span gas.**

For calibration, it is necessary to input the concentration of the span gas in ppm. For every range, the span gas concentration can be changed. After pressing F1 in the setup menu, a screen appears in which changes can be made. Select with the cursor buttons the range to change. The selected field turns black. To change parameters, switch to parameter input mode by pressing the Enter key. The input cursor (horizontal bar under the first character) then appears in the active edit field (black background). The cursor can be positioned with the right and left cursor keys, and the value displayed (number or letter) can be changed with the up and down cursor keys or entered directly. Every input has to be concluded by pressing the Enter key again. Then the input cursor disappears and a new range can be selected. The changes are saved by leaving the screen by pressing "Main" or "Back." At the right side of the screen, the range limits of the 4 ranges are displayed. They cannot be changed in this screen.

**Figure 7-18 Change Span Gas Settings**

## 7.5.2. Calibration Settings

**Figure 7-19 Change Auto Calibration Settings**

In the calibration settings menu, times and, deviations can be changed.

### 7.5.3. F1 Times

There are four times (in seconds) for auto calibration that can be changed. Purge, measuring, calibration and verifying time. Changes are made and saved as above.

### 7.5.4. F2 Measuring Deviations

During auto calibration, the measured value is only saved if it is within a certain time within an upper and a lower limit. These two limits format a working window. In the setup menu, the deviation is in percent.

### 7.5.5. F3 Deviations

Here you can change absolute and relative deviation in percent. After auto calibration, it is checked to assure the deviations are within this limit. If the deviations are not in this limit, a warning "Deviation error!" appears.

### 7.5.6. F3 Range Limits

There are 4 different ranges. The user can define the upper range limits in ppm.


Figure 7-20 Change Range Limits

### 7.5.7. F1 Range 1-4 (Change Upper Range Limits)

In this menu, the upper range limits can be changed. The new settings are saved by pressing MAIN or BACK. The auto range limits are automatically adapted. This means that if the upper range limit of range 1 for example has reached 90% of the upper range limit in the auto range mode, it is switched automatically to the second range.


Figure 7-21 Change Upper Range Limits

### 7.5.8. F2 Change Auto Range Limits

Although the auto range limits are adapted automatically, it is possible to define them manually. Up means the value when the next higher range is selected in auto range mode, down the value when the next lower range is selected.


Figure 7-22 Change Auto Range Limits

### 7.5.9. F4 Alarms

Error reports are always displayed in the lowest line of the screen. There are two pressures, three temperatures, one concentration, and two voltages with alarm limits that can be defined. The user can define the range limits and, if exceeded, will display an error-message.


Figure 7-23 Alarm screen I


Figure 7-24 Alarm screen II

Set Temperature Alarms Set Concentration, Pressure and Voltage Alarms

### 7.5.10. F5 Password

After turning on the analyzer, you are in access level 1. To change the access level or to change the passwords, press F5 (Setup) in the main menu and Press F5 (Password) again. The following screen appears:

Figure 7-25 Enter / Change Password

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To change access level, press F1. The following screen appears:

**Figure 7-26 Access Level Screen**

F1 to F4 selects an access level. Move the cursor to the access level to be modified. You must enter the correct password for the access level desired. The passwords for the various operation levels consist of three numbers that must be entered on the numeric keypad. If the code word is incorrect, you are asked to re-enter the codeword.

**IMPORTANT TIP:** *When a new analyzer is powered up, it defaults to access level 1 (User). To operate ALL parameters and gain complete access, select F4. Press the Enter key twice and enter 444.*

### 7.5.12. F2 Change Password

The passwords can only be changed, if you are in access level 4. After F2, enter your new 3 digit passwords.

**IMPORTANT TIP:** *You MUST remember and record this new password. If this is lost, you will need to consult the factory for the default password!*

### 7.5.13. F3 Reset Passwords

The passwords can only be changed, if you are in access level 4. Reset passwords will revert to the factory defaults.

7.5.14. F6 Linearization

Pressing F6 on the Setup screen brings up the Linearization screen. The analyzer can be linearized by a polynomial with 5 coefficients. By pressing F1, these 5 coefficients can be changed for each range. By pressing F2, the raw value can be displayed. This is the value before linearization and offset span correction. There are two values on the screen: The value at the top is the linearized, offset-span-corrected value, and the other value is the raw-value.


Figure 7-27 Linearization Screen


Figure 7-28 Coefficients Range Select


Figure 7-29 Change Coefficients


Figure 7-30 Linearized and raw data

## 7.6. F7 System Settings

This screen allows all the system settings to be displayed and modified.


Figure 7-31 System Setup Screen

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This brings up the clock time set screen; auto cal and auto cal enable screens.


Figure 7-32 Clock and Timing Setup Screen

Figure 7-33 Clock set screen

The current time may be set by using the cursor to highlight the entry and using the numeric keys to change the values.

F2 brings up the auto cal time set. As above, the date and times can be set by using the cursor to highlight the entry and using the numeric keys to change the values. F3 Sets autocalibration ranges.

```

Setup RTC          K3:SMAN
F1 Set time
F2 Set auto-calibration time
F3 Select calibration range
F4 Autocalibration on/off

Autocalibration is enabled

F10 Show time

Setup scheduler for auto
calibration

Mon Nov 17 11:04:34 2003

```

```

schedule          K1:SMAN
Starttime : 2: 0: 0 on Wed
Date      : 11. Jun 2003
Every : 1 hour(s).
F1 - change to weekly
F2 - change to daily
F3 - change to hourly

press MAIN or BACK to exit
F1, F2, F3 to save changes

Mon Jan 01 01:09:16 2003

```

Figure 7-34 Set Auto Cal Timing

```

Range selection K1:STBY
Enter Range [0..4] : 1
Enter Channel : 1

Mon Jan 01 01:09:21 2003

```

Figure 7-35 Set Auto Cal Ranges

```

Setup RTC          K3:SMAN
F1 Set time
F2 Set auto-calibration time
F3 Select calibration range
F4 Autocalibration on/off

Autocalibration is enabled

F10 Show time

enable/disable
autocalibration

Mon Nov 17 11:04:43 2003

```

Figure 7-36 F4 Toggles Auto Cal ON of OFF

### 7.6.2. System Setup F2 Displays TCP/IP Address

```

TCP/IP setup      K1: SMGA
-----
IP-address: 192.000.000.220
netmask   : 255.255.255.000
Port      : 0
Gateway   : 000.000.000.000
WinIfPort: 2000

HWaddress  : 00.E0.4B.01.9D.F9

enter IP-Address
take effect after reboot

Mon Nov 17 11:05:02 2003

```

Figure 7-37 TCP/IP Address

### 7.6.3. Systems Setup F3 Displays Output Signal Assignments (Used to Adjust Analog Output Channels)

```

assignment      K1: SARA
-----
Output  Signal      Error
on AOut
  1      off
  2      off
  3      off
  4      off

Mon Nov 17 11:05:11 2003

```

Figure 7-38 Output Assignments

### 7.6.4. System Setup F4 Displays Output Ranges (Used to Adjust Scale of Analog Output Channels)

```

output ranges   K3: SMGA
-----
Output  lower  Limit  upper  Mode
                upper
  1      0      0      0      mA
  2      0      0      0      mA
  3      0      0      0      mA
  4      0      0      0      mA

enter 0 in both fields to
use the default range.

Mon Nov 17 11:05:27 2003

```

Figure 7-39 Output Ranges

### 7.6.5. F5 Turns Status Line On or Off

The status line displays the AK Protocol action on the top line of the display.

K3: SMGA

Figure 7-40 Status line


Figure 7-41 Status line on/off

### 7.6.6. F8 Measure Settings

This screen allows several of the system settings to be displayed and modified.


Figure 7-42 Measure setup

### 7.6.7. F1 Set Lowpass filter


Figure 7-43 Lowpass filter

### 7.6.8. F2 Purge Time

F2 on the Menu Settings screen sets the purge time before continuing with a zero or span calibration.


Figure 7-44 Purge Time

### 7.6.9. F3 Set Temperature Compensation


Figure 7-45 T-compensation

### 7.6.10. F4 Pressure Compensation

Figure 7-46 P-Compensation

### 7.6.11. F3 Low Pass Filter Time Constant

F3 on the Menu Settings screen allows the software time constant to be set between 1 and 60 seconds. This is very useful in eliminating noise when measuring low-level concentrations.


Figure 7-47 Low pass filter time constant

### 7.6.12. F10 Displays the Current Analyzer and Software Versions

This displays the analyzer's information, including the factory recommended air and sample pressure settings.


Figure 7-48 Analyzer Information


Figure 7-49 Software Version

## 7.7. F7 Remote / Manual Control

The analyzer can be remote-controlled either by a master computer or via contact closures. The TCP/IP and serial communication fully corresponds to the specifications of the AK protocol. To change remote/manual control, press F6 in the main menu. This toggles between remote and manual control.




Figure 7-50 F6 Remote/Manual

### Main Menu (User Level 4)

**7.8. F8 Standby****Figure 7-51 F7 Standby****Main Menu (User Level 4)**


In Standby mode, pump is turned off and the solenoids are closed. The CAI logo is displayed.

**Figure 7-52 Standby Mode**

## 8. Analyzer Components

### 8.1. Rear Panel

The following details the rear panel connections:

**Figure 32: Rear Panel**

- 1) **Sample Gas Inlet:** Feeds sample gas to the analyzer. ¼ Inch Tube.
- 2) **Sample Gas Bypass Outlet (Vent):** Exhaust for sample. ¼ Inch Tube.
- 3) **TCP/IP Connection:** Connect Network Connector.
- 4) 28 Pin output terminal strip/main 1 (standard).
- 5) 28 Pin output terminal strip/main 2 (optional).
- 6) 28 Pin output terminal strip/main 3 (optional).
- 7) **Power Entry Module:** Power connection, power switch, fuse compartment (2 Amp) **With Rear Panel Power ON/OFF Switch:** Turns ON/OFF line power to instrument
- 8) **Serial Connector:** Connect Serial Connector

## 8.2. Rear Panel Connectors

28 Pin <u>Main</u> Connector Assignments:						28 Pin <u>Auxiliary</u> Connector Assignments:			
Signal Type	Main 1 Analog		Main 2 Optional Analog		Main 3 Optional Analog		Signal Type	601, 602, 603 Analog	
	Pin #		Pin #		Pin #		Spare	Pin #	
A Output	1	GND (analog)	1	GND (analog)	1	GND (analog)	A Input	1	GND (analog)
A Output	2	Channel 1	2	Channel 1	2	Channel 1	A Input	2	External Analog 1
A Output	3	Channel 2	3	Channel 2	3	Channel 2	A Input	3	External Analog 2
A Output	4	Channel 3	4	Channel 3	4	Channel 3	A Input	4	Spare analog
A Output	5	Aux	5	Aux	5	Aux	A Input	5	Spare analog
		<b>Digital</b>		<b>Digital</b>		<b>Digital</b>	<b>Alarms</b>		<b>Digital</b>
D Output	6	GND (Digital)		GND (Digital)		GND (Digital)	D Output	6	GND (Alarm)
D Output	7	Sense AutoRange		Sense AutoRange		Sense AutoRange	D Output	7	General Alarm
D Output *	8	Sense Range 4		Sense Range 4		Sense Range 4	D Output	8	Ch. 1 Conc. 1 Alarm
D Output *	9	Sense Range 3		Sense Range 3		Sense Range 3	D Output	9	Ch. 1 Conc. 2 Alarm
D Output *	10	Sense Range 2		Sense Range 2		Sense Range 2	D Output	10	Ch. 2 Conc. 1 Alarm
D Output *	11	Sense Range 1		Sense Range 1		Sense Range 1	D Output	11	Ch. 2 Conc. 2 Alarm
		* Hi = True							
D Input	12	Set AutoRange	12	Set AutoRange	12	Set AutoRange	D Output	12	Ch. 3 Conc. 1 Alarm
D Input	13	Control Range 1	13	Control Range 1	13	Control Range 1	D Output	13	Ch. 3 Conc. 2 Alarm
D Input	14	Control Range 2	14	Control Range 2	14	Control Range 2	D Output	14	Reserved
D Input	15	Control Range 3	15	Control Range 3	15	Control Range 3	D Output	15	GND (Alarm)
D Input	16	Control Range 4	16	Control Range 4	16	Control Range 4	D Output	16	Calibration Alarm 1
D Input	17	Auto Cal	17	Auto Cal	17	Auto Cal	D Output	17	Calibration Alarm 2
D Input	18	Calibrate	18	Calibrate	18	Calibrate	D Output	18	Calibration Alarm 2
D Input	19	Zero	19	Zero	19	Zero	D Output	19	Reserved
D Input	20	Span	20	Span	20	Span	D Output	20	Spare
D Input	21	Pump	21	Pump	21	Pump	D Output	21	Spare
D Output	22	Zero Gas Flow	22	Zero Gas Flow	22	Zero Gas Flow	D Output	22	Spare
D Output	23	Span Gas Flow	23	Span Gas Flow	23	Span Gas Flow			Spare
D Output	24	Sample Gas Flow	24	Sample Gas Flow	24	Sample Gas Flow			Spare
D Output	25	Local/Remote	25	Local/Remote	25	Local/Remote			Spare
D Output	26	Read Cal Mode	26	Read Cal Mode	26	Read Cal Mode			Spare
D Output	27	Reserved	27	Reserved	27	Reserved			Spare
D Output	28	Reserved	28	Reserved	28	Reserved			Spare

### 8.3. Digital Outputs

#### 8.3.1. RS-232 (Standard 9 Pin DIN Connector)

Pin	Function
1	DCD Carrier Detect
2	RxD Receive Data
3	TxD Transmit Data
4	DTR Data Terminal Ready
5	Ground
6	DSR Data Set Ready
7	RTS Ready to Send
8	CTS Clear to Send
9	RI Ring Indicator


#### 8.3.2. TCP/IP (8 Pin RJ-47 Connector)

Pin	Function
1	TDX+
2	TDX-
3	RXD+
4	Open
5	Open
6	RXD-
7	LNLED
8	LNLED

**IMPORTANT TIP:** For direct connect to a PC a crossover cable is required. Connection to a hub requires a straight cable.

#### 8.4. Internal Component Locations

**Figure 8-1 Internal components**

- 1) Signal conditioning PCB.
- 2) NDIR detector.
- 3) Sample cell.
- 4) IR source and chopper motor assembly.
- 5) DC power supply module.

### 8.5. Main Electronics Board (Potentiometers)



Figure 8-2 Main Electronic Board Potentiometers

1	RP2	Detector voltage adjust
2	RP3	Fine zero adjust
3	RP4	Not used (Fully CW)
4	RP5	Coarse zero adjust
5	RP6	Span adjust

**NOTE:** Potentiometers are clearly labeled on both sides of the Wave PCB.

**8.6. NDIR Detector Assembly****Figure 8-3 NDIR detector assembly**

- 1) Detector Assembly.
- 2) Sample cell assembly.
- 3) Light source and chopper motor assembly
- 4) Sample inlet.
- 5) Sample outlet.
- 6) Temperature sensing element (to J2 of Wave PCB)

## 9. Operation

### 9.1. External Wiring

Make sure that the external wires have been connected as described in Section 3 Installation.

### 9.2. External Piping

Review Section 3, 3.7 and 3.8

### 9.3. Operation & Calibration

#### 9.3.1. Power On:

Turn on the power switch (located on the rear panel). The LCD display should illuminate. Allow the instrument to warm up for approximately one hour. It is preferable, but not essential, that zero gas flow through the instrument at a rate of about 1 Liter/minute (LPM).

**Note: To achieve final stability, the analyzer may require an additional warm-up period of up to four hours (depending on variables in the analyzer's environment).**

#### 9.3.2. Zero/San Adjustment:

Follow the "Quick Start Guide" at the front of this manual.

**Note: Span gas concentration should not be less than 80% of the range to be spanned.**

#### 9.3.3. Start-Up and Routine Maintenance:

Prepare and check the sample system. Adjust the flow of sample gas to about 1 L/min. The instrument should show a meter indication. The infrared gas analyzer is designed for extended operation and may be left switched on continuously.

## 10. MAINTENANCE

### Warning

All replacement parts must be as supplied and/or specified by California Analytical Instruments. Failure to use specified parts may reduce the safety features of the instrument or create a hazardous condition.

#### 10.1. Zero and Span Calibration

The zero and span levels should be checked and/or calibrated daily (or as often as required.)

#### 10.2. Cleaning of the Optical Bench Measuring Cell (Infrared Analyzers Only)

Dust or water droplets entering the measuring cell may cause drift due to contamination. When it is impossible to adjust the zero level with the zero control mounted on the front panel, check the measuring cell for contamination. If contamination is present, check the sampling system, especially the filters, to eliminate the source of contamination. Periodic maintenance is generally not required. Cleaning is accomplished by use of a cleaning agent (such as isopropyl alcohol or household glass cleaner) and a non-abrasive, lint-free cloth or tissue.

#### 10.3. Optical Bench Configuration

Infrared analyzers may be configured with three types of optical benches: with a pipe cell; with a block cell; or with both a pipe cell and a block cell (see Table 10-1))

**Table 10-1 Optical Bench Configuration**

Optical Bench Type	Illustration	Paragraph
Pipe Cell (Cell length: 64 mm, 125 mm, or 250 mm).	Figure 10-1 This figure illustrates the pipe cell	10.3.1 describes the disassembly cleaning and reassembly
Block Cell (Cell length: 4 mm, 8 mm, 16 mm or 32 mm)	Figure 10-2 This figure illustrates the block cell.	10.3.2 describes the disassembly cleaning and reassembly
Combination	Figure 10-3 This illustrates the combination assembly.	10.3.3 describes the disassembly cleaning and reassembly


Figure 10-1 Optical Bench with Pipe Sample Cell

Item	Description	Item	Description
1-3	Screw (M4)	12	Pipe sample cell
4	Base	13	O-Ring
5	Infrared source unit	14	Window
6-7	Screw (M4)	15	Detector
8	Connector	16	Bridge circuit board
9	Sealed tubes	17	Optical filter (if installed)
10	Support	18	Wave PCB
11	Clamp		

**Note:** Pipe cell foil liner not shown

**10.4. Removal of Pipe Cell**

(Figure 10-1)

- 1) Discontinue the sample gas flow. When it contains harmful gas, purge the measuring cell sufficiently with zero gas.
- 2) Turn the power switch to OFF.
- 3) After loosening the retaining screws on the sides of the top cover, lift off top cover and locate the pipe cell.
- 4) Disconnect the sample gas inlet and outlet tubes from the measuring cell.
- 5) After loosening (do not remove) the two screws (No. 1 in Figure 10-1) which are used for securing the infrared light source unit (No. 5 in Figure 10-1) to the base plate, shift the infrared light source away from the pipe cell (No. 12 in Figure 10-1) to form a gap.
- 6) After removing the pipe cell retainer screws (No. 7 in Figure 10-1), remove the retaining clamps (No. 11 in Figure 10-1).
- 7) Carefully remove the cell from the optical bench and remove both windows (right-hand threaded) (No. 14 in Figure 5-1).
- 8) At this time, inspect the O-Ring (No. 13 in Figure 10-1) for signs of deterioration. Replace if necessary.
- 9) The  $\text{CaF}_2$  window is bonded to the window holder. Inspect and clean the windows as necessary using cotton swabs and a suitable cleaning solution.

***Alcohol or an alcohol-based glass cleaner is a suitable cleaning solution. A soft cloth or tissue that will not deposit lint should be used to clean the liner & windows.***

- 10) The pipe cells contain a reflective metal foil liner (not shown in Figure 10-1) to enhance the light energy through put in the cell. Normally it is not necessary to remove the liner for cleaning; however, the liner should be removed if the cell is subjected to grossly excessive moisture. If necessary, clean both sides of the liner and the inside of the pipe cell. If the liner has been subjected to a corrosive substance, it should be replaced. In either case, ensure that the small gas holes in the liner are aligned with the gas fittings at both ends of the pipe cell before reassembling the windows.
- 11) The pipe cell can be re-assembled by following the reverse of the disassembly procedures. In re-assembly, reserve gaps of approximately 0.5 mm between the infrared light source unit and cell and between the cell and detector, respectively. Larger gaps are undesirable.


Figure 10-2 Optical Bench with Block Sample Cell

Item	Description	Item	Description
1-3	Screw (M4)	11	Block sample cell
4	Base	12	O-Ring
5	Infrared source unit	13	Window
6	Screw (M4)	14	Detector
7	Screw (Round or Pan Head)	15	Bridge circuit board
8	Connector	16	Optical filter (if installed)
9	Sealed tubes	17	Wave PCB
10	Screw (Flat Head)		

**10.5. Removal of Block Cell**

(Figure 10-2)

- 1) Discontinue the sample gas flow. When it contains harmful gas, purge the measuring cell sufficiently with zero gas.
- 2) Turn the power switch to OFF.
- 3) After loosening the retaining screws on the sides of the top cover, lift off the top cover and locate the block cell.
- 4) Disconnect the sample gas inlet and outlet tubes from the measuring cell.
- 5) Disconnect the detector output-cable-connector from the main circuit board.
- 6) Remove the two screws (No. 7 in Figure 10-2) attaching the detector to the infrared light source assembly and remove the detector from the optical bench. The cell is removed together with the detector as a unit.
- 7) While holding the detector in the palm of your hand, remove the two flat-head screws (No. 10 in figure 10-2) which fix the cell to the detector.

***Note: The cell window (No. 10 in figure 10-2) is loose and is only retained by the clamping action between the detector and the block sample cell. Take care not to drop the window when separating the block cell from the detector.***

- 8) Clean the cell interior and CaF<sub>2</sub> windows using a soft cloth or tissue (see Note on previous page). Inspect the O-Ring for flatness or deterioration and replace if necessary.
- 9) The block cell can be re-assembled by following the reverse of the disassembly procedures. Note the orientation of the loose window and O-ring during disassembly.


Figure 10-3 Optical Bench with Pipe and Block Type Sample Cells

Item	Description	Item	Description
1-3	Screw (M4)	14	Window
4	Base	15	Detector
5	Infrared source unit	16	Bridge circuit board
6	Screw (M4)	17	Optical filter (if installed)
7	Screw	18	Screw (Flat Head)
8	Connector	19	Screw (Round or Flat Head)
9	Sealed tubes	20	Block Sample Cell
10	Support	21	O-Ring
11	Clamp	22	Window
12	Pipe sample cell	23	Wave PCB
13	O-Ring		

### 10.6. Disassembly of Combination Pipe & Block Type Cells

(Figure 10-3)

- 1) Discontinue the sample gas flow. When it contains harmful gas, purge the measuring cell sufficiently with zero gas.
- 2) Turn OFF the power switch.
- 3) After loosening the retaining screws on the sides of the top cover, lift off top cover and locate the block cell.
- 4) Disconnect the inlet and outlet tube from both measuring cells.
- 5) Disconnect the wires from the two terminal screws of the infrared source unit and unplug the 2-pin connector (No. 8 in figure 10-3) from the chopper motor.
- 6) Disconnect the block cell's detector output-cable-connector from the main circuit board.
- 7) Remove the two screws (No. 1 in figure 10-3) that secure the infrared source unit to the base.
- 8) Remove the sub-assembly consisting of the source unit, block cell, and detector.
- 9) Remove the two screws (No. 19 in figure 10-3) that secure the detector to the infrared source then separate the detector and block cell from the infrared source unit.
- 10) While holding the detector in the palm of your hand, remove the two flat-head screws (No. 18 in figure 10-3) which fix the cell to the detector.

**Note: the cell window (No. 10 in figure 10-3) is loose and is only retained by the clamping action between the detector and the block ample cell Take care not to drop the window when separating the block cell from the detector.**

- 11) Clean the cell interior and CaF<sub>2</sub> windows using a soft cloth or tissue (see Note on previous page). Inspect the O-Ring for flatness or deterioration and replace if necessary.
- 12) The block cell can be re-assembled by following the reverse of the disassembly procedures. Note the orientation of the loose window and O-ring during disassembly.

**Note: Before re-installing the infrared source unit, block cell and detector sub-assembly you should first remove and clean the pipe cell.**

- 13) Remove the pipe cell retaining screws (No. 7 in figure 10-3)
- 14) Remove the retaining clamps (No. 11 in figure 10-3).
- 15) Lift the pipe cell out of its two supports (No. 10. in figure 10-3) and remove both windows (right-hand threaded) (No. 14 in figure 10-3).
- 16) At this time, inspect the O-Ring (No. 13 in Figure 10-3) for signs of deterioration. Replace if necessary.

- 17) The CaF<sub>2</sub> window is bonded to the window holder. Inspect and clean the windows as necessary using cotton swabs and a suitable cleaning solution.

***Alcohol or an alcohol-based glass cleaner is a suitable cleaning solution. A soft cloth or tissue that will not deposit lint should be used to clean the liner & windows.***

- 18) The pipe cells contain a reflective metal foil liner (not shown in Figure 10-3) to enhance the light energy through put in the cell. Normally it is not necessary to remove the liner for cleaning; however, the liner should be removed if the cell is subjected to grossly excessive moisture. If necessary, clean both sides of the liner and the inside of the pipe cell. If the liner has been subjected to a corrosive substance, it should be replaced. In either case, ensure that the small gas holes in the liner are aligned with the gas fittings at both ends of the pipe cell before reassembling the windows.
- 19) The pipe cell can be re-assembled and installed by placing it in its supports and securing with its retaining brackets.

***Note: Before tightening the retaining brackets, be certain that a gap of approximately 0.5-1.0 mm exists between the pipe cell and the pipe cell detector.***

- 20) Reinstall the infrared source unit, block cell and detector sub-assembly by performing steps e, f, g, h, and I in reverse, being sure to leave a slight gap of approximately 0.5-1.0 mm between the back of the block cell and the front pipe-cell window.

## 11. ADJUSTMENTS CHECKS AND REPAIRS

### 11.1. Adjustment of Detector Voltage (NDIR'S only)

**Note:** Adjustment is required if detector or NDIR amplifier board is replaced. **Important:** Turn RP2 to its maximum counter clockwise position before initial power up of a replacement detector or NDIR.

The detector operating voltage is specified on the label attached to the side of the detector below the Serialization Number and Type designation. On the NDIR amplifier board, the voltage supplied to the detector can be measured between TP-1 (+) and TP-9 (common). Adjust RP2 to achieve the voltage specified on the detector to within 0.001 VDC.

**WARNING:** The detector may be damaged if the applied voltage is excessive. For this reason, it is recommended to adjust RP2 to its maximum counter clockwise position before initial power up of a replacement detector or NDIR amplifier board.

### 11.2. Coarse Zero Adjustment

With zero gas still flowing connect a DC voltmeter to TP7 (+) and TP9 (common). Switch JP8 to the on position and leave JP6 and JP7 in the off position. The target DC voltage for TP7 is 0mV. Adjust RP5 first and use RP3 for the final adjustment.

### 11.3. Span Gain Adjustment

With a full-scale span gas flowing connect a DC voltmeter to TP8 (+) and TP9 (common). The target voltage for TP8 with a full-scale span gas flowing is 8.00VDC. Turn RP4 fully clockwise (12 turns). To achieve this voltage select a resistor value, or combination of values to install into R30, R35, R36, and / or R37. Set the corresponding Jumper (JP9-JP12) to the on position for each resistor used to set the gain. The value needed may range between 1K and 100K ohms. A smaller resistor value increases the gain and achieves a larger DC voltage at TP8 when span gas is flowing. RP6 can be used to make a small gain adjustment when the target voltage is close. When finished setting the gain adjustment, flow zero gas again and repeat the previous step to achieve a good zero at TP8.

**11.4. Check and Repair Detector**

(No. 15 in Figure 10-1 and Figure 10-3, No. 14 in Figure 10-2)

**11.4.1. Problem:**

Microflow sensor broken, bridge resistor defective or gas leak in detector.

**11.4.2. Symptom:**

Zero adjustment impossible.

**11.4.3. Check and/or replace:**

- 1) The microflow sensor and bridge resistors (No. 16 in Figure 10-1) are normal when DC voltages of about 1.5 to 2V are measured between bridge circuit board terminals numbered 1-3 and between terminals numbered 2-3, respectively, and the difference between these voltages is a few tenths of a volt. These two voltage readings can be measured from the top of the NDIR amplifier board between ground (TP9) and pin 1 of J3 and pin 3 of J3 respectively when R7 is removed. J3 is not labeled on the top of the circuit board as it is mounted on the bottom side of the circuit board. The pins are accessible on the top of the board between U3 and JP1. Pin 1 is designated with a square pad.
- 2) Connect an oscilloscope between the NDIR amplifier board check terminals T<sub>Pp</sub> (common) and TP3. When the NDIR amplifier board, infrared light source unit, and chopper are normal, but ac waveform (approximately 10 Hz) is not observable at TP3, gas is leaking in the detector and the detector should be replaced.
- 3) When check in item 1 (above) indicates abnormal voltage, turn the power switch to OFF, and disassemble the NDIR amplifier board from the top of the detector to access the bridge circuit board mounted on the detector.
- 4) Check the microflow sensor for its resistance. Measure resistance values between terminals 1-3 and between terminals 2-3, respectively, on the bridge circuit board. When the resistance values are about 25 to 50 ohms, the microflow sensor is normal, but the bridge resistor is defective. If the resistance is infinite, the microflow sensor is faulty.
- 5) Replacement: Replace the detector with a new one (refer to 10-1). If the instrument is using the pipe cell, the detector is attached to the optical bench from below. Remove the four screws (No. 3 of Figures 10-1 and 10-3) that secure the base (No. 4 of Figures 10-1 and 10-3) to access the detectors retaining screws (No. 2 of figures 10-1 and 10-3).
- 6) After the detector has been replaced, adjust the detector voltage to the specified value (refer to section 11.1 of this manual).
- 7) Adjust zero and span.

**11.5. Check and Repair Infrared Light Source Unit**

(No. 5 in Figure 10-1)

**11.5.1. Problem:**

A faulty infrared light source or a leaky gas seal.

**11.5.2. Symptom:**

The unit reads off scale or the output is unstable.

**11.5.3. Check and/or replace:**

- 1) After turning the power switch to OFF and disconnecting the lead wires from the two terminal screws, measure the resistance between the two terminals. The resistance should normally be about 38 ohms. If the resistance is infinite, the infrared light source is faulty (Output drifts in the negative direction as resistance decreases).
- 2) When indicator output drifts due to influence from atmosphere in spite of normal operations of the detector and main circuit board, gas may be leaking into the infrared light source unit.

**Note:** *In the case of a low-concentration CO<sub>2</sub> analyzer, the indicator output may fluctuate due to atmospheric CO<sub>2</sub> penetrating the gaps in the optical bench. This would be normal and not necessarily indicative of a problem. Purge the analyzer case with N<sub>2</sub> gas.*

- 3) Replacement: After disconnecting wires from the two terminal screws and motor connector, remove the two screws that are used to attach the infrared light source unit to the optical bench. The light source assembly can be replaced, referring to Figure 10-1 or Figure 10-2.
- 4) After replacement of the infrared light source unit, adjust the zero level and span.

**11.6. Check and Replace Chopper**

(See Figure 2-1)

**11.6.1. Problem:**

Rotation abnormal.

**11.6.2. Symptom:**

Indicator output unstable or unresponsive.

**11.6.3. Check and/or replace:**

- 1) With the power switch turned ON, listen for a frictional noise from the chopper blade. If noise can be heard, it is necessary to adjust the chopper blade so that no contact is made with other parts. Remove the infrared light source assembly and detach the protective cover for access to the chopper blade. Take care not to damage the blade, as it is made of thin material. No adjustment is required as long as the output is normal.

- 2) If the motor shaft does not rotate after energizing the instrument, disconnect the power supply connector from the motor and check to see if AC 100 Volts is supplied to the connector on the power supply side. When power is supplied but the motor shaft does not rotate, check the shaft and blade sector for an obstruction. When the motor does not rotate, and there is no abnormal contact on the shaft or sector blade, the motor itself is defective.
- 3) Replacement: When the motor is defective, the infrared light source assembly must be replaced.

### **11.7. Check and Repair Measuring Cell**

Detector Window and Infrared Light Source Window

#### **11.7.1. Problem:**

Cell and window badly contaminated.

#### **11.7.2. Symptom:**

Zero adjustment impossible.

#### **11.7.3. Check and or replace:**

After removing the measuring cell, check the cell and windows for contamination. If contaminated, remove contaminant with a soft cloth and alcohol. Take care not to damage the windows since they are fragile. For details, refer to section 10 of this manual.

### **11.8. Check and Repair Tubing Trouble**

#### **11.8.1. Problem:**

Tubing loosened, disconnected, contaminated, or restricted.

#### **11.8.2. Symptom**

Indicator output unstable or response is too slow.

#### **11.8.3. Check and/or replace:**

- 1) When tubing is disconnected or loosened, firmly reconnect it.
- 2) When tubing is contaminated or restricted, disconnect it and blow out contaminants with high-pressure air, or replace the tubing.


### **Caution**

**Cracked, broken or frayed tubing may cause a hazardous condition. Replacement tubing must be as specified by California Analytical Instruments to prevent possible damage to operating personnel or equipment.**

### 11.9. Check and Repair Amplifier Circuit

- 1) Connect an oscilloscope across check terminals TP3 (signal) and TP9 (common) and observe the ac waveform. While zero gas is flowing, the amplitude of the waveform should be between 2.48Vpp and 2.8Vpp as shown below. This can be adjusted by following the procedure described under section 12.2 before making this adjustment it is recommended that a check for contamination in the cell and cell windows be done first, as this can lead to a reduction in amplitude of the ac signal. **Note: Check cell first for contamination and clean if necessary.**

**Figure 11-1 Amplifier circuit ac wave form**


- 2) If an ac waveform cannot be observed in step a) above, observe the ac waveform across terminals TP4 (signal) and TP9 (common). If an ac waveform of 10 Hz and “peak-to-peak” amplitude of 147mV to 164mV is observed, ac amplifier U1 is normal and AC amplifier U7 is not functioning.
- 3) If an AC Waveform is not observed across terminals TP4 (signal) and TP9 (common), check the detector per section 12.4 of this manual.

## 12. Communication Master Computer / Analyzer (AK Protocol)

### 12.1. Serial Interface and AK-Commands

The serial interface enables remote control of the analyzer by a master computer. It is implemented as an RS232 V24 interface and meets all requirements of the AK protocol.

A 9-pin male connector at the back of the unit is used to connect a master computer with the following pin assignment:

Pin 2 = Rxd (receive)

Pin 3 = Txd (transmit)

Pin 5 = Gnd (ground)

Figure 12-1 Serial Interface

### 12.2. Interface Specifications

Speed:	9600 bps
Character Length	1 start bit
	8 data bits
	1 stop bit
Parity:	none
Handshake	no

### 12.3. Protocol Description

#### 12.3.1. Instruction command

	Character	Explanation
1 <sup>st</sup> Byte	STX	ASCII code 02
2 <sup>nd</sup> Byte	Don't Care	Any ASCII code
3 <sup>rd</sup> Byte	Function Code 1	AK instruction e.g.: ASTF
4 <sup>th</sup> Byte	Function Code 2	
5 <sup>th</sup> Byte	Function Code 3	
6 <sup>th</sup> Byte	Function Code 4	
7 <sup>th</sup> Byte	Blank	
8 <sup>th</sup> Byte	K	
9 <sup>th</sup> Byte	0	
10 <sup>th</sup> Byte	Blank	
	D	AK instruction parameters, length is variable
	A	
	T	
	A	
n <sup>th</sup> Byte	ETX	ASCII code 03

**Table 12-1 Structure of an instruction command**

#### 12.3.2. Acknowledgement command

	Character	Explanation
1 <sup>st</sup> Byte	STX	ASCII code 02
2 <sup>nd</sup> Byte	Don't Care	Any ASCII code
3 <sup>rd</sup> Byte	Function Code 1	Echo of the AK instruction command
4 <sup>th</sup> Byte	Function Code 2	
5 <sup>th</sup> Byte	Function Code 3	
6 <sup>th</sup> Byte	Function Code 4	
7 <sup>th</sup> Byte	Blank	
8 <sup>th</sup> Byte	K	
9 <sup>th</sup> Byte	0	
10 <sup>th</sup> Byte	Blank	
	D	AK acknowledgement parameters, length is variable
	A	
	T	
	A	
n <sup>th</sup> Byte	ETX	ASCII code 03

**Table 12-2 Structure of an acknowledgement command**

### 12.3.3. Data Description

Each command begins with STX (Start of Text) in the first byte. The “don't care” byte can be any ASCII character. Generally, a blank or an underscore ( ) is used for readability reasons. The four function bytes represent the AK command. A blank comes next, followed by K and the channel number. The analyzer is a single-channel device, and because of that the channel number is usually 0. For delimiting the command parameters from the channel number, another blank follows. This may be followed by command parameters with a variable length. Every command ends with the ETX (End of Text) character. The error status byte in the acknowledgment command signals if internal errors in the analyzer occurred. It is zero when no error appeared, and it is unequal zero when one or more errors occurred. Every time a change in the errors happens the error status byte is incremented by one, no matter if one or several errors disappear or are added. If it had the value 10, it would be reset to 1. The error status byte does not indicate the real number of errors. If the analyzer does not have errors, the error status byte contains the value 0.

In general, AK commands are subdivided into three classes:

- Control commands (Sxxx)
- Inquiry commands (Axxx)
- Configuration commands (Exxx)

### 12.3.4. Error Handling

It might happen that an unknown instruction is sent, that the analyzer is busy with a function which is not the desired one, or that an error occurred in the command parameters. Table 1-4 summarizes all errors that can appear upon any master instruction.

Analyzers Acknowledgement	Explanation
???? f1	Analyzer does not know the instruction sent.
xxxx <sup>2</sup> f BS	Analyzer is busy with another function.
xxxx f SE	Syntax error within command parameters or incomplete command.
xxxx f NA	Requested function or data are not available.
xxxx f DF	Data error: The kind or number of given parameters are not valid.
xxxx f OF	Offline: Analyzer is offline, i.e. analyzer is in local Mode. Only inquiry commands and SREM (set analyzer in remote mode) are allowed.

**Table 12-3 Acknowledgment response in case of error**

<sup>1</sup> f stands for the error status byte.

<sup>2</sup> xxxx stands for the function code of the sent master command

**12.4. Scan Commands****12.4.1. AKON:** Measured concentration value

Command	Response	Description
_AKON_K0	_AKON_s_z.z_y.y_x.x_w.w	Measured concentration value of all channels is responded
_AKON_Km	_AKON_s_z.z_t	t = Timestamp (1/10 sec) Measured concentration value of channel m is responded t = Timestamp (1/10 sec)

**12.4.2. AEMB:** Get measuring range

Command	Response	Description
_AEMB_K0	_AEMB_s_Mn_Mn_Mn	Current measuring range of all channels is responded
_AEMB_Km	_AKON_s_Mn	Current measuring range of channel m is responded

**12.4.3. AMBE:** Measuring range limit

Command	Response	Description
_AMBE_K0	_AMBE_s_M1_w.w _M2_x.x _M3_y.y _M4_z.z	All existing measuring range limits of channel m are responded
_AMBE_K0_Mn	_AMBE_s_Mn_z.z	Range limit of Range Mn is responded

**12.4.4. AKAK:** Calibration gas concentrations

Command	Response	Description
_AKAK_Km	_AKAK_s_M1_w.w _M2_x.x _M3_y.y _M4_z.z	All existing calibration gas values are responded for selected channel m
_AKAK_Km_Mn	AKAK_s_Mn_z.z	For selected channel m calibration gas value of Range Mn is responded

**12.4.5. AMBU:** Upper and lower range switchover values for auto range

Command	Response	Description
_AMBU_Km	_AMBU_s_M1_w.w_W.W _M2_x.x_X.X _M3_y.y_Y.Y _M4_z.z_Z.Z	Lower and upper range switchover value of auto range are responded for channel m
_AMBU_Km_Mn	_AMBU_s_Mn_w.w_W.W	Lower and upper range switchover value of auto range are responded for channel m range n

**12.4.6. ASTZ: Normal device status**

Command	Response	Description
_ASTZ_K0	_ASTZ_s_K1_State1_State2_State3 _K2_State1_State2_State3 _K3_State1_State2_State3	Respond device status for all channels
_ASTZ_Km	_ASTZ_s_State1_State2_State3	Respond device status only for channel m

**Possible states:**

State 1	State 2	State 3
SREM: remote SMAN: manual	STBY: standby SPAU: pause SMGA: measuring gas SNGA: zero gas SEGA: end gas SATK SNGA: zero gas during auto cal SATK SEGA: end gas during auto cal	SARE: auto range on SARA: auto range off

**12.4.7. ASTF: Error status**

Command	Response	Description
_ASTF_K0	_ASTF_s_f1_f2_f3_..._f10	Current error numbers of all are responded

**Errors:**

1	Channel 1 Flow Failure	8	Channel 1 not calibrated
2	Channel 2 Flow Failure	9	Channel 2 not calibrated
3	Channel 3 Flow Failure	10	Channel 3 not calibrated
4	External Analog 1 Failure	11-13	Ch1...3: Low concentration warning
5	External Analog 2 Failure	14-16	Ch1...3: High concentration warning
6	Pressure Failure	17-19	Ch1...3: Temperature failure
7	Temperature Failure	20-22	Ch1...3: EPC Voltage failure

**12.4.8. AKEN: Device identification**

Command	Response	Description
_AKEN_K0	_AKEN_s_devicename	Device identification is responded
_AKEN_K1	_AKEN_s_model	Device model
_AKEN_K2	_AKEN_s_serialno	Device serial number
_AKEN_K3	_AKEN_s_samplepressure	Suggested input sample pressure

**12.4.9. ARMU: Raw value**

Command	Response	Description
_ARMU_K0	_ARMU_s_z.z_y.y_x.x_t	Raw value before linearization and offset-span-correction is responded for all channels t = Timestamp (1/10 sec) Raw value before linearization and offset-span-correction is responded for channel m t = Timestamp (1/10 sec)

**12.4.10. ATEM: Temperatures**

Command	Response	Description
_ATEM_K0	_ATEM_s_z.z	Device temperature in degrees Celsius is responded
_ATEM_Km	_ATEM_s_z.z	Detector temperature of channel m is returned in z.z

**12.4.11. ADRU: Pressures/ Valve voltage**

Command	Response	Description
_ADRU_K0	_ADRU_s_z.z	Pressure in is responded
_ADRU_Km	_ADRU_s_z.z	EPC voltage of channel m is returned in z.z

**12.4.12. ADUF: Flows**

Command	Response	Description
_ADUF_K0	_ADUF_s_z.z_y.y_x.x	Sample gas flow of all channels is responded
_ADUF_Km	_ADUF_s_z.z	Sample gas flow of channel m is responded

**12.4.13. AGRD: Polynom coefficients**

Command	Response	Description
_AGRD_Km_Mn	_AGRD_s_Mn_a0_a1_a2_a3_a4	Polynomial coefficients of channel m range Mn are responded

**12.4.14. AANG: Deviation from zero point after autocalibration**

Command	Response	Description
_AANG_Km	_AANG_s_M1_z.z_da_dr _M2_z.z_da_dr _M3_z.z_da_dr _M4_z.z_da_dr	Verifying deviations from zero point after auto calibration. Values: measured value (z.z), absolute dev (da), relative dev (dr)

**12.4.15. AAEG: Deviation from end point after autocalibration**

Command	Response	Description
_AAEG_Km	_AANG_s_M1_z.z_da_dr _M2_z.z_da_dr _M3_z.z_da_dr _M4_z.z_da_dr	Verifying deviation from end point after auto calibration Values: measured value (z.z), absolute dev (da), relative dev (dr)

**12.4.16. AFDA: Purge and Autocalibration times**

Command	Response	Description
_AFDA_Km_SATK	_AFDA_s_z_y_x_w	Auto calibration times of channel m: z: Purge time y: Calibration time x: Total calibration time w: Verify time (z, y, x, w in seconds)
_AFDA_K0_SSPL	_AFDA_s_z	Purge time will be responded

**12.4.17. APAR: Request Autocalibration tolerance values**

Command	Response	Description
_APAR_Km_SATK	_APAR_s_z.z_y.y_x.x_w.w	Autocalibration tolerance value (%): z.z: Range 1 y.y: Range 2 x.x: Range 3 w.w: Range 4

**12.4.18. AKAL: Deviations from calibration**

Command	Response	Description
_AKAL_Km_	_AKAL_s_M1_z.z_y.y_x.x_w.w _M2_z.z_y.y_x.x_w.w _M3_z.z_y.y_x.x_w.w _M4_z.z_y.y_x.x_w.w	Deviation: z.z: Zero gas relative last calibration y.y: Zero gas factory calibration. x.x: Span gas relative last calibration. w.w: Span gas factory calibration

**12.4.19. ASYZ: Respond System Time**

Command	Response	Description
_ASYZ_K0_	_ASYZ_s_yymmdd_hhmmss	System time: yymmdd: year, month, day (each 2 characters wide, no spaces) hhmmss: hour, minute, second (each 2 characters wide, no spaces)

**12.4.20. AT90: Respond Lowpass filter time**

Command	Response	Description
_AT90_K0_	_AT90_s_t	Respond low pass filter time t=filter time in seconds

**12.4.21. ADAL: Diagnostic alarm limits**

Command	Response	Description
_ADAL_K0	_ADAL_s_a1.min_a1.max_... _a16.min_a16.max	All alarms limits are responded
_ADAL_K0_x	_ADAL_s_x.min_x.max	Alarm limits of x

**Alarm Limits:**

1	Sample gas flow channel 1	7	Temperature
2	Sample gas flow channel 2	8	Sample concentration channel 1
3	Sample gas flow channel 3	9	Sample concentration channel 2
4	External input 1	10	Sample concentration channel 3
5	External input 2	11-13	Temperature channel 1...3
6	Barometric -Pressure	14-16	EPC voltage channel 1...3

**12.4.22. ATCP: Query TCP/IP settings**

Command	Response	Description
_ATCP_K0	_ADAL_s_zzz.zzz.zzz.zzz _yyy.yyy.yyy.yyy _xxxx	zzz: TCP/IP Address yyy: TCP/IP subnet mask xxxx: TCP/IP port

**12.4.23. AVER:** Query Software version

Command	Response	Description
_AVER_K0	_AVER_s_3MAIN_z_3USER_y_OS MSR_x	z: Main version x.xxx.b_dd.mm.yyyy y: User version x.xxx.b_dd.mm.yyyy x: OSMSR version x.xxx_dd.mm.yyyy

**12.4.24. AH2O:** QueryH<sub>2</sub>O correction parameter

Command	Response	Description
_AH2O_Km	_AH2O_s_z.z_y.y_x.x	z.z: Dry – voltage of A in with no water present y.y: 1st order coefficient x.x.: 2nd order coefficient

**12.4.25. ACO2:** Query CO<sub>2</sub> correction parameter

Command	Response	Description
_ACO2_Km	_ACO2_s_z.z_y.y_x.x_w.w	z.z: Offset – voltage of A in with no CO <sub>2</sub> present y.y: Min A in – if A in is below this value no CO <sub>2</sub> correction will be done. x.x: 1st order coefficient w.w.: 2nd order coefficient

**12.4.26. AUPD:** Query UDP data streaming parameter

Command	Response	Description
_AUDP_K0	_AUDP_s_<UDPPort>_ <DataFrequency>_ [<Mode>] _ [ <UDP_IP>]	Port: port for open the UDP connection DataFrequency: Frequency for transmit the data in Hz Mode: A: ASCII Mode UDP_IP: Alternative IP address for open the UDP connection when it should not use the IP of connected TCP/IP client

**12.5. Control commands****12.5.1. SRES: Reset**

Command	Response	Description
_SRES_K0	_SRES_s	Reset

**12.5.2. SPAU: Pause**

Command	Response	Description
_SPAU_K0	_SPAU_s	Pause mode

**12.5.3. STBY: Standby**

Command	Response	Description
_STBY_K0	_STBY_s	Standby mode for all channels
_STBY_Km	_STBY_s	Standby mode for channel m

**12.5.4. SNGA: Open valve for zero gas calibration**

Command	Response	Description
_SNGA_K0	_SNGA_s	Open all three zero gas valves
_SNGA_Km	_SNGA_s	Open valve for zero gas calibration of actual measuring range
_SNGA_Km_Mn	_SNGA_s	Open valve for zero gas calibration of range Mn

**12.5.5. SEGA: Open valve for end gas calibration**

Command	Response	Description
_SEGA_K0	_SEGA_s	Open all three end gas valves
_SEGA_Km	_SEGA_s	Open valve for end gas calibration of actual measuring range
_SEGA_Km_Mn	_SEGA_s	Open valve for end gas calibration of range Mn

**12.5.6. SSPL: Purge Analyzer with zero gas**

Command	Response	Description
_SSPL_K0	_SSPL_s	Open zero gas valve and purge all channels

**12.5.7. SATK: Start automatic calibration**

Command	Response	Description
_SATK_Km	_SATK_s	Start automatic calibration with selected range of channel m
_SATK_Km_Mn	_SATK_s	Start automatic calibration for channel m, Range n

**12.5.8. SEMB: Set measuring range**

Command	Response	Description
_SEMB_Km_Mn	_SEMB_s	Set measuring range Auto range will be disabled

**12.5.9. SARE:** Auto range on

Command	Response	Description
_SARE_K0	_SARE_s	Set auto range on for all channels
_SARE_Km	_SARE_s	Set auto range on for channel m

**12.5.10. SARA:** Auto range off

Command	Response	Description
_SARA_K0	_SARA_s	Set auto range off for all channels
_SARA_Km	_SARE_s	Set auto range off

**12.5.11. SREM:** Remote mode for AK-commands

Command	Response	Description
_SREM_K0	_SREM_s	Set device in remote mode

**12.5.12. SMAN:** Manual control to control device manually

Command	Response	Description
_SMAN_K0	_SMAN_s	Set device in manual mode

**12.5.13. SMGA:** Start measuring

Command	Response	Description
_SMGA_K0	_SMGA_s	Start measuring Open all sample valves
_SMGA_Km	_SMGA_s	Open sample valve of channel m

**12.5.14. SNKA:** Saves measured value as new offset.

Command	Response	Description
_SNKA_K0	_SNKA_s	Saves measured value of actual range for each channel as new offset if zero valve is opened
_SNKA_Km	_SNKA_s	Saves measured value of actual range as new offset if zero valve is opened

**12.5.15. SEKA:** Saves measured value as new span value

Command	Response	Description
_SEKA_K0	_SEKA_s	Saves new span values for each channel if span valve is opened
_SEKA_Km	_SEKA_s	Saves measured value of actual range as new span value if span valve is opened

**12.5.16. SUDP:** Start /Stop UDP data streaming

Command	Response	Description
_SUDP_K0_ON	_SUDP_s	Start Data streaming via the UDP channel. You need to configure the channel before with EUDP command
_SUDP_K0_OFF	_SUDP_s	Stop streaming via the UDP channel

## 12.6. Settings

**12.6.1. EKAK:** The four span gas concentration values are set

Command	Response	Description
_EKAK_Km_M1_w.w_M2_x.x_M3_y.y_M4_z.z	_EKAK_s	Set end gas values for channel m

**12.6.2. EMBE:** The four measuring range end values are set

Command	Response	Description
_EMBE_Km_M1_w.w_M2_x.x_M3_y.y_M4_z.z	_EMBE_s	Set range limits

**12.6.3. EMBU:** The upper and the lower range switchover for auto range are set

Command	Response	Description
_EMBU_Km_M1_w.w_W.W_M2_x.x_X.X_M3_y.y_Y.Y_M4_z.z_Z.Z	_EMBU_s	Set lower and upper range switchover limits

**12.6.4. EKEN:** Set new device identification and information

Command	Response	Description
_EKEN_K0_new device-name	_EKEN_s	Set new device identification Maximum length of device name is 40 characters

**NOTE:** To change device identification, you must first rename the device to “RESET”. Now a name up to 40 characters can be given.

**NOTE:** The device name must not have any blanks between characters, e.g. “CAI CLD” is not allowed. You can use underscores, e.g.. “CAI\_CLD”.

**12.6.5. EGRD:** Set polynom coefficients

Command	Response	Description
_EGRD_Km_Mn_a0_a1_a2_a3_a4	_EGRD_s	Set polynomial coefficients of range Mn on channel m

**12.6.6. EFDA:** Set autocalibration and purge times

Command	Response	Description
_EFDA_Km_SATK_z_y_x_w	_EFDA_s	Set auto cal. times for channel m: z: Purge time y: Calibration time x: Total calibration time w: Verify time (z, y, x, w in seconds)
_EFDA_K0_SSPL_z	_EFDA_s	Set analyzer purge time to z seconds

**12.6.7. EPAR:** Set autocalibration tolerance values

Command	Response	Description
_EPAR_Km_SATK_z.z_y.y_x.x_w.w	_EPAR_s	Autocalibration tolerance value (%): z.z= Range 1 y.y= Range 2 x.x= Range 3 w.w= Range 4

**12.6.8. ESYZ: Set System Time**

Command	Response	Description
_ESYZ_K0_yymmdd_hhmmss	_ESYA_s	Set system time: yymmdd: year, month, day (each 2 characters wide, no spaces) hhmmss: hour, minutes, seconds (each 2 characters, no spaces)

**12.6.9. ET90: Set Lowpass Filter Time**

Command	Response	Description
_ET90_K0_t	_ET90_s	Set lowpass filter time: t= filter time in seconds

**12.6.10. EDAL: Diagnostic alarm limits**

Command	Response	Description
_EDAL_K0_a1.min_a1.mas_..._a12max	_EDAL_s	Set all alarm limits
_EDAL_K0_x x.min xmax	_EDAL_s	Set alarm limits of x

**Alarm Limits:**

1	Flow of channel 1	7	Temperature
2	Flow of channel 2	8	Sample concentration channel 1
3	Flow of channel 3	9	Sample concentration channel 2
4	External analog in 1	10	Sample concentration channel 3
5	External analog in 2	11-13	Temperature alarm limits channel 1...3
6	Pressure	14-16	EPC voltage alarm limits channel 1...3

**12.6.11. ETCP: Set TCP/IP Parameters**

Command	Response	Description
_ETCP_K0_zzz.zzz.zzz.zzz _yyy.yyy.yyy.yyy _xxxx	_EDAL_s	zzz= TCP/IP address yyy= TCP/IP subnet mask xxxx= TCP/IP port All changes take effect after next power on cycle

**12.6.12. EH2O Set H<sub>2</sub>O correction parameters**

Command	Response	Description
_EH2O_Km_z.z_y.y_x.x	_EH2O_s	z.z: dry y.y: 1st order coefficient x.x: 2nd order coefficient

**12.6.13. ECO2 Set CO<sub>2</sub> correction parameters**

Command	Response	Description
_EH2O_Km_z.z_y.y_x.x	_EH2O_s	z.z: dry y.y: 1st order coefficient x.x: 2nd order coefficient

**12.6.14. EUDP Set UDP Data streaming parameters**

Command	Response	Description
_EUDP_K0_<UDPPort>_<DataFrequency>_[<Mode>]_[<UDP_IP>]	_EUDP_s	<p>Configure an UDP channel for data streaming of the measuring values via Ethernet UDP. Port: port for open the UDP connection</p> <p>DataFrequency: Frequency for transmit the data in Hz</p> <p>Mode: A: ASCII Mode (optional) UDP_IP: Alternative IP address for open the UDP connection when it should not use the IP of connected TCP/IP client (optional)</p>

**12.6.15. Format of the streaming Data via UDP:**

ASCII Mode:

The measuring values will be sent with ASCII signs. The format is:

<Sequence number>\_x.x\_y.y\_z.z


The sequence number will be incremented with every data packet, which is sent.

The measuring values x.x, y.y and z.z will be sent like in AKON K0 telegram

**12.7. Abbreviations used**

Km	: K' + channel Number (→K1....K3)
Mn	: Measuring range number
M1 .. M4	: Measuring Range 1 .. 4
w.w..zz	: Numerical value
W.W...Z.Z	: Numerical value
T	: Numeric integer value
x	: Number
a0 .. a4	: Polynomial coefficients
s	: Status

# 13. Appendix 1 Electrical Block Diagram

CALIFORNIA ANALYTICAL INSTRUMENTS	
600 - HCID	
MAIN	
Doc Number	110731-1
Rev	1
Sheet	1 of 18

## 14. Appendix 2- Starting With Serial Number U06081

### 1.0 INTRODUCTION

The Model 600 NDIR Series of instruments starting with Serial Number U06081 have several new Hardware and Software features.

The Hardware includes the use of a new memory system, isolation of the analog output signals and 15 relays that are used to buffer the many new digital output signals that are now available

The available digital signals consist of a SERVICE group, that can be used to externally monitor a number of conditions for preventative maintenance and diagnostics.

A second STATUS group, is provided to define the operation of the instrument such as Spanning, Zeroing, Calibrating and the current Range(1, 2, 3, 4, AUTO)

Many of the various signals are duplicated because an NDIR instrument can consist of up to three different channels.

The Software includes modifications to existing functions, changes to the Measurement screen, additional Short-Cut Keys and several New Functions that are listed as follows:

- **MEASUREMENT**

<b>Over Range</b>	<b>888888</b>
<b>Diagnostics</b>	<b>F3</b>
<b>Zero</b>	<b>F5</b>
<b>Span</b>	<b>F6</b>
<b>Standby</b>	<b>F7</b>
<b>Range Limits</b>	<b>F8</b>
<b>Span Values</b>	<b>F9</b>
<b>Outputs</b>	<b>F10</b>

Note: The operator can use these Short-Cut Keys or continue to use existing procedures.

- **NEW FUNCTIONS**

<b>Auto Startup</b>	<b>F5, F7, F7</b>
<b>ALARMS</b>	<b>F5, F7, (Use F6 to toggle ON/OFF)</b>
<b>Offsets&amp; Gains</b>	<b>F4, F3, F5</b>
<b>D/A Calibration</b>	<b>F5, F7, F8</b>
<b>Save Data Archiving Time</b>	<b>F5, F7, F1, F5 (Use ENTER to change recording time)</b>
<b>User Digital Outputs</b>	<b>F5, F9</b>

- **Modifications**

<b>Saved/Not good</b>	<b>F4, F2, F1 or F2 (To flow Zero or Span Gas)</b>
<b>Re-Set Calibration Values</b>	<b>F4, F5</b>


## 2.0 OPERATION OF MEASUREMENT KEYS:

Note: USE the F1 & F2 Keys to view the complete list of menu items.

**2.1 Diagnostics** Use **F3** to toggle between MEASUREMENT an DIAGNOSTIC

**2.2 Zero:** Select the required channel and range then press **F5**.

Note: For instruments with an internal Zero Solenoid select Calibration By Valves. (**F5, F2, F4**)

Zero Gas will be enabled and the observed results can be use to evaluate instrument performance


Press **F3** to toggle to the Diagnostic screen for additional information

Press **F1** to save the value and complete a ZERO calibration for this channel and range

Press **MAIN, F1** to return to the MEASUREMENT screen to select other channels and ranges and repeat the process or **F2** to return to the Manual Calibration screen

**2.3 Span:** Select the required Channel and Range then press **F6**.

Note: For instruments with an internal Span Solenoid select Calibration By Valves (**F5, F2,F4**).

Span Gas for Channel1 and Range 1 will be enabled and the observed results can be used to evaluate instrument performance.


Press **F3** to toggle to the Diagnostic screen for additional information

Press **F1** to save the new value and complete the SPAN calibration for this Channel and Range.

Note that the span gas value used for this channel and range is highlighted and can be changed if necessary. Use the Enter key and the numeric keys as required.

Press **MAIN**, **F1** to return to the MEASUREMENT screen to select other Channels and Ranges and repeat the process or **F2** to return to the Manual Calibration screen.


#### 2.4 Range Limits: (F8);

The Channels and Ranges are factory defined and application specific. Consult California Analytical if any changes are required.

Use the ← → keys to observe the other Channels.

#### 2.5 Span Values: F9


Use the ← → keys to select the desired Channel and the ↑ to select the Range

Note that the span gas value used for this channel and range is highlighted and can be changed if necessary.

Use the **Enter** key and the **NUMERIC** keys as required.

## 2.6 Outputs: F10

The screenshot shows a terminal window titled 'assignment'. It contains a table with two columns: 'Output' and 'Signal'. The 'Signal' column has a black highlight over the word 'Channel' in the first row. At the bottom of the screen, there is a timestamp: 'Thu Jun 28 09:25:18 2007'.

Output	Signal
1	Channel 1
2	Channel 2
3	Channel 3
4	Pressure

Thu Jun 28 09:25:18 2007

Use this screen to define the signals and their location that will be monitored by a remote reordering device.

### 3.0 NEW FUNCTIONS

#### 3.1 Auto Start Up: (F5, F7, F7)

```
Auto Startup |
-----
•Auto Startup      : On
•Wait for [min]   : 1
•Access Level     : 1
•Remote/Manual    : Manual

•Calibrations      :   Ch1  Ch2  Ch3
•Startrange       :   1    1    1
                   :   1    1    1
-----
MAIN SAVE
BACK SAVE
-----
Mon Jan 01 01:40:07 2001 |
```

All key analyzer parameters are stored in a secure memory location and retained when power is removed

In the event of an unexpected power failure it may be desirable to change some parameters until an operator can resume control.

This screen may be used to establish several desirable special instrument start-up parameters that define how the analyzer recovers from loss of AC power

When enabled this screen will define the following:

**Wait:** The time delay in minutes before proceeding  
If **zero** is selected the instrument will not start until warnings are cleared


**Access Level:** The final access level

**Remote/Manual:** The final operating MODE

**Calibrations:** The number of attempts to complete a successful calibration as required in the operator defined Deviation Tables.  
If calibration is not successful the instrument will continue reporting results using the last completed calibration.  
The analyzer can be configured use the previous calibration by selecting zero Calibrations.

**Starting Range:** When all defined actions are completed the analyzer will return to the Measurement Screen and to the range specified.

### 3.2 Alarms On/Off: (F5, F7)

Use **F6** to toggle ON/OFF

The instrument has an extensive library of error messages that will aid in the identification of various anomalous events and are displayed at the bottom of the screen. These messages will assist in Diagnostics and indicating the need for preventative maintenance.

This screen provides an **option** to disable these messages during initial start-up or as may be desired for a particular application.

### 3.3 Offset & Gain: (F4, F3, F5)

Offset and Gain		
Zero/ Span applied correction		
	Offset	Gain
Range1:	0.00	1.00
Range2:	0.00	1.00
Range3:	0.00	1.00
Range4:	0.00	1.00
[CO2]	CO	02
← to select channel		
CAI Tel. (800) 959 0959		

This screen can be used to provide an additional means to display calibration deviations.

Use ← → keys to change channel.

The OFFSET is the value stored during zero calibration.

The GAIN is the value stored during span gas calibration using the operator defined calibration gas.

An increasing or decreasing change to the OFFSET or GAIN when used in conjunction with "Deviations" will provide insight to both short and long term changes to system performance.

### 3.4 Calibrate the Digital Outputs

F5, F7, F8 D/A Calibration

F5, F7, F3 Output Assignment

- **Overview**

The 600 CLD Series of instruments are designed to provide analog outputs that can be configured as 0-10v, 4-20 ma, or 0-20 ma.

With this version the outputs can also be configured to include an additional 1.0 volt and 5.0 volt output and a **calibration** capability.

The outputs can be calibrated to exactly match the results obtained on a PLC, Recorder, Data Logger or other remote recording device that may be connected to the analyzer.

The operator will first select the OUTPUT ASSIGNMENT screen and choose the output that is to be calibrated. All outputs of interest may be selected. When calibration is completed, the operator will return the outputs to their original assignment

The D-A CALIBRATION screen will be then be used to complete the calibration procedure.

This screen provides a section that is used to record the zero signal corrections (zero offset) and a second area to record the 100% signal corrections (Gain) for each of the four output signals that may be defined to develop a voltage or current signal..

Since this is a Digital to Analog conversion, the calibration will require the completion of a simple "trial and error" procedure. The operator will observe the results of a "zero or full scale (Gain) signal generated by the analyzer to the remote recording device and select a correction factor The operator will save this value and then observe the results on attached the remote recording device.

The process of selection and saving for "zero" and "span" will be repeated until a satisfactory calibration is achieved. For 0-1V, 0-5V, 0-10V and a 0-20 ma outputs the Offset and Gain values are independent and do not interact... With the 4-20 ma output, the "Offset (zero)" and "Gain (span)" values interact and may require a few more trials.


The following is a table of typical values:

OUTPUT	OFFSET	GAIN
--------	--------	------

0-20 ma	0.000	0.927
4-20 ma	1.820	0.740
0-1 V	1.300	0.820
0-5 V	1.100	0.820
0-10 V	1.050	0.820

### Procedure

3.4.1 From the Main Menu press **F5,F7,F3**,to obtain following screen:

**F5, F7, F3**

3.4.2 Use the  $\uparrow$  to highlight the outputs that require calibration.

Note: In the above example only Output 1 will be calibrated  
Record the name of these signals, they will be restored.

3.4.3 Press enter to provide access to all the menu of signals that are available. (Real Time, THC, CH<sub>4</sub>, Calibration, Sample pressure, etc.)

3.4.4 Select Calibration and press **ENTER** to complete the selection

Note: Any or all of the four outputs can be selected for calibration  
This screen will not be used again until calibration has been completed.

3.4.5 Press **BACK** to return to the SYSTEM SETUP screen (**F5, F7** from the main menu)

3.4.6 Press **F8** to obtain the following screen:  
Press **Enter** to activate highlighted area

output scaling			☰
Calibrate	D/A	Outputs	
Output	Offset	Gain	mA ↕
1	0.0000	1.2210	
2	1.8000	0.7000	
3	1.8000	0.7000	
4	1.8000	2.7000	
**Set Output-F3-Calibrate**			
MAIN / BACK SAVE			
FS F1: 0% FS: 100%			
CAI Tel. (800) 959 0959			

**F5, F7, F8**

- 3.4.7** Use the ↑ to select the desired output press **ENTER**
- 3.4.8** Press F1 to select a ZERO signal and observe the results on the remote device
- 3.4.9** Change the offset value press **BACK** to save the new value.
- 3.4.10** Press **F8** to return to the D-A Calibration screen and note the results on the remote device.
- 3.4.11** Repeat steps 8.0 thru 10.0 until a satisfactory ZERO calibration is achieved.
- 3.4.12** Complete steps 8.0 thru 10.0 for each of the remaining outputs that require calibration.
- 3.4.13** Press **F5** to produce a full scale (100%) signal
- 3.4.14** Use the arrow keys to position the cursor at the require GAIN value
- 3.4.15** Observe the results on the remote device and make a correction to the GAIN value for the output of interest. Press **BACK** to save this new value
- 3.4.16** Press **F8** to return to the D-A calibration screen
- 3.4.17** Observe the results on the remote device and repeat steps change the GAIN value by repeating steps 3.4.14 thru 3.4.16 as needed for each output.
- 3.4.18** Return to the OUTPUT Assignment screen **F5, F7, F3** from the main menu and change the output signals from CALIBRATE to their original values as defined in step 3.2.2.

### 3.5 Save Data Archiving Time (F5, F7, F1, F5)

Setup RTC		☰
F1	Set time	
F2	Set autocalibration time	
F3	Select calibration range	
F4	Autocalibration on/off	
F5	Save Data Archiving Time	
	Archive Time(secs)	10
F10	Show time	
Archive Interval		
0 = Off F5:SAVE		
CAI Tel. (800) 959 0959		☒

Use ENTER to change recording time

### 3.6 User Digital Outputs

#### Overview

The 600 NDIR Series of instruments have 15 solid state, optically coupled, isolated relays that can be programmed by the operator to indicate the status of numerous signals.

The following is a list of digital signals that are available.

- **SERVICE**

Signal	Displayed Message	MIN	MAX	Definition
F11	Flow 1 ,Check	0.5	3.5	Flow 1
F12	Flow 2 ,Check	0.5	3.5	Flow 2
F13	Flow 3 ,Check	0.5	3.5	Flow 3
E1	Extern Analog 1 ,Check			External Analog 1
E2	Extern Analog 2 ,Check			External Analog 2
P	Pressure ,Check	10	16	Pressure
T	Temperature ,Check	20	50	Temperature
1NC	Channel 1 is not calibrated			Temperature Ch.1
2NC	Channel 2 is not calibrated			Temperature Ch.2
3NC	Channel 3 is not calibrated			Temperature Ch.3
1LoC	Ch1: Low conc. Warning	20	60	EPC Ch.1
2LoC	Ch2: Low conc. Warning	20	60	EPC Ch.2
3LoC	Ch3: Low conc. Warning	20	60	EPC Ch.3
1HiC	Ch1: High conc. Warning	2	8	
2HiC	Ch2: High conc. Warning	2	8	
3HiC	Ch3: High conc. Warning	2	8	
1DT	Ch1: Temperature !			
2DT	Ch2: Temperature !			
3DT	Ch3: Temperature !			
1EV	Ch1: EPC ,Check			
2EV	Ch2: EPC ,Check			
3EV	Ch3: EPC ,Check			
1OR	Ch1: Range overflow			
2OR	Ch2: Range overflow			
3OR	Ch3: Range overflow			
1AU	Ch1: ADC Under Range			
2AU	Ch2: ADC Under Range			
3AU	Ch3: ADC Under Range			
1AO	Ch1: ADC Over Range			
2AO	Ch2: ADC Over Range			
3AO	Ch3: ADC Over Range			
Off				


- **STATUS**

1R2	1 Range 2
1R3	1 Range 3
1R4	1 Range 4
1C	1 In Calibrate
1Z	1 In Zero
1S	1 In Span
1Sa	1 In Sample
2AR	2 Auto Range
2R1	2 Range 1
2R2	2 Range 2
2R3	2 Range 3
2R4	2 Range 4
2C	2 In Calibrate
2Z	2 In Zero
2S	2 In Span
2Sa	2 In Sample
3AR	3 Auto Range
3R1	3 Range 1
3R2	3 Range 2
3R3	3 Range 3
3R4	3 Range 4
3C	3 In Calibrate
3Z	3 In Zero
3S	3 In Span
3Sa	3 In Sample

- **STATUS**

1R2	1 Range 2
1R3	1 Range 3
1R4	1 Range 4
1C	1 In Calibrate
1Z	1 In Zero
1S	1 In Span
1Sa	1 In Sample
2AR	2 Auto Range
2R1	2 Range 1
2R2	2 Range 2
2R3	2 Range 3
2R4	2 Range 4
2C	2 In Calibrate
2Z	2 In Zero
2S	2 In Span
2Sa	2 In Sample
3AR	3 Auto Range
3R1	3 Range 1
3R2	3 Range 2
3R3	3 Range 3
3R4	3 Range 4
3C	3 In Calibrate
3Z	3 In Zero
3S	3 In Span
3Sa	3 In Sample


Typical Relay


These contacts(3, 4) will drive continuously up to 500 MA using a customer voltage supply that does not to exceed 60 VDC.

• **Operation**

- Use F5, F9 to select the first eight outputs
- Use the ↓ to select the desired output
- Press ENTER and use ↑ to select desired item
- NOTE the 600 FID can provide 35 digital output signals



Press **F1** to observe the remaining seven outputs  
 Program as desired per the above

**4.0 Changes to existing functions**


**4.1 Saved or NOT GOOD**

During Manual Calibration the following screens will be displayed to indicate the instruments response to the value of the zero or span gas using the amount that the operator defined in the deviation table.

The above is shown using Zero/SPAN Gas

From Measurement use: **F5 or F6**  
 From Main Menu use: **F4, F2, F1 or F2**


**F5, F2, F3**

Note: This screen is used by the operator to define the maximum acceptable limits of the Zero and Span gas for both Manual and Automatic Calibrating.


#### 4.2 Reset Calibration Values

When the re-set calibrations value function is used all recorded deviations

will be set to zero


**F4, F5**

**F4 F3, F4**

(Used to observe Auto Cal Results)

**F4, F3, F2**

(Used to observe Manual Cal results)

The above are the new deviations after the operator elects to re-set the calibration values

# MELISSA


**TECHNICAL NOTE**

**96.1**



## ***APPENDIX 13***

This document is confidential property of the MELISSA partners and shall not be used, duplicated, modified or transmitted without their authorization

Memorandum of Understanding ESTEC 4000 100 293/10/NL/PA

	A	B	C	D	E	F	G	H	I	J	K	L	M	Kg
<b>Tecnoplus</b>	439	108	261	200	216	241	253,8	88	9	243	196	G1"	G1"	10,5

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**Temperatura ambiente**  
de -10° C a 50° C

**Maximum operating pressure**  
6 bar

**Inlet charge**  
Maximum 0,5 bars lower than the adjusted starting pressure of the set.

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The pump can operate at any point inside the indicated area. The characteristic curves depend on the delivery pressure.

By way of example, the curves are shown for delivery pressures of 150, 200 and 300 kPa. The operating limit curve corresponds to the maximum rotating speed.

230 V 50 Hz	230/400 V 50 Hz	A			P1 (kW)		kW	HP	µF
		1~ 230 V	230 V	3~ 400 V	1~	3~			
<b>Tecnoplus</b>	-	3.6	-	-	0.75	-	0.55	0.75	12

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