



# MELiSSA Pilot Plant



Document Identification : COO3 – WP94.6 – Test Protocol for ramp-up of compartment 1	Type	Chrono	Issue	Page : 1 / 17
	TN	94.61	0	



Universitat Autònoma  
de Barcelona

## *TECHNICAL NOTE 94.61*

### **Call Off Order 3 – COMPARTMENT I Additional Characterization**

#### **Work Package 94.6**

#### **Test Protocol – Ramp-up phase**

Prepared by/Préparé par      A. Fossen, E. Peiro, F. Godia  
Reference/Référence          MELiSSA Pilot Plant Frame Contract 19445/05/NL/CP  
Issue/Edition                  0  
Revision/Révision              0  
Date of issue/Date d'édition   8/2/2011  
Status/Statut                  Final



# MELiSSA Pilot Plant



Document Identification : COO3 – WP94.6 – Test Protocol for ramp-up of compartment 1	Type	Chrono	Issue	Page : 2 / 17
	TN	94.61	0	

## APPROVAL

Title Titre	Test Protocol for ramp up of compartment 1	Issue Edition	0	Revision Révision	0
----------------	--	------------------	---	----------------------	---

Prepared by Auteur	A. Fossen	Date Date	25/01/11
-----------------------	-----------	--------------	----------

Checked by Verifié par	E. Peiro	Date Date	06/02/2011
---------------------------	----------	--------------	------------

Approved by Approuvé par	F.Gòdia	Date Date	08/02/2011
-----------------------------	---------	--------------	------------

Approved by customer Approuvé par le client	B. Lamaze	Date Date	8/2/2011
--	-----------	--------------	----------

## CHANGE LOG

Issue/Edition	Revision/Révision	Status/Statut	Date/Date
0	0	Final	8/2/11

## Distribution List

Name/Nom	Company/Société	Quantity/Quantité
Brigitte LAMAZE	ESA	2 hardcopies + electronic version
Technical Team	MPP	1 copy in QP 0007
Overall Manager	MPP	1 copy in QC1119
Quality Manager	MPP	The original signed copy in QC1004



# MELiSSA Pilot Plant



Document Identification : COO3 – WP94.6 – Test Protocol for ramp-up of compartment 1	Type	Chrono	Issue	Page : 3 / 17
	TN	94.61	0	

## Table of contents

1. SCOPE .....	4
2. APPENDIX : MPP-OP-10-1004(0), 13 PAGES.....	4



Document Identification : COO3 – WP94.6 – Test Protocol for ramp-up of compartment 1	Type	Chrono	Issue	Page : 4 / 17
	TN	94.61	0	

## 1. Scope

In the frame of Call Off Order 3, the objective of this document is to summarize the protocol to follow in order to grow the inoculum of compartment 1 up to 100L in the bioreactor and reach the adequate dry weight of 55 g/L to 75 g/L before starting the 10 days hydraulic residence time.

As this growth had to be carried out after the loss of the C1 inoculum, and that an MPP operating procedure was rewritten for that same purpose, it was agreed with ESA to only provide said operating procedure MPP-OP-10-1004 in appendix of the present protocol, for implementation.

The recovery of the lost inoculum is envisaged in two steps : first in the Influent Tank that will be prepared for the semi continuous culture of the recovered inoculum, and once the maximum volume is reached in the Influent Tank, a second semi continuous culture will be started in C1 Bioreactor.

## 2. Appendix :

MPP-OP-10-1004(0), 13 pages



# MELiSSA Pilot Plant








Universitat Autònoma de Barcelona

Document Identification : C1 Inoculum recovery in Influent Tank Semi-continuous Operating Procedure	Type	Chrono	Issue	Page : 1 / 13
	MPP-OP	10-1004	(0)	

## C1 Inoculum recovery in Influent Tank Semi-continuous Operating Procedure

Approval Loop :

Date	Issue	Prepared by (visa):	Checked by (visa):	Approved by (visa):
24/11/2010	0	V.GARCIA Analysis technician  R.MOYANO Maintenance technician  E.PEIRO Technical Manager 	E.PEIRO  Technical Manager	A. FOSSEN  Quality Manager

Distribution List of controlled copies :

Company	Name or position	Number of copies / location
MPP	MPP Technical Team	1 copy in QP 0007
	MPP Overall Manager	1 copy in QC 1119
	MPP Quality Manager	the original signed copy

This document is confidential property of the MELISSA partners and shall not be used, duplicated, modified or transmitted without their prior authorization

Memorandum of Understanding ESTEC 4000 100 293/10/NL/PA



# MELISSA Pilot Plant

Document Identification : C1 Inoculum recovery in Influent Tank Semi-continuous Operating Procedure	Type	Chrono	Issue	Page : 2 / 13
	MPP-OP	10-1004	(0)	

### Change log :

Date	Issue	Reason of the change	Modified paragraphs
24/11/2010	0	Creation	

### Table of content:

1. Scope.....	3
2. Reference documents .....	3
3. General description of C1 inoculum maintenance process .....	4
4. Personal protective equipment .....	6
5. Reagents .....	6
6. Equipment.....	6
7. Operating procedure .....	6
8. Labelling .....	7
9. Disposal of C1 reactor content, samples and supernatant.....	8
10. Follow up (MPP-REC-10-1008) filling out instructions .....	9
11. Liquid Samples handling .....	9
12. Gas analysis .....	9
13. Liquid samples analysis .....	9
14. CI Software modification for Broth startup in the Influent.....	11
15. Process Diagram.....	13





# MELISSA Pilot Plant

Document Identification : C1 Inoculum recovery in Influent Tank Semi-continuous Operating Procedure	Type	Chrono	Issue	Page : 3 / 13
	MPP-OP	10-1004	(0)	

## 1. Scope

This procedure describes the semi-continuous operating procedure of the C1 inoculum to be cultured in the Influent Tank during the recovery and growth start-up phase previous to the transfer to the bioreactor, after the loss of the bioreactor content on 18<sup>th</sup> Oct 2010 (MPP-NCR-10-1006).

As previously agreed with UBP and ESA (see MPP-MOM-10-1003(0)-AF-20101119), the strategy for the recovery and growth start-up period in order to have the CI inoculum ready for the CI characterisation phase will be as follows:

Step 1 : Restarting 50 L reactor with the available volume of back-up Inoculum (maintained at 4°C)

1. Transfer the inoculum with initial feed (50:50) in the influent tank (working vol. 50L)
2. Feed with only beet (which is easily degradable), until reaching a nominal biomass concentration (close to the one expected for starting continuous operation of the reactor).
3. Progressively, feed with the mix of the 3 higher plants. The frequency / quantity should be adapted to the response of the reactor (probably following the CO<sub>2</sub> gas production is sufficient to know the status of the reactor). The objective is to fit a weekly feeding quantity as it is foreseen in the test plan.
4. When the objective of the nominal operation of the reactor with the feeding protocol expected in the test plan is achieved, start the 100L reactor.

Step 2 : Restarting the 100 L reactor

1. Transfer the culture broth from the influent tank into the 100 L bioreactor
2. Adapt the ramp in the frequency / quantity of feed (plant composition identical to the one used in the influent tank) as function of the response of the reactor knowing that the objective is to reach the Nominal 10 days RT as depicted in the test plan.

The procedure includes a general description of the process, the follow up record sheets (MPP-REC-10-1008) filling out instructions and the analyses to be performed.

## 2. Reference documents

- |     |  |   |
|-----|--|---|
| RD1 | CI Software Description and Procedures, September 2010 | SHERPA ENGINEERING, CI Software Description and Procedures. |
|-----|--|---|



# MELiSSA Pilot Plant

UAB

Universitat Autònoma  
de Barcelona

Document Identification : C1 Inoculum recovery in Influent Tank Semi-continuous Operating Procedure	Type	Chrono	Issue	Page : 4 / 13
	MPP-OP	10-1004	(0)	

RD2 NTE-MCI-HB-012, NTE, CI HMI SW User Manual for the MELISSA CS CI.  
Issue 1.0, May 2009.

### **3. General description of C1 inoculum maintenance process**

C1 backup inoculum is maintained in the influent tank of pilot scale C1 located in C1 room. It is a 60L volume tank with 50L working volume. Some hardware modifications were done on the influent tank to adapt it to C1 culture. In process diagram MPP-PID-10-1002-A1, hardware modifications are shown in different colour. This diagram is shown below in paragraph "15. Process Diagram."

At present, influent tank is equipped with mechanical agitation and a recirculation loop that it is used to mix the content and to introduce the feeding by the "getting cake process" ([RD1], paragraph 3.7, Procedure 7).

Temperature is controlled by means of the Temperature Control loop and dedicated hardware of the bioreactor, but connected to the influent tank jacket. For this purpose, the following changes were performed:

- The thermal control piping connected to the bioreactor jacket was mounted on the influent jacket

The pH is also controlled by the pH Control loop of the bioreactor, by means of one of its pH probes connected to a pH transmitter/controller (both Mettler Toledo). pH is maintained at 5.3 ( $\pm 0.2$ ) by the addition of acid (HCl, 3M) and base (NaOH, 3M) by the two peristaltic pumps installed in the bioreactor skid. For this purpose, the following changes were performed:

- The bioreactor pH probe was mounted on the influent tank directly, and the acid/base addition piping was mounted on the top of the influent tank.
- The influent blender was linked to the pH addition control, to be sure that the blender is ON when is needed to add acid/base.

The tank headspace pressure is maintained at 120 mBar (to avoid the shut down of recirculation pump, which stops above 135 mbar, approximately) by the addition of a nitrogen flow of 100mL/min by a mass flow controller calibrated with Air with a working range of 20 – 1000 mL/min (Bronkhorst Hi-Tech).





# MELISSA Pilot Plant

**UAB**

Universitat Autònoma  
de Barcelona

Document Identification : C1Inoculum recovery in Influent Tank Semi-continuous Operating Procedure	Type	Chrono	Issue	Page : 5 / 13
	MPP-OP	10-1004	(0)	

Gas analysis: the piping for the analysis gas loop was mounted in the free ports of the influent tank, to have the values of CO<sub>2</sub> and Methane. This point does not affect to the PLC program.

In the Human-Machine Interface (HMI) located in the control room of the MPP, the temperature, pH and pressure control, as well as the agitation and recirculation pump, can be actuated; other parameter values can be visually checked [RD2].

Further details about the software control configuration are given in the Section 14 of this document.

C1inoculum is maintained in semi-continuous mode. The tank is fed once per week with a solid mixture of maintenance feed composed only by beet in the first stage (MPP-OP-10-1005) and then progressively incorporating the complete maintenance feed (MPP-OP-09-1001) Culture nominal volume is 50 L. Solids are returned to the reactor by a centrifugation step and separation of supernatant. Residence time is 20 days and the corresponding load of substrate is 0,58 gDW/(day.L).



# MELISSA Pilot Plant

Document Identification : C1 Inoculum recovery in Influent Tank Semi-continuous Operating Procedure	Type	Chrono	Issue	Page : 6 / 13
	MPP-OP	10-1004	(0)	

## **4. Personal protective equipment**

- 4.1. Lab coat
- 4.2. Gloves
- 4.3. Goggles
- 4.4. Full face shield
- 4.5. Half-face organic vapour mask

## **5. Reagents**

- 5.1. Bleach

## **6. Equipment**

- 6.1. C1 Influent Tank of 55L working volume
- 6.2. Steel funnel (the funnel connected to valve H3V\_1001\_01 can be used)
- 6.3. Disposal bottle of 10L volume
- 6.4. Graduated plastic bucket of 20L volume
- 6.5. Plastic rod
- 6.6. Infrared gas analyser GA94 (Geotechnical Instruments)

## **7. Operating procedure**

### **7.1. Inoculum back up and feeding 1 day procedure**

#### 7.1.1. On Fridays:

- Defrost one C1 beet feed bag of 2,53-2,57Kg (if only beet is fed, see MPP-OP-10-1005) or 4,1Kg (if the whole vegetable mix is used, see MPP-OP-09-1001) storing it at 3 – 5 ° C for next Monday feeding operation.

#### 7.1.2. On Mondays:

1. Dispose a sample of 80 – 100 mL of reactor content. Take another sample of 80 – 100 mL of reactor content. Label the sample as indicated below and name it BEFORE FEEDING (see 8).



# MELISSA Pilot Plant



Universitat Autònoma de Barcelona

Document Identification : C1Inoculum recovery in Influent Tank Semi-continuous Operating Procedure	Type	Chrono	Issue	Page : 7 / 13
	MPP-OP	10-1004	(0)	

2. Take a sample of 9L of reactor content in order to centrifugate the sample and keep the pellet.
3. Take a sample of supernatant of 80 – 100 mL volume and label it as indicated below. Name it FILTRATE (see 8). Dispose the supernatant as indicated in “9. Disposal of C1 reactor content, samples and supernatant.”
4. In case that only beet is used, check the particle size of the defrosted feed bag of 2.53-2.57 Kg of red beet, by sieving it through a 2 mm mesh sieve. Use only feed with particle size lower than this size to feed the CI culture.
5. Mix in a graduated vessel the defrosted 2.53-2.57 Kg feed plus the pellet obtained in the centrifugation with distilled water, until a final volume of 9 L ( $\pm 0.2$ ).
6. Transfer the feed mixed in the WPU in a graduated plastic bucket of 20L of volume
7. Take a sample of feed of 80 – 100 mL volume and label it as indicated below. Name it FEED (see 8).
8. Open one port of the top of the influent tank. Use the full face visor until step 11.
9. Introduce the content of feed into the influent tank by means of the graduated plastic bucket in three times. Fill the plastic bucket with a maximum of 6L approx. of the content of the bottle each time.
10. Close the top lid of the influent tank.
11. After 1 hour, dispose a sample of 80 – 100 mL of reactor content. Take another sample of 80 – 100 mL of reactor content. Label the sample as indicated below and name it AFTER FEEDING.
12. Fill out the follow up record sheet MPP-REC-10-1008.

## 8. Labelling

### 8.1. Samples labelling

- Use the DYMO or a pen to label the samples as indicated below:

MPP C1	{	RC BF	(REACTOR CONTENT BEFORE FEEDING)	DATE	(YY/MM/DD)
		RC AF	(REACTOR CONTENT AFTER FEEDING)		
		FEED # ...			
		FILTRATE			

This document is confidential property of the MELISSA partners and shall not be used, duplicated, modified or transmitted without their prior authorization





# MELISSA Pilot Plant

Document Identification : C1Inoculum recovery in Influent Tank Semi-continuous Operating Procedure	Type	Chrono	Issue	Page : 8 / 13
	MPP-OP	10-1004	(0)	

- Examples:

MPP C1 RC BF 090126
MPP C1 RC AF 090126
MPP C1 FEED #12 090126
MPP C1 FILTRATE 090126

## 8.2. Inoculum back ups labelling

MELISSA C1 INOCULUM BACK UP + DATE (YY/MM/DD)

## 9. Disposal of C1 reactor content, samples and supernatant

To dispose big volumes (>100 mL) use the organic vapour mask.

9.1. Add to the disposal volume a quantity equal to 10% of the disposal volume.

- For example: To dispose 9L of waste use 1L of bleach.

9.2. Close the disposal vessel, shake it gently until complete mixing and throw the content to the liquid waste collector of the MPP.



# MELISSA Pilot Plant

**UAB**

Universitat Autònoma  
de Barcelona

Document Identification : C1 Inoculum recovery in Influent Tank Semi-continuous Operating Procedure	Type	Chrono	Issue	Page : 9 / 13
	MPP-OP	10-1004	(0)	

## **10. Follow up (MPP-REC-10-1008) filling out instructions**

10.1. Fill out the record sheet MPP-REC-10-1008 every day and when it is indicated in this procedure. Fill it out also after taking any action (extra sampling, hardware modification or fixing) or when some parameter is out of the range (although any action is done).

## **11. Liquid Samples handling**

11.1. Analyse pH of all samples at the moment of taking them and note down the result in the follow up record sheet MPP-REC-09-1003.

11.2. Store the samples at 3 – 5 °C to be further analysed as it is indicated in paragraph “13. Liquid samples analysis”.

## **12. Gas analysis**

12.1. Analyse influent tank headspace gas with the infrared gas analyser GA94 (Geotechnical Instruments) once per week, just before feeding.

12.2. Connect the gas analyser inlet to the valve HV\_1052\_02. (If the gas sampling tubing presents condensates, clean it before the measurement).

12.3. Connect to the gas analyser outlet the H<sub>2</sub> probe and connect the probe outlet to the valve HV\_1052\_03.

12.4. Open the valves HV\_1052\_02 and HV\_1052\_03 and start the gas analyser.

12.5. With the analyser pump in OFF mode wait 5 – 10 minutes (until a stable measure) and note down the results in the follow up record sheet (MPP-REC-10-1008)

12.6. Repeat the measurement with the H<sub>2</sub>S probe and note down the H<sub>2</sub>S result in the same follow up sheet.

12.7. Close the valves HV\_1052\_02 and HV\_1052\_03 and disconnect the gas analyser.

12.8. Clean the H<sub>2</sub>S and H<sub>2</sub> probes: Connect one probe and start the analyser with the pump in ON mode in a clean atmosphere during 5 – 10 minutes. Repeat the process with the other probe.

## **13. Liquid samples analysis**

### **13.1. Dr Lange analysis**

This document is confidential property of the MELISSA partners and shall not be used, duplicated, modified or transmitted without their prior authorization





# MELISSA Pilot Plant

Document Identification : C1 Inoculum recovery in Influent Tank Semi-continuous Operating Procedure	Type	Chrono	Issue	Page : 10 / 13
	MPP-OP	10-1004	(0)	

- Every week analyse Feed, Reactor Content and Filtrate samples (as it is shown in Table 1) by ammonium (MPP-QCP-07-00011), total nitrogen (MPP-QCP-07-00020), COD total (MPP-QCP-07-00005) and COD soluble (MPP-QCP-07-00004) using Dr. Lange kits.
- Record the results in MPP-REC-09-1003.

Table 1. Analyses performed on C1 samples

Sample Type	Analysis parameter			
	Total fraction	Soluble fraction		
	COD	COD	NH <sub>4</sub> <sup>+</sup>	N-total
Feed	X	X	X	X
Reactor Content	X	X	X	X
Filtrate	X			

### 13.2. Dry weight analysis

- Every week analyse the dry weight (MPP-QCP-07-00007), ashes (MPP-QCP-07-00007) and electroconductivity (MPP-QCP-07-00008) of Feed, Reactor Content and Filtrate samples.
- Record the results in the sheet MPP-REC-09-1005 and MPP-REC-09-1003.

### 13.3. VFA analysis

- Every week, when analysing the parameters by Dr. Lange kits, store 1mL of each sample soluble fraction at -78/-80°C in GC vials previously labelled for VFA analysis to be performed later on.

### 13.4. Cell count

- Every two weeks, perform the aerobic and anaerobic count of Reactor content, and perform the same count of each batch of Feed (MPP-QCP-09-1001).



# MELiSSA Pilot Plant

Document Identification : C1Inoculum recovery in Influent Tank Semi-continuous Operating Procedure	Type	Chrono	Issue	Page : 11 / 13
	MPP-OP	10-1004	(0)	

## 13.5. Elemental composition and Minerals analysis

- Every 1 month take and freeze-dry 100 – 150 mL of Reactor content (MPP-QCP-07-0010).
- Following the annual calendar previously agreed, send the Feed (MPP-OP-09-1001) and Reactor content freeze-dried samples to the “Servei d’Anàlisi Química” (SAQ) of the UAB, located in the Biosciences Faculty.
- Fill out one application for all the samples for minerals analysis. Fill out one application for each sample for elemental composition analysis.

## **14. CI Software modification for C1 inoculum recovery in the Influent**

The following list covers the modifications performed by SHERPA on CI running software version (V01\_01\_CI) in order to provide the adequate conditions for the process control already described in the Section 3 of this document. The temporary version is called “Temp\_CI”.

- 1- In order to have the control of the temperature linked to the bioreactor heater, the PLC address TT\_1008\_01 is now connected to the TT\_1002\_01 PLC input (PLC@:300014). The sensor failure tag (TT\_1008\_01\_ERR) is also connected to the TT\_1002\_01\_ERR output.
- 2- TT\_1002\_01 is connected to the TT\_1008\_01 PLC input (PLC@:300015). Alarm previously linked to the bioreactor level is now implemented on the Influent level. In case of Very Low Level in the influent tank, the pH, the blender and the temperature control loops are triggered to OFF mode (the very low level threshold of influent tank is 10 litres).
- 3- As only one pH probe is used for the control inside the Influent (AT\_1011\_02), the automatic probe switches in case of failure are removed. If AT\_1011\_02 fails, the pH control Loop is triggered to OFF mode.
- 4- If the pH is in auto mode when the Influent Blender stops, the tag “CL1011\_Blender\_Stopped\_A” prevents any injection. It ensures that Base and Acid injections are not done if the broth is not mixed.



# MELISSA Pilot Plant

Document Identification : C1Inoculum recovery in Influent Tank Semi-continuous Operating Procedure	Type	Chrono	Issue	Page : 12 / 13
	MPP-OP	10-1004	(0)	

- 5- The valve SV\_1003\_01 is opened when the pressure reaches -50 mbar (very low threshold) and 120 mbar (Very high threshold).
- 6- The pH Very Low threshold alarm (which triggers to OFF mode the pH control Loop) is changed from 5 to 4.




# MELISSA Pilot Plant



Document Identification :		Type	Chrono	Issue
C1 Inoculum recovery in Influent Tank		MPP-OP	10-1004	(0)
Semi-continuous Operating Procedure				Page : 13 / 13

## 15. Process Diagram

Process diagram file:

**MPP-PID-10-1002-A1**



A1	RM	17/1/2010	Checked by	Approval date	Approval by	Reason for change
	Author	Date				
 <b>J.T.C. Consultant</b> <small>of the Institut de Recerca i Innovació Tecnològica (IRIT)</small>						
						
Your document is confidential property of the MELISSA partners and shall not be used, duplicated, modified or transmitted without their prior authorization.						
Adaptation of Influent tank to be the Bioreactor Adaptation of Influent tank to be the Bioreactor						
Scale: 1/1 Project Title: C1 Inoculum recovery in Influent Tank Document Code: MPP-PID-10-1002-A1						
Provisional CI						

