Document Identification : COO3 – WP94.2 – Sampling and Analysis Protocols Issue 2 - Test with 10 days hydraulic residence time	Type	Chrono	Issue	Page : 1 / 47
	TN	94.22	0	

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TECHNICAL NOTE 94.22

Call Off Order 3 – COMPARTMENT I Additional Characterization

Work Package 94.2

Sampling and Analysis Protocols – Issue 2 - Test with 10 days hydraulic residence time

Prepared by/Préparé par	A. Fossen, E. Peiro, F. Godia
Reference/Référence	MELiSSA Pilot Plant Frame Contract 19445/05/NL/CP
Issue/Edition	0
Revision/Révision	0
Date of issue/Date d'édition	13/10/11
Status/Statut	Final

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
Document Identification : COO3 – WP94.2 – Sampling and Analysis Protocols Issue 2 - Test with 10 days hydraulic residence time	Type TN	Chrono 94.22	Issue 0	Page : 2 / 47
--	------------	-----------------	------------	---------------

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Distribution List

Name/Nom	Company/Société	Quantity/Quantité
Brigitte LAMAZE	ESA	2 hardcopies + electronic version
Technical Team	MPP	1 copy in QP 0007
Overall Manager	MPP	1 copy in QC1119
Quality Manager	MPP	The original signed copy in QC1004

MELiSSA Pilot Plant

Document Identification : COO3 – WP94.2 – Sampling and Analysis Protocols Issue 2 - Test with 10 days hydraulic residence time	Type	Chrono	Issue	Page : 3 / 47
	TN	94.22	0	

Table of contents

1. INTRODUCTION	5
2. APPLICABLE AND REFERENCE DOCUMENTS	5
2.1. Applicable documents.....	5
2.....	5
2.2. Reference documents.....	7
3. ACRONYMS AND DEFINITIONS.....	7
4. RESPONSIBILITIES.....	7
5. TEST ITEMS	8
5.1. Description (PID, technical drawings, user manual).....	8
5.2. Hazards induced by test item and safety measures to be taken	8
5.3. Instructions for operation.....	8
5.4. Instructions for maintenance.....	8
6. TEST LOGIC.....	9
6.1. Objectives of the tests	9
6.2. Approach followed	9
6.3. Applicable requirements	10
6.4. Recall of test sequence.....	10
6.5. Success/failure criteria for the 10 days HRT.....	10
7. DATA COLLECTION PLAN – SAMPLING PLAN	11
7.1. Uncertainty acceptance level	11
7.2. Measurement plan.....	11
7.3. Sampling techniques	12
7.4. Sampling ports location	12
7.5. Sampling recipients.....	16
7.6. Use of the volumes extracted during the bleeding.....	16
7.7. Conditioning of the samples	16
7.8. Labelling of samples.....	18
7.9. Conservation of the samples :.....	18
7.10. Sampling and analyses frequency on the liquid-solid phase.....	19
7.11. Sampling and Analysis frequency on the gas phase	26
7.12. Analyses	26
8. RESOURCES SPECIFICATION FOR THE TESTS.....	28
8.1. Personnel: staff qualification and training needs.....	28
8.2. Hardware: instruments, specific part, hardware for software operation.....	28
8.3. Software : verification of software, backup needs.....	28
8.4. Facilities : environmental needs, test conditions, interfaces needs, utilities needs	28

MELiSSA Pilot Plant

Document Identification : COO3 – WP94.2 – Sampling and Analysis Protocols Issue 2 - Test with 10 days hydraulic residence time	Type	Chrono	Issue	Page : 4 / 47
	TN	94.22	0	

8.5. Test conditions.....	28
9. MEASUREMENT AND DATA SAMPLING	28
9.1. Data logfile	28
9.2. Reporting of status for a test.....	28
9.3. Deviations and non conformances	29
9.4. Record for the performance of samplings and analysis.....	29
10. APPENDIX 1 – APPROVED MPP PROTOCOLS FOR ANALYSES ON SAMPLES	30
11. APPENDIX 2 - EXTERNAL PROTOCOLS IMPLEMENTED FOR ANALYSES	31
11.1. Inorganic elemental composition (minerals) SAQEAt0001_00	31
11.2. C,H,N,S organic elemental composition : SAQEAE0001_00_CHNS	31
11.3. Fibers content : not available in English, added in Catalan language.....	32
12. APPENDIX 3 – EXAMPLE OF MPP-REC-10-1001 FOR THE FOLLOW-UP OF THE COMPARTMENT 1	39
13. APPENDIX 4 – EXAMPLE OF RECORD FOR DRY WEIGHT AND ASHES : MPP-REC-09- 1005 40	
14. APPENDIX 5 – EXAMPLE OF RECORD FOR ELECTROCONDUCTIVITY, NH4, N TOTAL, COD TOTAL, COD SOLUBLE : MPP-REC-09-1004 AND MPP-REC-09-1003.....	41
15. APPENDIX 6 – EXAMPLE OF ANALYSIS REPORT FOR THE CHNS ANALYSIS	43
16. APPENDIX 7 – EXAMPLE OF REPORT FOR MINERALS ANALYSIS	44
17. APPENDIX 8 – EXAMPLE OF THE DATABASE USED FOR C1 RESULTS OF ANALYSIS AND OPERATIONAL PARAMETERS	46

Document Identification : COO3 – WP94.2 – Sampling and Analysis Protocols Issue 2 - Test with 10 days hydraulic residence time	Type	Chrono	Issue	Page : 5 / 47
	TN	94.22	0	

1. Introduction

The objective of this document is to summarize the sampling and analysis strategy to be implemented on Compartment 1 for the test at 10 days HRT.

This document is defining a first set of samplings and analyses that might evolve depending on the feedback of the C1 bioreactor operation and on the results of said samplings and analyses.

2. Applicable and reference documents

2.1. Applicable documents

Ref.	Title	Reference	Issue	Date
AD1	MPP Proposal for Call Off Order 3 – C1 additional characterization	OFR-ESA-03/07-UAB	1	30/11/07
AD2	MPP Quality Manual	MPP-QA-07-0001	2	
AD3	MPP Rules for Good Laboratory Practices	MPP-QA-07-0003	0	
AD4	Test Plan for C1 additional characterization tests	TN94.5	0	
AD5	PID of Compartment 1	MPP-PID-10-1001	B3	5/10/2011
AD6	MPP Operation Manual for C1	MPP-OP-12-1001	0	9/2/2012
AD7	C1 Acceptance Review Datapackage including HMI and PLC software user manuals	DP94.1	1	October 11
AD8	MPP Maintenance Manual for C1	MPP-UM-11-1001	0	
AD9	Test Protocol for nominal operation with 10 days HRT	TN94.62		
AD10	Analytical Procedure for Ashes determination	MPP-QCP-07-0001	(1)	16/10/2007
AD11	Analytical Procedure for COD soluble determination	MPP-QCP-07-0004	(1)	16/10/2007
AD12	Analytical Procedure for COD total determination	MPP-QCP-07-0005	(1)	16/10/2007
AD13	Analytical Procedure for Dry matter determination	MPP-QCP-07-0007	(1)	16/10/2007

MELISSA Pilot Plant

Document Identification : COO3 – WP94.2 – Sampling and Analysis Protocols Issue 2 - Test with 10 days hydraulic residence time	Type	Chrono	Issue	Page : 6 / 47
	TN	94.22	0	

AD14	Analytical Procedure for EC determination	MPP-QCP-07-0008	(1)	16/10/2007
AD15	Analytical Procedure for lyophilization	MPP-QCP-07-0010	(1)	16/08/2007
AD16	Analytical Procedure for NH4-N determination	MPP-QCP-07-0011	(1)	16/10/2007
AD17	Analytical Procedure for NO2-N determination	MPP-QCP-07-0012	(1)	16/10/2007
AD18	Analytical Procedure for NO3-N determination	MPP-QCP-07-0013	(1)	16/10/2007
AD19	Analytical Procedure for O.D. determination	MPP-QCP-07-0014	(1)	16/10/2007
AD20	Analytical Procedure for pH determination	MPP-QCP-07-0015	(1)	16/10/2007
AD21	Analytical Procedure for Protein determination	MPP-QCP-07-0018	(1)	16/10/2007
AD22	Analytical Procedure for Total nitrogen determination	MPP-QCP-07-0020	(1)	16/10/2007
AD23	Quality control procedure for counting bacteria from C1	MPP-QCP-09-1001 0	0	20/10/2009
AD24	Measurement procedure for pH	MPP-QCP-10-1001	(0)	22/02/2010
AD25	Measurement procedure for gas mass flow measurement in C1	MPP-QCP-10-1002	(0)	22/02/2010
AD26	Quality control procedure for Alkalinity determination	MPP-QCP-11-0001	(0)	06/02/2012
AD27	Quality control procedure for the measurement of VFAs by HPLC	MPP-QCP-11-0002	(0)	06/02/2012
AD28	Analytical Procedure for Protein determination from C1	MPP-QCP-12-1001	(0)	01/03/2012
AD29	Procedure for the control of sterility of C1 filtrate line and effluent tank	MPP-QCP-12-1002	(0)	18/6/12
AD30	Quality Control Procedure for total nitrogen determination (total broth)	MPP-QCP-12-0001	(0)	20/6/12

Document Identification : COO3 – WP94.2 – Sampling and Analysis Protocols Issue 2 - Test with 10 days hydraulic residence time	Type	Chrono	Issue	Page : 7 / 47
	TN	94.22	0	

2.2. Reference documents

Ref.	Title	Reference	Issue	Date
RD1	TN 94.11 Compartment I Integration in MPP	TN 94.11	0	13.02.09
RD2	HAZOP on Compartment 1	MPP-TN-08-1001	0	01/09/2008
RD3	Gas Chromatograph User Manual	MPP-UM-09-0009	1	23.10.06
RD4	Portable Gas Analyzer User Manual	MPP-UM-09-0012	0	
RD5	TN 83.7 Expertise of level 0 control loops on the 100 L pilot reactor	TN 83.7	1	23.10.06
RD6	Minutes of meeting MPP/UBP on C1 characterization	MPP-MOM-08-1007	0	16.04.2008
RD7	EPAS EWC User Manual	User Manual	1	12.06.07
RD8	HPLC Ultimate 3000 User Manual	MPP-UM-11-0002	0	25.11.2011

3. Acronyms and definitions

BR: bioreactor
 CI : compartment I
 CST : capillary suction time
 FU: Filtration unit
 GL: Gas loop
 HMI: human interface
 HPLC High Pressure Liquid Chromatography
 HRT: hydraulic residence time, equivalent to liquid residence time
 ICP-MS : Induced Coupled Plasma Mass Spectrometry
 MELiSSA: Micro-Ecological Life Support System Alternative
 SFC: Sequential function chart
 TAR test acceptance review
 TRR test readiness review
 UAB: Universitat Autònoma de Barcelona
 VFA: volatile fatty acids

4. Responsibilities

For the whole call-off order 3, concerning the characterization tests on C1, the responsibilities of the different organizations are detailed in AD4.

Document Identification : COO3 – WP94.2 – Sampling and Analysis Protocols Issue 2 - Test with 10 days hydraulic residence time	Type	Chrono	Issue	Page : 8 / 47
	TN	94.22	0	

5. Test items

5.1. Description (PID, technical drawings, user manual)

The compartment 1 was delivered in the MPP and installed as described in RD1.

It consists of 3 subunits or modules that are described on the PID AD5 and in the User Manual AD6, namely :

- The bioreactor and influent tank skid
- The gas loop skid
- The filtration unit skid

The system is operated automatically from a programmable logical controller (PLC) as described in AD7.

5.2. Hazards induced by test item and safety measures to be taken

As explained in the hazard and operability study carried out on compartment 1 (cf. RD2), the main hazards induced by the operation of compartment 1 are:

- pressure (gas: up to 3 barg, liquid: up to 5 barg)
- temperature (steam sterilization)
- chemical (acid/base for pH control)
- biological (biohazard level 2 as a maximum when using faeces for the feeding of C1)
- flammable gases (H₂, CH₄) ;

The adequate individual protection measures shall be taken by the operators in order to limit the exposure to these hazards. As detailed in AD4, these measures include :

- wearing of a labcoat
- wearing of safety goggles
- wearing of gloves when manipulating materials or equipments
- respect of the user and maintenance instructions, in particular the respect of the confined and anaerobic conditions in the bioreactor

5.3. Instructions for operation

See AD6 and AD7, and RD7

5.4. Instructions for maintenance

See AD6 and AD7, AD8 and RD7

Document Identification : COO3 – WP94.2 – Sampling and Analysis Protocols Issue 2 - Test with 10 days hydraulic residence time	Type	Chrono	Issue	Page : 9 / 47
	TN	94.22	0	

6. Test logic

6.1. Objectives of the tests

The objective of the C1 characterization tests is to collect as many data as possible on steady state periods in order to provide the parameters necessary for the understanding of the C1 process behaviour and for the construction of a knowledge model.

For phase 3, the objective is to reach the steady state operating at 10 days HRT and collect all the main parameters characterizing the process at that hydraulic residence time. These parameters will then allow to determine if the success criteria fixed in AD4 have been reached.

6.2. Approach followed

The approach followed for phase 3 of the characterization tests is to operate the C1 reactor in continuous mode in order to reach the targeted suspended matter concentration, and then to reach the steady state at 10 days HRT. For that steady state, a mass balance can be calculated on the bioreactor, as per the following equation

Solid&liquid feed input -> reactor content (liquid+solid) + gas output + filtrate output + reactor samplings/bleedings

This mass balance can be drawn at overall level and for the main chemical elements (by order of precedence C balance, then N balance, then O balance and H balance), and be related to the operational parameters of the C1 unit, which will provide a first set of equations for the knowledge model.

The different parameters to be recorded during the tests have been grouped in three categories by order of priority, as follows :

Priority 1 (high priority): all the data necessary for the characterization of compartment 1 and the long term operation of C1 in the MPP integrated loop including operating parameters measured online (like the pH, the temperature, the pressure, the gas composition in CO₂ and CH₄), and parameters measured offline (like the sterility checks of the filtrate output, the VFAs, the dry matter, the COD, the pH, the electroconductivity, the bacterial counts)

Priority 2 (medium priority): all the data necessary for computing mass balance on C1 bioreactor as per the hereabove equation (total and soluble nitrogen, soluble COD, ammonium, organic elemental composition, ashes, gas composition in H₂, H₂S and O₂)

Priority 3 (low priority) : the remaining parameters used to refine the models later on (particles size, capillary suction time, proteins in total and soluble fractions, alkalinity, mineral elemental composition, gas contaminants)

The steady state indicators of the bioreactor are measured more frequently than once per HRT in order to detect the establishment of the steady state, these are:

- Dry matter content
- Total Chemical Oxygen Demand (COD_{tot})
- CO₂ production rate
- VFA production rate and, for information, the ratio between the various VFAs compounds

When the steady state is reached, all the parameters of priority 1 and 2 are measured twice, separated by one HRT

Document Identification : COO3 – WP94.2 – Sampling and Analysis Protocols Issue 2 - Test with 10 days hydraulic residence time	Type	Chrono	Issue	Page : 10 / 47
	TN	94.22	0	

6.3. *Applicable requirements*

See the test protocol AD9.

6.4. *Recall of test sequence*

The characterization tests sequence can be summarized as follows :

Phase 1 : maintenance of the inoculum, 20 days HRT

Phase 2 : ramp-up of the culture in the C1 bioreactor up to continuous conditions, HRT evolving from 20 days up to 10 days, reaching a dry matter content between 40g/L and 70g/L

Phase 3 : 10 days liquid residence time test

Phase 4 : 7 days liquid residence time test

Phase 5 : 13 days liquid residence time test

Phase 6 : 5 days liquid residence time test

6.5. *Success/failure criteria for the 10 days HRT*

The characterization tests are considered successful if the specified conditions have been maintained and the expected data have been collected as per the sampling plan.

As per AD9, the steady state has been successfully reached if the following parameters are detected stable within the following deviations around their average on a period of 3 HRTs:

- Dry matter content : 15%
- Total Chemical Oxygen Demand (COD_{tot}) : 15%
- CO₂ production rate : 20%
- VFA production rate : 20%

The degree of closure of the mass balance is also considered as a success criterion for the sampling and analyses activities. The ratio of measured/calculated output mass by the measured/calculated input masses on the bioreactor should be higher than 90%.

Similar success criteria on C mass balance and N mass balance closures are defined and set to 80%.

For more details see the test protocol AD9.

Document Identification : COO3 – WP94.2 – Sampling and Analysis Protocols Issue 2 - Test with 10 days hydraulic residence time	Type	Chrono	Issue	Page : 11 / 47
	TN	94.22	0	

7. Data collection plan – Sampling plan

7.1. Uncertainty acceptance level

The uncertainty budget has not been exhaustively assessed for all the measurement techniques to be implemented.

A general approach is to accept on all biological samples an uncertainty of 10% due to the natural variety present in the sample.

For three measurement techniques, the uncertainty was assessed, and the budget was calculated :

- pH see AD24 for the procedure details
- gas mass flow , see AD25 for the procedure details
- VFAs by Gas Chromatography , see AD26 for the procedure details

The calculated expanded uncertainties with a level of confidence of 95% are respectively ± 0.065 pH unit for pH and 2% for CH₄ gas mass flow.

For the VFAs measurement using the gas chromatography, the method reached an expanded uncertainty of 44% to 69% with a level of confidence of 95%, which led to the abandon of this technique for the alternative of VFAs measurements by HPLC that improved the repeatability and reduced uncertainties.

7.2. Measurement plan

The measurement plan as discussed among the partners of call off order 3 includes the following parameters (cf. RD6)

Phase	Physical or chemical or biological parameter
Liquid-solid phase	total liquid flow or volumes
	Dry matter
	ashes
	sample volume and weight
	CHONS total
	Minerals: P, Ca, Mg, Na, K, Si, S, Fe, Al, Ba, Cr, Cu, Mn, Ni, Sr, Zn, Mo, Ti, Be, V, Co, As, Se, Pd, Pb, Cd, Sn, Sb, W, Hg
	VFAs
	NH ₄ ⁺
	COD soluble
	COD total
	N total
	N soluble
	aerobic count
	anaerobic count
	EC
	pH
	Temperature
	Proteins (total and soluble)
	Fibers
	alkalinity
CST	
turbidity	

Document Identification : COO3 – WP94.2 – Sampling and Analysis Protocols Issue 2 - Test with 10 days hydraulic residence time	Type	Chrono	Issue	Page : 12 / 47
	TN	94.22	0	

Gas phase	Sterility of filtrate solution
	particles size
	total mass gas flow or volumes
	CO2 concentration
	CH4 concentration
	H2 concentration
	H2S concentration
	O2 concentration
	Sample volume
	Pressure
	temperature
	gas contaminants

7.3. Sampling techniques

The sampling techniques are used to recover gas and liquid/solid samples.

For gas samples, a dedicated circuit allows to continuously circulate, dry out and analyze the biogas for CO2 and CH4 assaying. A bypass line also allows to force the biogas from C1 bioreactor to a portable analyzer in order to make further assays (CH4, CO2, but also O2, H2 and H2S).

For liquid samples, various ports allow to bleed through manual valves the content of the tanks or circuits.

No continuous sampling of liquid/solid phase is planned.

When making a sampling, the first bled mL are discarded in order to take a sample that be representative of the sampling point. The discarded quantity is weighed. The sampled quantity is weighed as well.

The bleeding is a particular sampling of the liquid-solid phase that intends to regulate the dry matter content in a continuously agitated tank with perfusion. For clarity purposes, the bleeding is distinguished from the other samplings. The bled quantity is also weighed.

All samples are taken in aerobic conditions.

The samplings made on the filtrate circuits and tank, ie downstream the microfiltration membranes are made in sterile conditions, with a previous sterilization, in order to preserve the sterility of the filtrate circuits and to collect a non contaminated sample (see AD29).

7.4. Sampling ports location

For the liquid-solid and gas samplings, various ports can be used.

7.4.1. Liquid sampling port on Influent Tank

One sampling port is available on the lower part of the influent tank, controlled by HV_1000_07.


Document Identification : COO3 – WP94.2 – Sampling and Analysis Protocols Issue 2 - Test with 10 days hydraulic residence time	Type	Chrono	Issue	Page : 13 / 47
	TN	94.22	0	

Picture 1- influent tank sampling port

7.4.2. Liquid sampling port on Bioreactor Tank

One sampling port is used to make the bleedings and samplings of liquid-solid phase on the bioreactor : the lower side port controlled by HV_1007_02.

Picture 2 – bioreactor sampling port

7.4.3. Liquid sampling port on Effluent Tank

The steam sterilizable port HV_1204_02 can be used to take filtrate sample from the effluent tank VSL2_1204_01 and to proceed to harvest.

Another sterilizable port HV_1204_01 can be used to take smaller quantities of filtrate sample from the effluent tank VSL2_1204_01. It is sterilized applying an ethanol wick on the septum before introducing the syringe into the port.

For more details on the steps to follow for sterilization of these ports, see AD29.

Document Identification : COO3 – WP94.2 – Sampling and Analysis Protocols Issue 2 - Test with 10 days hydraulic residence time	Type	Chrono	Issue	Page : 14 / 47
	TN	94.22	0	


Picture 3 – effluent tank harvesting port

7.4.1. Liquid sampling port on filtrate line

The filtrate line downstream the ultrafiltration membranes and upstream the peristaltic pump PP_1202_01 is equipped with a steam sterilizable sampling port HV_1210_03, as shown in the following PID excerpt.


For more details on the steps to follow for sterilization of this port, see AD29.

Picture 4 - filtrate line sterilizable sampling port

Document Identification :
 COO3 – WP94.2 – Sampling and Analysis Protocols
 Issue 2 - Test with 10 days hydraulic residence time


Type	Chrono	Issue	Page : 15 / 47
TN	94.22	0	

Picture 5 – sampling ports on effluent tank and filtrate line

7.4.2. Gas sampling ports on the bioreactor

A bypass line allows to force the biogas from C1 bioreactor to a portable analyzer in order to make further assays (CH₄, CO₂, but also O₂, H₂ and H₂S) and to return it to the bioreactor through HV .

Picture 6 – gas sampling ports for portable gas analyzer

Document Identification : COO3 – WP94.2 – Sampling and Analysis Protocols Issue 2 - Test with 10 days hydraulic residence time	Type	Chrono	Issue	Page : 16 / 47
	TN	94.22	0	

7.5. Sampling recipients

For the bleeding, a clean 2L bottle is used to recover the necessary volume. This bottle can then be used as a dispenser to extract several samples if needed for further analyses.

For the analyses, the samples are collected in a clean plastic recipient, with a maximum volume of 100mL and a screwable lid. The exact volume recovered into this recipient can be lower than 100mL but should be traced every time a sample is taken, preferably recording the mass of the sample.

Then depending on the analyses, specific containers can be used that are adapted to the requirements of the analysis technique :

- for HPLC, 1,5mL single use glass vials are used
- for lyophilized samples, a 50mL single use Falcon tube is used

Some backup samples are taken in case analyses need to be repeated ; they are kept in the following recipients for conservation in frozen state at -20°C :

- for HPLC samples for VFA analysis, one 2mL plastic Eppendorf is used for the deproteinized fraction of the sample
- for soluble fraction : two single use 2mL plastic Eppendorf are used
- for total broth : two single use 2mL plastic Eppendorf are used
- for lyophilized samples : the remaining quantity inside the 50mL single use Falcon tube

7.6. Use of the volumes extracted during the bleeding

The quantity of liquid-solid phase recovered while making the bleeding can if needed be used as a sample to perform further analyses on the bioreactor content.

In that case, the volume of the bled bioreactor content has to be sufficient to provide the volume necessary for samples preparation.

7.7. Conditioning of the samples

Important remark : before taking any aliquot from a sample to perform the mentioned analyses, the full sample should be previously agitated so that it is homogeneous.

The following analyses are made on the raw liquid sample :

pH, EC, COD of total fraction, proteins of total fraction, dry weight, ashes, particle size, anaerobic and aerobic counting, nitrogen of total fraction, particles size.

The following analyses are made on a sample that has been previously centrifuged and filtered on a 0.22 micron filter :

Ammonium, nitrogen of soluble fraction, COD of soluble fraction, proteins of soluble fraction, alkalinity.

For the VFAs analyses, the samples are first treated to degrade their proteins, then centrifuged and filtered on a 0.22 micron filter.

The following analyses are made on a sample that was previously frozen at -80°C, lyophilized and homogenized by milling:

Minerals analysis, COHNS elemental analysis, fibers composition.

Document Identification : COO3 – WP94.2 – Sampling and Analysis Protocols Issue 2 - Test with 10 days hydraulic residence time	Type	Chrono	Issue	Page : 17 / 47
	TN	94.22	0	

The backup samples follow the same treatment as for the nominal analysis but they are then reserved for further analysis and stored in frozen state at -20°C.

The following table indicates the minimal volumes of sample needed to perform each analysis.

Item	Priority	Type of recipient	Preparation for analysis			Backup for later analysis		
			Raw sample	Filtered	Lyophilized sample	Volume/content	Type of recipient	
Bleeding volume	1	2L glass bottle	variable				1L to 2.5L	Glass bottle
pH/EC	1	100mL plastic container	20mL				N/A (fresh)	
Dry weight/ashes	1		60mL				?	
Total N	2		5mL				2x2mL broth	Eppendorf
Proteins total	3		2mL					
COD total	1		2mL					
Ammonium	2		5mL	3mL			2x2mL soluble fraction	Eppendorf
COD soluble	2							
Soluble nitrogen	2							
Soluble proteins	3		2mL	2mL				
Alkalinity	3		300mL	150mL			N/A	
VFA	1		5mL	4mL			2 mL deproteinized filtered broth	Eppendorf
CHONS	2		200mL		3g		10g-20g	50mL Falcon tube
Minerals	2				3g			
Fibres	3			300mL		5g		
TIC/TOC	2		12mL	10mL				
CST	3		TBD					
Anaerobic count	1		0,5mL				N/A(fresh)	
Aerobic count	1		0,5mL				N/A(fresh)	
Sterility check	1		10mL				N/A(fresh)	
Particle size	3		TBD					

Document Identification : COO3 – WP94.2 – Sampling and Analysis Protocols Issue 2 - Test with 10 days hydraulic residence time	Type	Chrono	Issue	Page : 18 / 47
	TN	94.22	0	

7.8. Labelling of samples

Samples are labelled in such a way that the label allows

- the identification of the originating process part (RC for bioreactor, ET for effluent tank, IT for influent tank, FL for filtration line)
- the identification of the date of the sampling
- the identification of the treatment implemented for the sample (RW raw sample, SF soluble fraction, DP deproteinized, LY lyophilized, DL dilution)

Additionally, the traceability to the record where the preparation steps of the sample are recorded is made on that record.

Example :

RC_dd.mm.yy_RW, ET_dd.mm.yy_DPDL20

The labelling is made by using an appropriate marker or by printing a sticker with the relevant information and protecting the inscriptions with adhesive tape.

7.9. Conservation of the samples :

In case an analysis cannot be performed on a fresh raw liquid sample, the sample has to be stored in the fridge at 4°C and analyzed within 8h.

In case an analysis cannot be performed on a filtered liquid sample, the sample has to be stored in the fridge at 4°C and analyzed within 48h ; for VFAs, the samples stored at -20°C can be analyzed within 2 months.

In case an analysis cannot be performed on a lyophilized sample, the sample has to be stored in dry conditions and analyzed within 6 months.

The samples that are taken as back-ups in order to perform or repeat some analyses are stored at -20°C and can be analyzed within 2 months.

The gas samples are not stored ; they should be analyzed on the spot.

Document Identification : COO3 – WP94.2 – Sampling and Analysis Protocols Issue 2 - Test with 10 days hydraulic residence time	Type	Chrono	Issue	Page : 19 / 47
	TN	94.22	0	

7.10. Sampling and analyses frequency on the liquid-solid phase

The following six tables summarize the circuit where to sample and the frequency at which the samples should be taken on the liquid-solid phase, when the bioreactor is operated with a HRT of 10days. The frequencies of operation are relative to the applicable HRT.

Remarks :

- the steady state indicators are identified by a yellow highlighting in the following tables
- the operations marked on week-ends can be performed indifferently on Saturdays on Sundays due to operating constraints
- for the transient phase analyses, the analyses that are not required for steady state monitoring might be adjusted to the contingencies of the laboratory
- some frequencies for priority 3 analyses still have to be decided
- analyses or preparations are marked in orange cells when they can be performed on a batch of several samples, collected at distinct dates
- the backup samples are all described as “storage in frozen or lyophilized state” in the preparations and samplings sections
- for the effluent tank, the CHONS analyses will be limited to TIC/TOC technique on the soluble fraction, while soluble nitrogen will allow N determination in the liquid sample
- for the filtrate line and effluent tank, as explained in detail in AD29, in order to validate the sterility of the filtrate line and of the effluent tank, the samplings will be made on a daily basis in order to make sure that sterility can be reached on a 4 weeks period, then the applicable frequency will be the one specified in the tables

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Table 3 - samplings and analyses on the C1 effluent tank in transient phase

Document Identification : COO3 – WP94.2 – Sampling and Analysis Protocols Issue 2 - Test with 10 days hydraulic residence time	Type	Chrono	Issue	Page : 23 / 47
	TN	94.22	0	

After establishment of the steady state phase :

C1 Process Operations		Steady State																				frequency per HRT in the planning	frequency specified per HRT
		Month n																					
		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20		
Bioreactor																							
samplings and preparations																							
1	Bleeding	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	8,5	10
2	Sampling for analysis of total fraction	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	6	2
3	Storage of total fraction in frozen state	•				•				•					•							2	2
4	Sampling for analysis of VFAs	•			•			•			•				•					•		3	2
5	Preparation of samples for VFA analysis	•			•			•			•				•					•		3	2
6	Storage of VFA samples in frozen state	•			•			•			•				•					•		3	2
7	Sampling for analysis of soluble fraction		•													•						1	1
8	Filtering of sample		•													•						1	1
9	Storage of filtered sample in frozen state		•													•						1	1
10	Lyophilization	•									•											1	1
11	storage in lyophilized state	•									•											1	1
analyses on total fraction																							
1	Dry weight	•		•				•		•				•		•						2,9	2
2	COD, total	•						•						•								1,4	2
3	Ashes	•								•												1,0	1
4	Nitrogen conc., total sample	•								•												1,0	1
5	Protein conc., total sample	•								•												1,0	1
6	pH (offline)		•					•							•							1,4	1
7	Electroconductivity		•					•							•							1,4	1
8	Microbial count (aerobic)			•													•					1,0	1
9	Microbial count (anaerobic)			•													•					1,0	1
10	Particle size			•													•					1,0 ?	
11	VIAMASS																					online	online
analyses on soluble fraction																							
12	VFA composition analysis (HPLC)							•												•		1,0	2
13	COD, soluble fraction		•												•							1,0	1
14	Nitrogen conc., soluble fraction		•												•							1,0	1
15	Ammonium conc., soluble fraction		•												•							1,0	1
16	Protein conc., soluble fraction		•												•							1,0	1
17	Alkalinity		•												•							1,0 ?	
analyses on lyophilized fraction																							
18	Fibres composition	•													•							1,0	1
19	Organic elements conc. (COHNS)	•													•							1,0	1
20	Inorganic elements conc.	•													•							1,0	1

Table 4 - samplings and analyses on the C1 bioreactor in steady state

C1 Process Operations		Steady State																				frequency per HRT in the planning	frequency specified per HRT
		Month n																					
		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20		
Effluent Tank																							
samplings and preparations																							
2	Sampling for analysis of total fraction	●	●		●	●		●	●		●	●		●	●		●	●		●	●	6,2	2
3	Storage of total fraction in frozen state	●			●			●			●			●			●			●		2,9	2
4	Sampling for analysis of VFAs	●			●			●			●			●			●			●		2,9	2
5	Preparation of samples for VFA analysis	●			●			●			●			●			●			●		2,9	2
6	Storage of VFA samples in frozen state	●			●			●			●			●			●			●		2,9	2
7	Sampling for analysis of soluble fraction		●														●					1,0	1
8	Filtering of sample		●														●					1,0	1
9	Storage of filtered sample in frozen state		●														●					1,0	1
10	Lyophilization																					0,0	N/A
11	storage in lyophilized state																					0,0	N/A
analyses on total fraction																							
1	Dry weight	●			●			●			●			●			●			●		2,9	2
2	COD, total	●			●			●			●			●			●			●		1,4	1,33
3	Ashes	●									●											1,0	1
4	Nitrogen conc., total sample	●									●											1,0	1
5	Protein conc., total sample	●									●											1,0	1
6	pH (offline)		●						●								●					1,4	1,33
7	Electroconductivity		●						●								●					1,4	1,33
24	Sterility control (aerobic)	●			●			●			●			●			●			●		2,9	3
25	Sterility control (anaerobic)	●			●			●			●			●			●			●		2,9	3
analyses on soluble fraction																							
12	VFA composition analysis (HPLC)							●												●		1,0	2
14	Nitrogen conc., soluble fraction		●														●					1,0	1
15	Ammonium conc., soluble fraction		●														●					1,0	1
16	Protein conc., soluble fraction		●														●					1,0	1
17	Alkalinity		●																			0,5	0,33
18	Fibres composition																					0,0	?
19	Organic elements conc. (COHNS)	●															●					1,0	1
20	Inorganic elements conc.	●															●					1,0	1
Filtrate Line																							
samplings and preparations																							
5	Sampling for analysis of total fraction	●			●			●			●			●			●			●		2,9	2
analyses on total fraction																							
3	COD, total	●			●			●			●			●			●			●		1,4	2
24	Sterility control (aerobic)	●			●			●			●			●			●			●		2,9	3
25	Sterility control (anaerobic)	●			●			●			●			●			●			●		2,9	3
10	Particle size																					0,0	?

Table 6 - samplings and analyses on the C1 effluent tank in steady state

Document Identification : COO3 – WP94.2 – Sampling and Analysis Protocols Issue 2 - Test with 10 days hydraulic residence time	Type	Chrono	Issue	Page : 27 / 47
	TN	94.22	0	

Concerning the timing to take into account for the obtention of the results, see the following table.

analyses on total fraction		Time between sampling and results
1	Dry weight	24h
2	COD, total	Immediate
3	Ashes	36h
4	Nitrogen conc., total sample	Immediate
5	Protein conc., total sample	Immediate
6	pH (offline)	Immediate
7	Electroconductivity	Immediate
8	Microbial count (aerobic)	3 to 7days
9	Microbial count (anaerobic)	3 days
10	Particle size	TBD
11	VIAMASS	immediate
analyses on soluble fraction		
12	VFA composition analysis (HPLC)	
13	COD, soluble fraction	immediate
14	Nitrogen conc., soluble fraction	immediate
15	Ammonium conc., soluble fraction	immediate
16	Protein conc., soluble fraction	1 day
17	Alkalinity	1 day
analyses on lyophilized fraction		
18	Fibres composition	1 month to 6 months
19	Organic elements conc. (COHNS)	1 month to 4 months
20	Inorganic elements conc.	1 month to 2 months
analyses on gas phase		
21	CO2 and CH4 production	Immediate
22	CO2, CH4, H2, H2S (offline)	Immediate
23	gas contaminants	2 months

Document Identification : COO3 – WP94.2 – Sampling and Analysis Protocols Issue 2 - Test with 10 days hydraulic residence time	Type	Chrono	Issue	Page : 28 / 47
	TN	94.22	0	

8. Resources specification for the tests

8.1. *Personnel: staff qualification and training needs*

The MPP technicians are qualified to operate the C1 compartment.

The MPP Analysis Technicians are qualified to perform the sampling operations and the MPP inhouse analyses (cf. appendix 1)

8.2. *Hardware: instruments, specific part, hardware for software operation*

C1 Hardware as described in AD6

The portable gas analyzer GA94 as described in RD4.

The gas chromatograph is used for VFA measurement. It is described in RD 3. Alternatively for VFA measurement, a HPLC is used. It is described in RD8.

The preparation of the samples is made with a lyophilizer and other equipments of common use in the Chemical Engineering Department..

8.3. *Software : verification of software, backup needs*

The software used is the Schneider Concept V2.6. for C1 control.

8.4. *Facilities : environmental needs, test conditions, interfaces needs, utilities needs*

All hardware involved in MPP utilities for C1 as specified in AD6 and AD5.

8.5. *Test conditions*

As specified by the test protocol AD9.

9. Measurement and data sampling

9.1. *Data logfile*

The samplings and analyses are performed routinely and are recorded in written on dedicated record sheet, internal or external to the MPP.

These raw data are then typed into the C1 database for analyses.

9.2. *Reporting of status for a test*

On a monthly basis, the Bioprocess Engineer or the Technical Manager reviews the raw data, checks the trends and spots the inconsistent values. An example of this database is given in appendix 8.

Document Identification : COO3 – WP94.2 – Sampling and Analysis Protocols Issue 2 - Test with 10 days hydraulic residence time	Type	Chrono	Issue	Page : 29 / 47
	TN	94.22	0	

At the end of the test phase (for example the 10 days HRT test), or at least every 3 months, a report is compiled with all the analyses results related to the same test phase and sent by the Technical Manager to the partners..

9.3. Deviations and non conformances

Records should be maintained for any deviation or non conformity from this Sampling and Analysis Protocol and included in the as-run procedures records or test batch records.

9.4. Record for the performance of samplings and analysis

9.4.1. Records of samplings

The samplings are recorded on the C1 follow-up sheet MPP-REC-10-1001 as per appendix 3.

9.4.2. Records of analyses

Various MPP records are used to trace the results of the analyses on C1 samples :

- For dry weight and ashes : MPP-REC-09-1005 (appendix 4)
- For Electroconductivity, NH₄, N total, Nsoluble, COD total, COD soluble : MPP-REC-09-1003 (appendix 5)
- For the CHNS analysis, the analysis report is given in appendix 6.
- For the minerals analysis, the report is given in appendix 7.

Document Identification : COO3 – WP94.2 – Sampling and Analysis Protocols Issue 2 - Test with 10 days hydraulic residence time	Type	Chrono	Issue	Page : 30 / 47
	TN	94.22	0	

10. Appendix 1 – approved MPP protocols for analyses on samples

Parameter analyzed	MPP reference	Principle of the analysis technique
electroconductivity	MPP-QCP-07-0008	potentiometric
pH	MPP-QCP-07-0015	potentiometric
	MPP-QCP-10-1001	
dry matter	MPP-QCP-07-0007	gravimetric (oven is maintained at 100°C in UAB)
ashes	MPP-QCP-07-0001	gravimetric (oven is held at 550°C in UAB)
VFA	MPP-QCP-10-0003	Gas chromatograph
NH4-N	MPP-QCP-07-0011	colorimetric
N total broth	MPP-QCP-12-0001	colorimetric
N soluble	MPP-QCP-07-0020	colorimetric
COD total	MPP-QCP-07-0005	colorimetric
COD soluble	MPP-QCP-07-0004	colorimetric
CO2,CH4	MPP-QCP-10-0002	Online gas analyzer
Protein (broth and soluble fraction)	MPP-QCP-12-1001	Lowry method adapted for Comp. I
H2, O2, H2S	MPP-UM-09-0012	Portable infrared gas analyzer
Anaerobic count	MPP-QCP-09-1001	Petri dish seeding, incubation and counting of colonies
Aerobic count	MPP-QCP-09-1001	Petri dish seeding, incubation and counting of colonies
Filtrate circuits sterility	MPP-QCP-12-1002	Filtering of a volume of liquid, Petri dish seeding, incubation and counting of colonies
VFAs with HPLC	MPP-QCP-11-0002	HPLC
Alkalinity	MPP-QCP-11-0001	Dosimetric method

Missing here, to be defined in an updated edition: particles size ; CST ; TIC/TOC

Document Identification : COO3 – WP94.2 – Sampling and Analysis Protocols Issue 2 - Test with 10 days hydraulic residence time	Type	Chrono	Issue	Page : 31 / 47
	TN	94.22	0	

11. Appendix 2 - External protocols implemented for analyses

11.1. Inorganic elemental composition (minerals) SAQEAt0001_00

SAQEAt0001_00_elements inorgànics_bior

11.2. C,H,N,S organic elemental composition : SAQAE0001_00_CHNS

SAQAE0001_00_CHNS_101208.pdf

Document Identification : COO3 – WP94.2 – Sampling and Analysis Protocols Issue 2 - Test with 10 days hydraulic residence time	Type	Chrono	Issue	Page : 32 / 47
	TN	94.22	0	

11.3. Fibers content : not available in English, added in Catalan language

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P4 = Pes del gresol i les cendres

OBSERVACIONS

Document Identification : COO3 – WP94.2 – Sampling and Analysis Protocols Issue 2 - Test with 10 days hydraulic residence time	Type	Chrono	Issue	Page : 33 / 47
	TN	94.22	0	

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Document Identification : COO3 – WP94.2 – Sampling and Analysis Protocols Issue 2 - Test with 10 days hydraulic residence time	Type	Chrono	Issue	Page : 34 / 47
	TN	94.22	0	

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MELISSA Pilot Plant

Document Identification : COO3 – WP94.2 – Sampling and Analysis Protocols Issue 2 - Test with 10 days hydraulic residence time	Type	Chrono	Issue	Page : 35 / 47
	TN	94.22	0	

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Document Identification : COO3 – WP94.2 – Sampling and Analysis Protocols Issue 2 - Test with 10 days hydraulic residence time	Type	Chrono	Issue	Page : 36 / 47
	TN	94.22	0	

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MELISSA Pilot Plant

Document Identification : COO3 – WP94.2 – Sampling and Analysis Protocols Issue 2 - Test with 10 days hydraulic residence time	Type	Chrono	Issue	Page : 37 / 47
	TN	94.22	0	

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MELiSSA Pilot Plant

Document Identification : COO3 – WP94.2 – Sampling and Analysis Protocols Issue 2 - Test with 10 days hydraulic residence time	Type	Chrono	Issue	Page : 38 / 47
	TN	94.22	0	

cover them. Leave them in a soak for 2-3', remove them and gently squeeze them to eliminate the excess of acetone.

Separate the bags and let them dry in the air. Complete the drying in an oven at 105°C. Remove the bags from the oven, introduce them in a desiccator until they reach the T^a environment and weigh them (P3). Incinerate the bag with the residue in a crucible previously weighed, cool again in a desiccator and weigh.

Cálculos. % Fibra Neutro detergente (libre de cenizas) =
$$\frac{(P4 - (P1 \times C2)) \times 100}{P2 \times MS}$$

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MELiSSA Pilot Plant

Document Identification : COO3 – WP94.2 – Sampling and Analysis Protocols Issue 2 - Test with 10 days hydraulic residence time	Type	Chrono	Issue	Page : 39 / 47
	TN	94.22	0	

12. Appendix 3 – example of MPP-REC-10-1001 for the follow-up of the Compartment 1

	MELISSA Pilot Plant			
Document Identification : <i>C1 in Bioreactor Follow up Record Sheet</i>	Type MPP-REC	Reference 10-1001 (0)	Chrono --	Page : 1 / 1

Analyst	
Date	
Hour	

Checked by:	
Initials	
Date	

HMI - C1 BIOREACTOR						
Emergency button	(on/off)		Level	LT_1010_01	(L)	
Agitation	BLE_1012_01	(on/off)	Headspace pressure	PT_1009_01	(mbar)	
On-line pH 2	AT_1011_02		Temperature	TT_1008_01	(°C)	
C1 ROOM GENERAL						
On-line pH 2	AT_1011_02		ACID bottle:	VSSL_1011_01		
Cooler	HX_1102_01	(on/off)	Observed Level		(mL)	
Condensates pump	PP_1102_01	(on/off)	HCl 3M added volume		(mL)	
Hot bath:	VSSL_1008_01		BASE bottle:	VSSL_1011_02		
Filled with water?	(yes/no)		Observed Level		(mL)	
Temperature	(°C)		NaOH 3M added volume		(mL)	
GAS LOOP						
N2 supply	HPCV_1052_01	(mbar)	Mass Flow Controller	FRC_1052_01	(%)	
GAS ANALYSIS:						
CH4	(%)		O2	(%)		
CO2	(%)		H2	(ppm)		
			H2S	(ppm)		
LIQUID LOOP						
Feeding: (On Mondays)			Decantation: (On Fridays)			
Volume of FILTRATE dumped	(L)		Volume of RC extracted (EFFLUENT)	(L)		
Weight of FEED	(Kg)		Samples: (total volume in mL)			
FEED Lot number	#		REACTOR CONTENT		pH off-line	
Volume of FEED + WATER	(L)		FILTRATE		pH off-line	
			FEED		pH off-line	
MELISSA MAINTENANCE FEED						
On Fridays:						
Bag of FEED stored at 3-5°C	(yes/no)		N° of bags left		#	
Weight of Bag	(Kg)					
FEED Lot number	#					
REMARKS						
.						

This document is confidential property of the MELiSSA partners and shall not be used, duplicated, modified or transmitted without their authorization Page :

Memorandum of Understanding ESTEC 4000 100 293/10/NL/PA

39/47

MELiSSA Pilot Plant

Document Identification : COO3 – WP94.2 – Sampling and Analysis Protocols Issue 2 - Test with 10 days hydraulic residence time	Type	Chrono	Issue	Page : 40 / 47
	TN	94.22	0	

13. Appendix 4 – example of record for dry weight and ashes : MPP-REC-09-1005

	MELiSSA Pilot Plant			
Document Identification : C1 samples Dry weight and ashes Analysis Record Sheet	Type MPP-REC	Reference 09-1005 (0)	Chrono --	Page : 1 / 1

Name:							
Date of analysis:							
Project name:						Ec:	
Sample name/date of sampling:							
Crucible nº	Volume (mL)	x1	x2	x3		DW	Ashes
	20						
	20						
Done by: (Initials)						Average:	
Project name:						Ec:	
Sample name/date of sampling:							
Crucible nº	Volume (mL)	x1	x2	x3		DW	Ashes
Done by: (Initials)						Average:	
Project name:						Ec:	
Sample name/date of sampling:							
Crucible nº	Volume (mL)	x1	x2	x3		DW	Ashes
Done by: (Initials)						Average:	
Project name:						Ec:	
Sample name/date of sampling:							
Crucible nº	Volume (mL)	x1	x2	x3		DW	Ashes
Done by: (Initials)						Average:	
Remarks:						Checked by:	
						Initials	Date

MELiSSA Pilot Plant

Document Identification : COO3 – WP94.2 – Sampling and Analysis Protocols Issue 2 - Test with 10 days hydraulic residence time	Type	Chrono	Issue	Page : 41 / 47
	TN	94.22	0	

14. Appendix 5 – example of record for Electroconductivity, NH4, N total, COD total, COD soluble : MPP-REC-09-1004 and MPP-REC-09-1003

		MELiSSA Pilot Plant		 de Barcelona	
Document Identification : C1 Filtrate samples Dr. Lange Analysis Record Sheet		Type MPP-REC	Reference 09-1004 (0)	Chrono ---	Page : 1 / 1
Name:		Project name:			
Date of analysis:		Sample name/date of sampling:			
Hour:					
Sample filtrated? Yes/No		No			
Parameter:		COD total (ppm)			
Equipment/ Dr Lange kit code:					
Initial time:					
End time					
	Dilution	Raw result	Corrected result		
1					
2					
3					
4					
Done by (initials):		Average:			
Sample filtrated? Yes/No					
Parameter:					
Equipment/ Dr Lange kit code:					
Initial time:					
End time					
	Dilution	Raw result	Corrected result		
1					
2					
3					
4					
Done by (initials):		Average:			
Sample filtrated? Yes/No					
Parameter:					
Equipment/ Dr Lange kit code:					
Initial time:					
End time					
	Dilution	Raw result	Corrected result		
1					
2					
3					
4					
Done by (initials):		Average:			
Sample filtrated? Yes/No					
Parameter:					
Equipment/ Dr Lange kit code:					
Initial time:					
End time					
	Dilution	Raw result	Corrected result		
1					
2					
3					
4					
Done by (initials):		Average:			
Remarks:			Checked by:		
			Initials		Date

MELiSSA Pilot Plant

Document Identification : COO3 – WP94.2 – Sampling and Analysis Protocols Issue 2 - Test with 10 days hydraulic residence time	Type	Chrono	Issue	Page : 42 / 47
	TN	94.22	0	

	MELiSSA Pilot Plant			
Document Identification :	Type	Reference	Chrono	Page : 1 / 1
Record for C1 samples analysis	MPP-REC	09-1003 (1)	--	

Analyst name:		Supervised by:	
Date of analysis:		Date:	

Sample name	Sample date	Tube number	EC (Sample not filtered)	COD total (Sample not filtered)				COD soluble				NH ₄ ⁺				N total				
				Kit Code:				Kit Code:				Kit Code:				Kit Code:				
				Kit Batch n°	Dilution	Raw result	Final Result	Kit Batch n°	Dilution	Raw result	Final Result	Kit Batch n°	Dilution	Raw result	Final Result	Kit Batch n°	Dilution	Raw result	Final Result	
MPPC1...	dd/mm/aa	#	mS/cm			mg/L	mg/L			mg/L	mg/L			mg/L	mg/L			mg/L	mg/L	

Analyses to be done per type of sample:

MPPC1BF	yes	yes	yes	yes
MPPC1AF	yes	yes	yes	yes
MPPC1 Filtrate	yes	yes	no	no
MPPC1Feed#...	yes	yes	yes	yes

Remarks:

MELiSSA Pilot Plant

Document Identification : COO3 – WP94.2 – Sampling and Analysis Protocols Issue 2 - Test with 10 days hydraulic residence time	Type	Chrono	Issue	Page : 43 / 47
	TN	94.22	0	

15. Appendix 6 – example of analysis report for the CHNS analysis

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

1 de 2

Document Identification : COO3 – WP94.2 – Sampling and Analysis Protocols Issue 2 - Test with 10 days hydraulic residence time	Type	Chrono	Issue	Page : 44 / 47
	TN	94.22	0	

16. Appendix 7 – example of report for minerals analysis

Muestra	MPP-C1-RC enero	MPP-C1-RC 042010
Codigo SAQ	10EAt153/001	10EAt153/004
	Contingut	Contingut
ug / g Be	< 0.5	< 0.5
% Na (p/p)	6.7	6.2
% Mg (p/p)	0.18	< 0.05
ug / g Al	251	80
% Si (p/p)	< 0.1	< 0.1
% P (p/p)	0.17	0.25
% S (p/p)	0.17	0.25
% K (p/p)	2.2	3.7
% Ca (p/p)	0.43	0.55
ug / g Ti	41	22
ug / g V	0.97	0.58
ug / g Cr	76	32
ug / g Mn	56	41
ug / g Fe	517	177
ug / g Co	0.63	0.55
ug / g Ni	129	35
ug / g Cu	16	15
ug / g Zn	131	146
ug / g As	< 0.5	< 0.5
ug / g Se	< 0.5	< 0.5
ug / g Sr	41	42
ug / g Mo	27	5.9
ug / g Pd	< 0.5	< 0.5
ug / g Cd	0.19	0.22
ug / g Sn	< 0.5	< 0.5
ug / g Sb	< 0.5	< 0.5
ug / g Ba	13	13

This document is confidential property of the MELiSSA partners and shall not be used, duplicated, modified or transmitted without their authorization Page :

MELiSSA Pilot Plant

Document Identification : COO3 – WP94.2 – Sampling and Analysis Protocols Issue 2 - Test with 10 days hydraulic residence time	Type	Chrono	Issue	Page : 45 / 47
	TN	94.22	0	

ug / g W	1.0	0.58
ug / g Pb	2.0	0.93
ug / g Hg	0.19	0.063

MELISSA Pilot Plant

Document Identification : COO3 – WP94.2 – Sampling and Analysis Protocols Issue 2 - Test with 10 days hydraulic residence time	Type	Chrono	Issue	Page : 46 / 47
	TN	94.22	0	

17. Appendix 8 – example of the database used for C1 results of analysis and operational parameters

Time			MPP RECORD SHEET			HMI - C1 BIOREACTOR & INFLUENT TANK -						
			Analyst	Reference	Chrono	Emergency button	Agitation	Level Alarm	Recirculation pump	Level	Headspace pressure	Temperature
Date	hour	day	(name)	Page number	(on/off)	(on/off)	(on/off)	(on/off)	(L)	(mbar)	(°C)	
11/03/2009	16:30	0.0	NM 09-1006 (0)	001	OFF	ON	ON	ON	55.8	-	45.0	
11/03/2009	17:30	0.0	NM 09-1006 (0)	002	OFF	ON	ON	ON	56.0	6.6	53.3	
11/03/2009	19:00	0.1	NM 09-1006 (0)	003	OFF	ON	ON	ON	56.0	3.9	55.8	
12/03/2009	10:20	0.7	NM 09-1006 (0)	004	OFF	ON	OFF	ON	54.5	71.1	57.0	
12/03/2009	12:00	0.8	NM 09-1006 (0)	005	OFF	ON	OFF	ON	54.4	50.4	54.7	
12/03/2009	20:00	1.1	NM 09-1006 (0)	006	OFF	ON	OFF	ON	54.3	65.5	55.0	
13/03/2009	19:30	2.1	NM 09-1006 (0)	007	OFF	ON	OFF	ON	54.3	55.1	55.2	
16/03/2009	11:15	4.8	NM 09-1006 (0)	008	OFF	ON	OFF	ON	43.9	86.9	55.2	
16/03/2009	16:21	5.0	NM 09-1006 (0)	009	OFF	ON	OFF	ON	43.6	98.4	55.5	
16/03/2009	17:00	5.0	NM 09-1006 (0)	010	OFF	ON	OFF	ON	-	28.7	-	
17/03/2009	11:00	5.8	NM 09-1006 (0)	011	OFF	ON	OFF	ON	53.3	92.1	50.6	
17/03/2009	18:00	6.1	NM 09-1006 (0)	012	OFF	ON	OFF	ON	53.4	79.9	54.7	
18/03/2009	13:00	6.9	NM 09-1006 (0)	013	OFF	ON	OFF	OFF	53.9	135.6	54.2	

Time			C1 ROOM GENERAL										GAS LOOP										
			HOT BATH					Acid					Base					N2 Supply HPCV_1003_03	Mass Flow Controller FRC_1003_01	CH4	CO2	O2	H2
Date	hour	day	Recirculation pump running	On-line pH	Level	Temperature hot bath	Bottle level	Addition of HCl 3M to the bottle	Addition to the reactor absolute	Addition to the reactor cumulative	Bottle level	Addition of NaOH 3M to the bottle	Addition to the reactor absolute	Addition to the reactor cumulative	(mbar)	(%)	(%)						
11/03/2009	16:30	0.0	YES	5.06	YES	63	1000	-	0	0	1000	-	0	0	-	0	0	-	-	-	-	-	-
11/03/2009	17:30	0.0	YES	5.47	NO	63	1000	-	0	0	500	-	500	500	-	-	-	-	-	-	-	-	-
11/03/2009	19:00	0.1	YES	5.38	-	-	1000	-	0	0	300	-	200	700	-	-	-	-	-	-	-	-	-
12/03/2009	10:20	0.7	YES	5.28	NO	60	1000	-	0	0	100	-	200	900	-	-	-	-	-	-	-	-	-
12/03/2009	12:00	0.8	YES	5.29	NO	60	1000	-	0	0	100	900	0	900	-	-	-	-	-	-	-	-	-
12/03/2009	20:00	1.1	YES	5.28	NO	60	1000	-	0	0	1000	-	0	900	-	-	0	0.4	0.5	45	3	-	
13/03/2009	19:30	2.1	YES	5.27	NO	60	1000	-	0	0	1000	-	0	900	-	-	0	1.2	2.2	>1000	-	-	
16/03/2009	11:15	4.8	YES	5.26	YES	60	1000	-	0	0	1000	-	0	900	-	-	0	2	4.2	>1000	-	-	
16/03/2009	16:21	5.0	YES	0	-	-	1000	-	0	0	1000	-	0	900	-	-	0	1.7	1.1	>1000	218	-	
16/03/2009	17:00	5.0	0	0	-	-	-	-	0	0	-	-	0	900	-	-	0	0.1	0.8	219	-	-	
17/03/2009	11:00	5.8	YES	5.35	YES	60	1000	-	0	0	1000	-	0	900	-	-	-	-	-	-	-	-	
17/03/2009	18:00	6.1	YES	5.37	NO	60	1000	-	0	0	1000	-	0	900	-	-	0	0.2	0.8	79	4	-	
18/03/2009	13:00	6.9	NO	5.33	NO	60	1000	-	0	0	1000	-	0	900	-	-	-	-	-	-	-	-	
19/03/2009	19:00	8.1	YES	5.30	YES	60	1000	-	0	0	800	-	200	1100	-	-	-	-	-	-	-	-	
19/03/2009	19:30	8.1	YES	5.30	YES	60	1000	-	0	0	800	-	0	1100	500	10	-	-	-	-	-	-	
20/03/2009	11:00	8.8	YES	5.28	YES	58.8	1000	-	0	0	800	-	0	1100	500	10	0	0.7	7.8	129	-	-	

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Document Identification : COO3 – WP94.2 – Sampling and Analysis Protocols Issue 2 - Test with 10 days hydraulic residence time	Type	Chrono	Issue	Page : 47 / 47
	TN	94.22	0	

Time			LIQUID LOOP											
			FEED			INFLUENT		REACTOR CONTENT		EFFLUENT		FILTRATE/SUPERNATANT		
			Lot	Weight	Total Volume	Sample volume	pH off-line	Total Volume Feed + Solids into C1	Sample volume	pH off-line	Volume	Volume not sedimented or disposed	Volume	Sample volume
Date	hour	day	#	(Kg)	(mL)	(mL)		(mL)		(mL)	(mL)	(mL)	(mL)	
11/03/2009	16:30	0.0	-	0	0	0	0	100	5.03	0		0	0	
11/03/2009	17:30	0.0	-	0	0	0	0	0		0		0	0	
11/03/2009	19:00	0.1	-	0	0	0	0	0		0		0	0	
12/03/2009	10:20	0.7	-	0	0	0	0	0		2000	2000	0	0	
12/03/2009	12:00	0.8	-	0	0	0	0	0		0		0	0	
12/03/2009	20:00	1.1	-	0	0	0	0	0		0		0	0	
13/03/2009	19:30	2.1	-	0	0	0	0	0		0		0	0	
16/03/2009	11:15	4.8	-	0	0	0	0	0		10000		0	0	
16/03/2009	16:21	5.0	-	0	0	0	0	0		0		0	0	
16/03/2009	17:00	5.0	-	0	0	0	0	0		0		0	0	
17/03/2009	11:00	5.8	14	3	7000	80	6.26	10000	5.17	0		7000	80	5.19
17/03/2009	18:00	6.1	-	0	0	0		150	5.26	0		0	0	
18/03/2009	13:00	6.9	-	0	0	0		0		0		0	0	
19/03/2009	19:00	8.1	-	0	0	0		0		0		0	0	
19/03/2009	19:30	8.1	-	0	0	0		0		0		0	0	
20/03/2009	11:00	8.8	-	0	0	0		0		17500		0	0	

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Codi S.A.Q.: 8AE-024, 8AE-067 i 8AE-068 respectivament.